


C300 Corpus 3G



*** Indicates a Required Field**

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* I) Knee to Foot: _____

Please send order/quote to fax#: +61 2 8883 2356

Email to: sales.australia@permobil.com

Prices effective January 16, 2012.

Permobil Australia Pty Ltd.
Unit 312/ 14-16 Lexington Dr,
Norwest Business Park
Bella Vista
NSW 2153
Australia
www.permobil.com



Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I105033-99-0	C300 Base Corpus 3G R-net - C300 PS2 <i>Price Includes: FWD C300 Power Base, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 8kmh, Weight Capacity = 136kg, STFH = 17.5".</i>	
<input type="checkbox"/> I10637	10kmh High Speed Upgrade Package - C300/K300 <i>Price Includes: ESP Module, High Speed 10kmh Motor Package, and 120A Power Module. Available with R-net Electronics only.</i>	

CHOOSE DRIVE TIRES:

Part Number	Description	Price AUD
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge

OPTIONAL DRIVE TIRES:

Part Number	Description	Price AUD
<input type="checkbox"/> I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels.</i>	

CHOOSE SHROUD COLOUR:

Part Number	Description	Price AUD
<input type="checkbox"/> I105761-99-0	C300/K300 Shroud - Cobalt Blue (New Style!)	No Charge
<input type="checkbox"/> I105762-99-0	C300/K300 Shroud - Purple Rain (New Style!)	No Charge
<input type="checkbox"/> I105763-99-0	C300/K300 Shroud - Onyx Black (New Style!)	No Charge
<input type="checkbox"/> I105764-99-0	C300/K300 Shroud - Phantom Gray (New Style!)	No Charge
<input type="checkbox"/> I105765-99-0	C300/K300 Shroud - Popstar Pink (New Style!)	No Charge
<input type="checkbox"/> I105766-99-0	C300/K300 Shroud - Volcano Red (New Style!)	No Charge

CHOOSE HUBCAP COLOUR:

Part Number	Description	Price AUD
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Color	No Charge
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver	No Charge

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price AUD
<input type="checkbox"/> I105138-99-0	Power Adjustable Seat Height - 8" Travel	
<input type="checkbox"/> I105421-99-0	Seat Tube Fixed	No Charge

MUST CHOOSE BATTERIES:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> IM34 SLDG	Batteries, Grp 34 (60Ah) Sealed LA/Gel, Installed	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10789	Seat Stabilizer Option	
<input type="checkbox"/> I0036-7002A	Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	
<input type="checkbox"/> I105715-99-0	LED Light Kit for C/K300 w/ R-net (and New Shroud)	



Base Options

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I102768-99-0	Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent & a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics).</i>	
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus 3G Seat <i>Accommodates Pulmonetics LTV Series and similar sized Laptop Ventilators. The tray is 13"T, adjustable in width from 8" - 13", and adjustable in depth from 2" - 6.5". Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option.</i>	
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray	
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	

Permolock

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a Licensed Vehicle Modifier.</i>	
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> ID51185	R-net Remote Color JS w/ Mono Jacks, No Lights	
<input type="checkbox"/> ID51188	R-net Remote Color JS w/ Toggle Switches, No Lghts	
<input type="checkbox"/> ID51189	R-net Remote Color JS w/ Spd Pot/Toggle, No Lights	
<input type="checkbox"/> ID51153	R-net Remote Color JS w/ Mono Jacks, Lights	
<input type="checkbox"/> ID51149	R-net Remote Color JS w/ Toggle Switches, Lights	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Spd Pot/Toggle, Lights	
<input type="checkbox"/> I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimble <i>Heavy-Duty modification enables gimble to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick Module	No Charge

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	



Joystick Options

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price AUD
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2	
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2	
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net <i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	
<input type="checkbox"/> I10793	Stealth Gatlin Mid Line Mount - 8 Position <i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	
<input type="checkbox"/> I10794	Stealth Gatlin Mid Line Mount - 2 Position <i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	
<input type="checkbox"/> I10803	Stealth Gatlin Bracket for R-net Joystick	
<input type="checkbox"/> I10804	Stealth Gatlin Bracket for Compact Joystick	
<input type="checkbox"/> I10805	Stealth Gatlin Bracket for Mini Joystick	
<input type="checkbox"/> I10806	Stealth Gatlin Bracket for Omni	
<input type="checkbox"/> I10802	Stealth Height Adjustable Mnt (ARM250-P) <i>Provides 1.75" of height adjustability to all Permobil joystick mounts and Stealth Gatlin Mid Line Mounts.</i>	
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus 3G	
<input type="checkbox"/> I105408-99-0	Attendant Control R-net - Corpus 3G	
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	

Seating System

CHOOSE SEAT MOUNTING PLATE TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I105056-99-0	Standard (0°) Tilt Mounting Plate <i>Range of Tilt will be 0° to +50°.</i>	No Charge



Seating System

CHOOSE SEAT MOUNTING PLATE TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I105068-99-0	5° Anterior Tilt Mounting Plate <i>Range of Tilt will be -5° to +45°.</i>	

CHOOSE SEATING SYSTEM:

Part Number	Description	Price AUD
PLEASE NOTE: <i>Weight Limit 136kg. For the Corpus 3G Seat on the C300, there is a half speed limitation when the Total Backrest Angle (TBA) > 150°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 45° Tilt + 100° Recline = 145° TBA (Full Speed) Example 2: 35° Tilt + 120° Recline = 155° TBA (Half Speed)</i>		
<input type="checkbox"/> I10647	C3G Seat w/ No Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>	
<input type="checkbox"/> I10648	C3G Seat w/ No Tilt & 175° Power Recline	
<input type="checkbox"/> I10649	C3G Seat w/ 50° DuraCore Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>	
<input type="checkbox"/> I10650	C3G Seat w/ 50° DuraCore Tilt & 175° Power Recline	

Seating System Options

<input checked="" type="checkbox"/> I10585	Seat Function Control Kit for R-net <i>Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	
--	---	--

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price AUD
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest	No Charge
<input type="checkbox"/> I105026-99-0	Mount ICS Switchbox on Right	No Charge
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	
<input type="checkbox"/> I105027-99-0	Mount ICS Switchbox on Left	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	
<input type="checkbox"/> I104769-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	
<input type="checkbox"/> I10750	Memory Seat Program, Corpus 3G w R-net	

Backrest Options



Backrest Options

CHOOSE FIXED OR SLIDING BACKREST KIT:

<input type="checkbox"/> Part Number	Description	Price AUD
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear	

CHOOSE BACKREST TYPE:

<input type="checkbox"/> Part Number	Description	Price AUD
BACKREST HEIGHT NOTE:		
<i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°.</i>		
<input type="checkbox"/> I105113-99-0	Corpus 3G Ergo Backrest - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	
<input type="checkbox"/> I105114-99-0	Corpus 3G Ergo Backrest - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	
<input type="checkbox"/> I105109-99-0	Corpus 3G Ergo Backrest - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	

PLEASE NOTE:

All taller backrests listed below include: lateral wedges, lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided.

<input type="checkbox"/> I105110-99-0	Corpus 3G Ergo Back -14"W x 23 - 28"T (C3GERGBCA)	
<input type="checkbox"/> I105111-99-0	Corpus 3G Ergo Back -16"W x 23 - 28"T (C3GERGBCA)	
<input type="checkbox"/> I105112-99-0	Corpus 3G Ergo Back -18"W x 23 - 28"T (C3GERGBCA)	
	Set Backrest at 23" _____	
	Set Backrest at 24" _____	
	Set Backrest at 25" _____	
	Set Backrest at 26" _____	
	Set Backrest at 27" _____	
	Set Backrest at 28" _____	

CHOOSE BACKREST CUSHION:

<input type="checkbox"/> Part Number	Description	Price AUD
<input type="checkbox"/> I10659	Leatherette Cushion - Corpus 3G Ergo Back	No Charge
<input type="checkbox"/> I10660	Black Mesh Cushion - Corpus 3G Ergo Back	
<input type="checkbox"/> I10661	Gray Mesh Cushion - Corpus 3G Ergo Back	
<input type="checkbox"/> I105117-99-0	Omit Backrest Cushion	No Charge

OPTIONAL ITEMS:

<input type="checkbox"/> Part Number	Description	Price AUD
<input type="checkbox"/> I10726	Tubular Backrest Frame - Corpus 3G (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Infinitely adjustable in width from 14" to 20". Height of canes is 21" - 23" measured from the seat pan to the top of the cane. Provides an infinitely adjustable 1"-4" gap between the backrest hinge and the cane. Depending on the aftermarket back, some seat depth may be lost when mounting this system. This is a custom option. Please also indicate the length of the UniTrack linkage disks below:</i> 1" UniTrack Link Kit, C3G Tubular Back: _____ 2" UniTrack Link Kit, C3G Tubular Back: _____ 3" UniTrack Link Kit, C3G Tubular Back: _____	
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus 3G <i>Custom Option.</i>	
<input type="checkbox"/> I10786	Tailored Backrest Cushion - Corpus 3G <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. You must also select a backrest cushion (I10659, I10660 or I10661) when ordering this custom option.</i>	



Seat Sizing

CHOOSE SEAT FRAME AND WIDTH SETTING:

Part Number	Description	
PLEASE NOTE:		
<i>Narrow Seat Frames can be adjusted to 17"W, 19"W, or 21"W without any additional parts.</i>		
<i>Wide Seat Frames can be adjusted to 19"W, 21"W, or 23"W without any additional parts.</i>		
<i>Please select both the Frame Kit and Width Setting below.</i>		
<input type="checkbox"/> I10673	Narrow Corpus 3G Seat Frame Kit (17", 19", 21")	No Charge
<input type="checkbox"/> I10674	Wide Corpus 3G Seat Frame Kit (19", 21", 23")	No Charge

Set Seat Width at 17": _____
 Set Seat Width at 19": _____
 Set Seat Width at 21": _____
 Set Seat Width at 23": _____

CHOOSE SEAT DEPTH:

Part Number	Description	
PLEASE NOTE:		
<i>Odd numbered seat depths can be made 1" shorter without any additional parts. For example, a 19"D Seat can also be set at 18"D.</i>		
<input type="checkbox"/> I105040-99-0	Seat Depth 14"	No Charge
<input type="checkbox"/> I105457-99-0	Seat Depth 15"	No Charge
<input type="checkbox"/> I105095-99-0	Seat Depth 16"	No Charge
<input type="checkbox"/> I105458-99-0	Seat Depth 17"	No Charge
<input type="checkbox"/> I105042-99-0	Seat Depth 18"	No Charge
<input type="checkbox"/> I105459-99-0	Seat Depth 19"	No Charge
<input type="checkbox"/> I105045-99-0	Seat Depth 20"	No Charge
<input type="checkbox"/> I105460-99-0	Seat Depth 21"	No Charge
<input type="checkbox"/> I105047-99-0	Seat Depth 22"	No Charge

CHOOSE SEAT CUSHION:

Part Number	Description	
NOTE:		
<i>For seat widths of 23" and seat depths of 22", a Corpus II Ergonomic Seat Cushion will be provided.</i>		
<input type="checkbox"/> I10651	Corpus 3G Ergo Seat Cush - Leatherette (C3GERGSCA)	
<input type="checkbox"/> I10652	Corpus 3G Ergo Seat - Black Mesh (C3GERGSCA)	
<input type="checkbox"/> I10653	Corpus 3G Ergo Seat Cush - Gray Mesh (C3GERGSCA)	
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge

Armrest Options

CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	
PLEASE NOTE:		
<i>Narrow Armbar Kit can be adjusted to 15"W, 17"W or 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted to 19"W, 21"W, or 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>		
<input type="checkbox"/> I105440-99-0	Narrow Armbar Kit, Adjustable 15", 17", 19"W	No Charge
<input type="checkbox"/> I105441-99-0	Wide Armbar Kit, Adjustable 19", 21", 23"W	No Charge

Set Armbar at 15"W Between Armrests: _____
 Set Armbar at 17"W Between Armrests: _____
 Set Armbar at 19"W Between Armrests: _____
 Set Armbar at 21"W Between Armrests: _____
 Set Armbar at 23"W Between Armrests: _____



Armrest Options

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	Price AUD
<input type="checkbox"/> I105081-99-0	Height Adjustable Amrests	No Charge
<input type="checkbox"/> I105069-99-0	Height Adjustable Amrests - NC <i>No charge when ordering power tilt and/or recline.</i>	

SET ARMREST HEIGHT AT: _____

Corpus 3G Armrests are infinitely adjustable from 7" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.

CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	Price AUD
ARMREST MEDIAL SWIVEL NOTE:		
<i>Every Corpus 3G Armrest can flip backwards for transfers. The armrests also include a feature that allows the armrest pads to swivel medially towards the center of the seating system. This swivel feature can be locked out by simply moving a bolt located under the armrest bracket. Please indicate below whether each armrest should have the swivel feature enabled or locked.</i>		
<input type="checkbox"/> I10645	Left Armrest Bracket - Medial Swivel Locked	No Charge
<input type="checkbox"/> I10646	Right Armrest Bracket - Medial Swivel Locked	No Charge
<input type="checkbox"/> I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge
<input type="checkbox"/> I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105100-99-0	4"x10" Corpus 3G Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105104-99-0	4"x13" Corpus 3G Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105105-99-0	4"x16" Corpus 3G Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105099-99-0	4"x10" Corpus 3G Armrest Pad Right - Black Mesh	No Charge
<input type="checkbox"/> I105103-99-0	4"x13" Corpus 3G Armrest Pad Right - Black Mesh	No Charge
<input type="checkbox"/> I105108-99-0	4"x16" Corpus 3G Armrest Pad Right - Black Mesh	No Charge
<input type="checkbox"/> I105199-99-0	4"x10" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge
<input type="checkbox"/> I105201-99-0	4"x13" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge
<input type="checkbox"/> I105203-99-0	4"x16" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge

CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105098-99-0	4"x10" Corpus 3G Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105102-99-0	4"x13" Corpus 3G Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105107-99-0	4"x16" Corpus 3G Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105097-99-0	4"x10" Corpus 3G Armrest Pad Left - Black Mesh	No Charge
<input type="checkbox"/> I105101-99-0	4"x13" Corpus 3G Armrest Pad Left - Black Mesh	No Charge
<input type="checkbox"/> I105106-99-0	4"x16" Corpus 3G Armrest Pad Left - Black Mesh	No Charge
<input type="checkbox"/> I105200-99-0	4"x10" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge
<input type="checkbox"/> I105202-99-0	4"x13" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge
<input type="checkbox"/> I105204-99-0	4"x16" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105463-99-0	Corpus 3G Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts.</i>	No Charge
<input type="checkbox"/> I105462-99-0	Corpus 3G Arm Trough Adpater 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts.</i>	
Install Arm Trough Adapters: _____		



Armrest Options

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I100671-99-0	Send Arm Trough Adapters as Parts: _____ Armrest Pouch-Right	
<input type="checkbox"/> I100672-99-0	Armrest Pouch-Left	
<input type="checkbox"/> I10705	Corpus 3G Armrest Side Support - Right Outer <i>Custom Option. Leatherette Pad (2"T x 8.5"L).</i>	
<input type="checkbox"/> I10706	Corpus 3G Armrest Side Support - Left Outer <i>Custom Option. Leatherette Pad (2"T x 8.5"L).</i>	
<input type="checkbox"/> I10707	Corpus 3G Armrest Side Support - Right Inner <i>Custom Option. Leatherette Pad (2"T x 8.5"L). Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	
<input type="checkbox"/> I10708	Corpus 3G Armrest Side Support - Left Inner <i>Custom Option. Leatherette Pad (2"T x 8.5"L). Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	

Legrest Options

CHOOSE LEGREST TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I105080-99-0	Corpus 3G Manual Adj Legrest Elevation	No Charge
<input type="checkbox"/> I105079-99-0	Power Elevating Center Mount Legrests - 170° <i>Includes mechanical articulation.</i>	

FOOTPLATE HEIGHT POSITION:

Part Number	Description	Price AUD
	SET FOOTPLATE HEIGHT AT: _____ <i>The height of the footplates are infinitely adjustable from 6" - 17". This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INLCUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments (ex: 12", 15", etc.), please indicate where you would like the footplates set. Calf Pads not available at footplate heights less than 9".</i>	

CHOOSE FOOTPLATE STYLE:

Part Number	Description	Price AUD
<input type="checkbox"/> I105076-99-0	2pc Corpus 3G Footplates - 9"D x 7.5"W	No Charge
<input type="checkbox"/> I105077-99-0	2pc Narrow Corpus 3G Footplates - 9"D x 6"W	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105252-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware <i>8"W x 8"T Pad.</i>	

Accessories

MUST CHOOSE POSITIONING BELT:

Part Number	Description	Price AUD
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge



Accessories

MUST CHOOSE POSITIONING BELT:

Part Number	Description	Price AUD
<input type="checkbox"/> I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	
<input type="checkbox"/> I10402	Bodypoint Non-Padded Hip Belt, UniTrack - Large	
<input type="checkbox"/> I10408	Bodypoint Padded Hip Belt UniTrack - Large	
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105291-99-0	Headrest Corpus 3G Back, 10"W x 5"T, Leatherette <i>Must also select PN I10552.</i>	
<input type="checkbox"/> I105290-99-0	Headrest Corpus 3G Back, 10"W x 5"T, Black Fabric <i>Must also select PN I10552.</i>	
<input type="checkbox"/> I105212-99-0	Headrest Corpus 3G Back, Ergo Shape, Black Fabric <i>Must also select PN I10552.</i>	
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	
<input type="checkbox"/> I100978-99-0	Abductor Kit - Leatherette <i>Must also select PN I10555.</i>	
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	
<input type="checkbox"/> I10783	Hand Heater Kit with Buddy Button <i>Kit includes Hand Heater (2"H x 1.5"W x 2.75"D), Buddy Button, Power Cable with XLR Plug and PERMOfix Mounting Adapter. Heater has three temp/fan speed settings. For best results, a hand cone (I10784 or I10785) is recommended to trap heat around the user's hand(s). To ensure best positioning for the user, heater will NOT be mounted when chair is received.</i>	
<input type="checkbox"/> I10784	One-Hand Cone Kit with PERMOfix Mounting Set <i>Clear plastic dome (5.5"W x 13.5"L x 8"H) with forearm opening at one end. Kit includes PERMOfix mounting hardware to attach cone to end of armrest. To ensure best positioning for the user, cone will NOT be mounted when wheelchair is received.</i>	
<input type="checkbox"/> I10785	Two-Hand Cone Kit with Mounting Pins <i>Clear plastic triangular dome (15"W x 12.5"L x 8"H) with two forearm openings. Use with Upper Extremity Support is recommended. Kit includes four pins to secure cone to tray surface. Pins remove without tools for easy transfers. To ensure the best positioning for the user, cone will NOT be mounted.</i>	



Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus 3G Seat	
<input type="checkbox"/> I105292-99-0	Push handles, Corpus 3G Seat	
<input type="checkbox"/> I103044-99-0	Cup/Cell Phone Holder	
<input type="checkbox"/> I105563-99-0	Medical Necess Bag Clips - Carabiner Style (C3G)	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	
<input type="checkbox"/> I105144-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i>	
<input type="checkbox"/> I105143-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i>	
<input type="checkbox"/> I105145-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i>	
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105221-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105222-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105220-99-0	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105262-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I105263-99-0	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	
<input type="checkbox"/> I105224-99-0	Oxygen holder	
<input type="checkbox"/> I105223-99-0	Crutch holder	

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10808	Stealth Ultra Headrest <i>Includes Stealth Occipital (500 series) and Sub Occipital (600 series) Pads, bilateral Facial Support Pads (#900) with Stealth Swing Away Hardware (SUS9-B) and Flip Down Headrest Mount (TWB480-FDM). Also includes Universal Headrest Adapter.</i>	
<input type="checkbox"/> I10800	Stealth Swing Away Egg Switch for Ultra Headrest <i>Includes Black Egg Switch and Stealth Ultra Swing Away Hardware (SUS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i>	
<input type="checkbox"/> I10801	Stealth Swing Away Facial Supprt Hdw for Ultra <i>Must also select support pad (I10810, I10811, I10812). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	
<input type="checkbox"/> I10810	Stealth #900 Support Pad (4.5" x 2") - Ultra	
<input type="checkbox"/> I10811	Stealth #910 Support Pad (3" x 2") - Ultra	
<input type="checkbox"/> I10812	Stealth #920 Support Pad (6.4" x 2") - Ultra	



Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10718	Stealth Comfort Plus Headrest - 8"W (CP480) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	
<input type="checkbox"/> I10719	Stealth Comfort Plus Headrest - 10"W (CP280) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	
<input type="checkbox"/> I10720	Stealth Comfort Plus Headrest - 14"W (CP180) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	
<input type="checkbox"/> I10721	Adj Stealth Headrest Hardware - Each (TWB480)	
<input type="checkbox"/> I10736	Stealth Swng Away Egg Sw for Comfort Plus Headrest <i>Includes Black Egg Switch and Stealth Comfort Plus Swing Away Hardware (CPS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i>	
<input type="checkbox"/> I10809	Stealth Swing Away Facial Supprt Hdw for Cmfrt Plus <i>Must also select support pad (I10813, I10814, I10815). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	
<input type="checkbox"/> I10813	Stealth #900 Support Pad (4.5" x 2") - Cmfrt Plus	
<input type="checkbox"/> I10814	Stealth #910 Support Pad (3" x 2") - Cmfrt Plus	
<input type="checkbox"/> I10815	Stealth #920 Support Pad (6.4" x 2") - Cmfrt Plus	
<input type="checkbox"/> I10778	Stealth Lateral Supports 5"D x 3"T (TWBSTL-CS3G35) <i>Must also select I10723. Contains the THIN Profile Stealth Mounting Hardware, which allows for flipping back the armrest for transfers. Works with Corpus 3G backrest shells only.</i>	
<input type="checkbox"/> I10735	Stealth Lateral Supports - 6"D x 4"T (TWBSTL-CS3G) <i>Must also select I10723. Contains the THIN Profile Stealth Mounting Hardware, which allows for flipping back the armrest for transfers. Works with Corpus 3G backrest shells only.</i>	
<input type="checkbox"/> I10723	Adj Removable Stealth Lateral Hardware (Pair)	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price AUD
<input type="checkbox"/> I104690-99-0	Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i>	
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i>	
<input type="checkbox"/> I10754	Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 10cm.</i>	



Alternative Drive Systems

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price AUD
<input type="checkbox"/> I10068	Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i>	
<input type="checkbox"/> I10795	Stealth Swing Away Chin Boom for Ultra Headrest <i>Includes Stealth Ultra Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest:</i> Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____	
<input type="checkbox"/> I10797	Stealth Sw Away Chin Boom for Comfort Plus Headrest <i>Includes Stealth Comfort Plus Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest:</i> Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____	
<input type="checkbox"/> I10807	Stealth Gatlin Eclipse Tray w/ Mini JS Mount <i>Designed to mount on Stealth Gatlin Mid Line Mounts. Includes bracket for mounting a Mini Joystick under tray. Through hole in tray surface provides joystick access. Please indicate side of tray for hole position:</i> Position Hole on Left: _____ Position Hole on Right: _____	

ALT JOYSTICK ACCESSORIES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g.</i> Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	

HEAD ARRAY SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10759	Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 14cm Wide x 14cm Tall. Requires Omni Input Module.</i>	
<input type="checkbox"/> I10639	Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below</i> Mount on Left: _____ Mount on Right: _____	
<input type="checkbox"/> I10760	Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 14cm Wide x 14cm Tall.</i>	
<input type="checkbox"/> I10761	Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i>	

SIP & PUFF SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	Price AUD
--------------------	--------------------	------------------



Alternative Drive Systems

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10763	VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i>	
<input type="checkbox"/> I10764	VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	
<input type="checkbox"/> I10765	VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	
<input type="checkbox"/> I10317	Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 30cm Wide x 24cm Long x 4cm Thick Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	
<input type="checkbox"/> I10638	RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	
<input type="checkbox"/> I10770	Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i>	

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price AUD
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	
<input type="checkbox"/> I10577	Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i>	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge
<input type="checkbox"/> I10079	Mount Seat Control Switchbox on Left	No Charge
<input type="checkbox"/> I10078	Mount Seat Control Switchbox on Right	No Charge

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	
<input type="checkbox"/> I10771	Touch Contact with CR2032 Battery Qty: _____	



Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10772	<p>Dome Switch, Backrest Mnt w Gooseneck Qty: _____</p> <p>Mount on Left: _____ Mount on Right: _____</p> <p>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</p>	

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group.