


K300 PS Jr



*** Indicates a Required Field**

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* I) Knee to Foot: _____

Please send order/quote to fax#: +61 2 8883 2356

Email to: sales.australia@permobil.com

Prices effective January 16, 2012.

Permobil Australia Pty Ltd.
Unit 312/ 14-16 Lexington Dr,
Norwest Business Park
Bella Vista
NSW 2153
Australia
www.permobil.com



Chair Model and Colors

		Price AUD
<input checked="" type="checkbox"/>	I103573-99-0 K300 Base PS - C300 PS0 <i>Price Includes: FWD C300 Power Base, Shock Absorbing Suspension System, R-net Controller, Tie-Down Hardware for Strap Systems, Anti-Tippers, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 8kmh, Weight Capacity = 75kg, STFH = 16.5".</i>	
<input type="checkbox"/>	I10637 10kmh High Speed Upgrade Package - C300/K300 <i>Price Includes: ESP Module, High Speed 10kmh Motor Package, and 120A Power Module. Available with R-net Electronics only.</i>	

CHOOSE DRIVE TIRES:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge
<input type="checkbox"/>	I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge

OPTIONAL DRIVE TIRES:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels.</i>	No Charge

CHOOSE SHROUD COLOUR:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I105761-99-0	C300/K300 Shroud - Cobalt Blue (New Style!)	No Charge
<input type="checkbox"/>	I105762-99-0	C300/K300 Shroud - Purple Rain (New Style!)	No Charge
<input type="checkbox"/>	I105763-99-0	C300/K300 Shroud - Onyx Black (New Style!)	No Charge
<input type="checkbox"/>	I105764-99-0	C300/K300 Shroud - Phantom Gray (New Style!)	No Charge
<input type="checkbox"/>	I105765-99-0	C300/K300 Shroud - Popstar Pink (New Style!)	No Charge
<input type="checkbox"/>	I105766-99-0	C300/K300 Shroud - Volcano Red (New Style!)	No Charge

CHOOSE HUBCAP COLOUR:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I10451	Match Hub Cap to Shroud Color	No Charge
<input type="checkbox"/>	I103537-99-0	Hub Cap Platinum Silver	No Charge

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I101995-99-0	Fixed Seat Tube for Tilt	No Charge
<input type="checkbox"/>	I103779-99-0	Power Adjustable Seat Height - 8" Travel <i>To assist in transfers and improve access for function.</i>	No Charge

MUST CHOOSE BATTERIES:

	Part Number	Description	Price AUD
<input checked="" type="checkbox"/>	IM34 SLDG	Batteries, Grp 34 (60Ah) Sealed LA/Gel, Installed	No Charge

OPTIONAL ITEMS:

	Part Number	Description	Price AUD
<input type="checkbox"/>	I0036-7002A	Enhanced Steering Perform Unit (ESP), R-net <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	No Charge
<input type="checkbox"/>	I105715-99-0	LED Light Kit for C/K300 w/ R-net (and New Shroud)	No Charge



Base Options

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I102768-99-0	Fixed Vent Shelf for C300 <i>If selecting this part, you MUST also pick: PN I0036-7002A - ESP Unit. Fits an LP10 Sized Vent & a U1 Battery. Must be ordered with an expandable controller (ex. R-net Electronics).</i>	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	

Permlock

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104680-99-0	Permlock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permlock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a Licensed Vehicle Modifier.</i>	
<input type="checkbox"/> I10608	Permlock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permlock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> ID51185	R-net Remote Color JS w/ Mono Jacks, No Lights	
<input type="checkbox"/> ID51188	R-net Remote Color JS w/ Toggle Switches, No Lghts	
<input type="checkbox"/> ID51153	R-net Remote Color JS w/ Mono Jacks, Lights	
<input type="checkbox"/> ID51149	R-net Remote Color JS w/ Toggle Switches, Lights	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Spd Pot/Toggle, Lights	
<input type="checkbox"/> I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimble <i>Heavy-Duty modification enables gimble to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick Module	No Charge

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price AUD
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	



Joystick Options

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price AUD
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2	
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2	
<input type="checkbox"/> I103472-99-0	Mount Standard Joystick on Backrest Shell	
<input type="checkbox"/> I10793	Stealth Gatlin Mid Line Mount - 8 Position <i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	
<input type="checkbox"/> I10794	Stealth Gatlin Mid Line Mount - 2 Position <i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	
<input type="checkbox"/> I10803	Stealth Gatlin Bracket for R-net Joystick	
<input type="checkbox"/> I10804	Stealth Gatlin Bracket for Compact Joystick	
<input type="checkbox"/> I10805	Stealth Gatlin Bracket for Mini Joystick	
<input type="checkbox"/> I10806	Stealth Gatlin Bracket for Omni	
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i>	
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I103712-99-0	PS Junior Seat Frame - Only <i>Weight limit 75kg.</i>	No Charge
<input type="checkbox"/> I103715-99-0	PS Jr Seat w/ 45° Manual Tilt <i>Must be ordered with Stroller Handles (I103811-99-0). Weight limit 75kg.</i>	
<input type="checkbox"/> I103716-99-0	PS Junior Seat Frame with 45° Power Tilt <i>Weight limit 75kg.</i>	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103811-99-0	Stroller Handle for PS Jr Seating System <i>Angle adjustable. This item Must be selected when ordering manual tilt (I103715-99-0).</i>	

Seating System Options



Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10585	Seat Function Control Kit for R-net <i>Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price AUD
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox	No Charge
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I103471-99-0	Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i>	No Charge
<input type="checkbox"/> I105026-99-0	Mount ICS Switchbox on Right	No Charge
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	
<input type="checkbox"/> I105027-99-0	Mount ICS Switchbox on Left	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	
<input type="checkbox"/> I104769-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I103717-99-0	PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Solid Back, Sling Back, or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10448	PS Jr Curved Metal Backrest <i>Backrest width will match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the backrest.</i>	No Charge
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge
<input type="checkbox"/> I104616-99-0	Backrest Adjustment Quick Lock, PS / PS Jr	No Charge

Seat Sizing



Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10432	PS Jr Seat Width 11"	No Charge
<input type="checkbox"/> I10433	PS Jr Seat Width 12"	No Charge
<input type="checkbox"/> I10434	PS Jr Seat Width 13"	No Charge
<input type="checkbox"/> I10435	PS Jr Seat Width 14"	No Charge
<input type="checkbox"/> I10436	PS Jr Seat Width 15"	No Charge
<input type="checkbox"/> I10437	PS Jr Seat Width 16"	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10438	PS Jr Seat Depth 10"	No Charge
<input type="checkbox"/> I10439	PS Jr Seat Depth 11"	No Charge
<input type="checkbox"/> I10440	PS Jr Seat Depth 12"	No Charge
<input type="checkbox"/> I10441	PS Jr Seat Depth 13"	No Charge
<input type="checkbox"/> I10442	PS Jr Seat Depth 14"	No Charge
<input type="checkbox"/> I10443	PS Jr Seat Depth 15"	No Charge
<input type="checkbox"/> I10444	PS Jr Seat Depth 16"	No Charge
<input type="checkbox"/> I10445	PS Jr Seat Depth 17"	No Charge
<input type="checkbox"/> I10446	PS Jr Seat Depth 18"	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10276	Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Not available with seat widths or depths greater than 14".</i>	
<input type="checkbox"/> I10038	Ergonomic Seat Cushion Leatherette <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	

Armrest Options

CHOOSE ARMREST STYLE:

Part Number	Description	Price AUD
<input type="checkbox"/> I10824	Flip Back, Height Adjustable Armrests <i>Mounted off the backrest canes. Includes: 45° of angle adjustability, rotational adjustment, 3" offset adjustment both medially and laterally, and height adjustment of 6" - 14" from the seat pan. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): _____</i>	
<input type="checkbox"/> I10825	Flip Back, Height Adjustable Armrests - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (6" - 14"): _____</i>	No Charge

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105100-99-0	4"x10" Corpus 3G Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105104-99-0	4"x13" Corpus 3G Armrest Pad Right - Leatherette	No Charge
<input type="checkbox"/> I105099-99-0	4"x10" Corpus 3G Armrest Pad Right - Black Mesh	No Charge
<input type="checkbox"/> I105103-99-0	4"x13" Corpus 3G Armrest Pad Right - Black Mesh	No Charge
<input type="checkbox"/> I105199-99-0	4"x10" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge
<input type="checkbox"/> I105201-99-0	4"x13" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge
<input type="checkbox"/> I104951-99-0	3.5"x8" Unitrack Armrest Pad Right (C)	No Charge
<input type="checkbox"/> I104950-99-0	3.5"x12.5" Unitrack Armrest Pad Right (G)	No Charge

CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105098-99-0	4"x10" Corpus 3G Armrest Pad Left - Leatherette	No Charge



Armrest Options

CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price AUD
<input type="checkbox"/> I105102-99-0	4"x13" Corpus 3G Armrest Pad Left - Leatherette	No Charge
<input type="checkbox"/> I105097-99-0	4"x10" Corpus 3G Armrest Pad Left - Black Mesh	No Charge
<input type="checkbox"/> I105101-99-0	4"x13" Corpus 3G Armrest Pad Left - Black Mesh	No Charge
<input type="checkbox"/> I105200-99-0	4"x10" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge
<input type="checkbox"/> I105202-99-0	4"x13" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge
<input type="checkbox"/> I104949-99-0	3.5"x12.5" Unitrack Armrest Pad Left (G)	No Charge
<input type="checkbox"/> I104952-99-0	3.5"x8" Unitrack Armrest Pad Left (C)	No Charge

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10586	70° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate.</i>	No Charge
<input type="checkbox"/> I10587	90° Swing Away Legrests w Multi Axis Adj Footplate <i>Includes 5"W x 8"D Multi Axis Footplate.</i>	No Charge
<input type="checkbox"/> I10823	Aeron Legrest Assembly for PS Jr Seat - Narrow <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes Narrow, Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104684-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i>	No Charge
<input type="checkbox"/> I104686-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i>	No Charge
<input type="checkbox"/> IFT202SR	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	No Charge

Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge
<input type="checkbox"/> I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	No Charge
<input type="checkbox"/> I10623	Bodypoint Padded Hip Belt UniTrack - Extra Small	No Charge
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	No Charge
<input type="checkbox"/> I100749-99-0	Miniflex Headrest Removable - Leatherette <i>Must also select PN I10552.</i>	No Charge
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	No Charge
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	No Charge
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	No Charge
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	No Charge
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	No Charge



Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	
<input type="checkbox"/> I101473-99-0	Medical Necessity Bag Clips - Carabiner Style	
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	
<input type="checkbox"/> I103798-99-0	After Market Seating System Adapter <i>Used to connect after market seating systems to the UniTrack Accessory Rail on the PS Jr Seat Frame. Includes Mounting Hardware.</i> Custom Swing Away Laterals Available Upon Request	

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10718	Stealth Comfort Plus Headrest - 8"W (CP480) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	
<input type="checkbox"/> I10721	Adj Stealth Headrest Hardware - Each (TWB480)	
<input type="checkbox"/> I10736	Stealth Swng Away Egg Sw for Comfort Plus Headrest <i>Includes Black Egg Switch and Stealth Comfort Plus Swing Away Hardware (CPS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i>	
<input type="checkbox"/> I10809	Stealth Swing Away Facial Supprt Hdw for Cmfrt Plus <i>Must also select support pad (I10813, I10814, I10815). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	
<input type="checkbox"/> I10813	Stealth #900 Support Pad (4.5" x 2") - Cmfrt Plus	
<input type="checkbox"/> I10814	Stealth #910 Support Pad (3" x 2") - Cmfrt Plus	
<input type="checkbox"/> I10815	Stealth #920 Support Pad (6.4" x 2") - Cmfrt Plus	

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price AUD
<input type="checkbox"/> I104690-99-0	Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	



Alternative Drive Systems

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price AUD
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i>	
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i>	
<input type="checkbox"/> I10754	Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 10cm.</i>	
<input type="checkbox"/> I10068	Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i>	
<input type="checkbox"/> I10797	Stealth Sw Away Chin Boom for Comfort Plus Headrest <i>Includes Stealth Comfort Plus Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest: Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____</i>	
<input type="checkbox"/> I10807	Stealth Gatlin Eclipse Tray w/ Mini JS Mount <i>Designed to mount on Stealth Gatlin Mid Line Mounts. Includes bracket for mounting a Mini Joystick under tray. Through hole in tray surface provides joystick access. Please indicate side of tray for hole position: Position Hole on Left: _____ Position Hole on Right: _____</i>	

ALT JOYSTICK ACCESSORIES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g. Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	

HEAD ARRAY SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10759	Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 14cm Wide x 14cm Tall. Requires Omni Input Module.</i>	
<input type="checkbox"/> I10639	Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I10760	Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 14cm Wide x 14cm Tall.</i>	
<input type="checkbox"/> I10761	Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i>	

SIP & PUFF SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	



Alternative Drive Systems

SIP & PUFF SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10763	VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i>	
<input type="checkbox"/> I10764	VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	
<input type="checkbox"/> I10765	VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	
<input type="checkbox"/> I10317	Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 30cm Wide x 24cm Long x 4cm Thick Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	
<input type="checkbox"/> I10638	RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	
<input type="checkbox"/> I10770	Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i>	

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price AUD
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	
<input type="checkbox"/> I10577	Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i>	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge
<input type="checkbox"/> I10079	Mount Seat Control Switchbox on Left	No Charge
<input type="checkbox"/> I10078	Mount Seat Control Switchbox on Right	No Charge

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
-------------	-------------	-----------



Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	
<input type="checkbox"/> I10771	Touch Contact with CR2032 Battery Qty: _____	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group.