


# K450 MX



**\* Indicates a Required Field**

## Dealer Information

\* Contact: \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

\* City: \_\_\_\_\_

\* State: \_\_\_\_\_

\* Phone#: \_\_\_\_\_

\* Fax#: \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_

\* Last Name: \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* Weight: \_\_\_\_\_

\* Height: \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_

E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\* I) Knee to Foot: \_\_\_\_\_

Please send order/quote to fax#: +61 2 8883 2356

Email to: [sales.australia@permobil.com](mailto:sales.australia@permobil.com)

Prices effective January 16, 2012.

Permobil Australia Pty Ltd.  
Unit 312/ 14-16 Lexington Dr,  
Norwest Business Park  
Bella Vista  
NSW 2153  
Australia  
[www.permobil.com](http://www.permobil.com)



**Chair Model and Colors**

		Price AUD
<input checked="" type="checkbox"/>	I104524-99-0 K450 Base w/ R-net 8kmh - K450 PS2 <i>Price Includes: K450 Power Base, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. Max. Speed = 8 km/h, Weight Capacity = 60 kg, STFH = 75 mm - 635 mm.</i>	

CHOOSE DRIVE TIRES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104525-99-0	Flat Free Drive Tires - K450	No Charge
<input type="checkbox"/> I10607	Pneumatic Drive Tires & Flat Free Casters - K450	No Charge

CHOOSE SHROUD COLOUR:

Part Number	Description	Price AUD
<input type="checkbox"/> I104534-99-0	Shroud for K450 - Chili Pepper Red	No Charge
<input type="checkbox"/> I104535-99-0	Shroud for K450 - Onyx Black	No Charge
<input type="checkbox"/> I104536-99-0	Shroud for K450 - Platinum Silver	No Charge
<input type="checkbox"/> I104537-99-0	Shroud for K450 - Pop Star Pink	No Charge
<input type="checkbox"/> I104538-99-0	Shroud for K450 - Titan Blue	No Charge
<input type="checkbox"/> I104539-99-0	Shroud for K450 - Extreme Yellow	No Charge

**Base Options**

MUST CHOOSE SEAT ELEVATOR:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10510	Power Adjustable Seat Height - 8" Travel <i>This item must be selected.</i>	

MUST CHOOSE BATTERIES:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> IM34 SLDG	Batteries, Grp 34 (60Ah) Sealed LA/Gel, Installed	

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104550-99-0	Anti-Tippers for K450	
<input type="checkbox"/> I104540-99-0	Light Kit - K450 R-net	
<input type="checkbox"/> I104553-99-0	WC 19 Transit Kit - K450 <i>This kit includes additional anchor points for tie down straps, seat stabilizing bars and hardware, and a WC-19 Approved integrated lap and shoulder belt. This kit has been approved according to ANSI/RESNA WC19. Cannot be ordered with a 10"W Seat Size (I104317-99-0). An ICS Switchbox must be installed on the chair in order to access WC19 mode.</i>	
<input type="checkbox"/> I104554-99-0	Fixed Ventilator Tray - K450 <i>Tray is 9.75"D x 13"W, and increases the overall length of a K450 by approximately 11". The ventilator tray is designed to accommodate an LP10 style ventilator and one U1 Battery. A custom charge will need to be applied if ordering vent tray with any of the following items: WC-19 Transit Kit (I104553-99-0), Attendant Control (ID50882), Push Handles (I104576-99-0), or Back Pack Clips (I10598).</i>	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	

**Permolock**

OPTIONAL ITEMS:

Part Number	Description	Price AUD
-------------	-------------	-----------



**Permolock**

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a Licensed Vehicle Modifier.</i>	
<input type="checkbox"/> I104652-99-0	Permolock C3 Dual Locking Pin Kit - K450 <i>This kit includes two pins that are mounted underneath the base of the wheelchair. The pins are fixed and do not push into the wheelchair base.</i>	
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	

**Joystick Options**

**CHOOSE JOYSTICK TYPE:**

Part Number	Description	Price AUD
<input type="checkbox"/> ID51188	R-net Remote Color JS w/ Toggle Switches, No Lghts	
<input type="checkbox"/> ID51153	R-net Remote Color JS w/ Mono Jacks, Lights	
<input type="checkbox"/> ID51149	R-net Remote Color JS w/ Toggle Switches, Lights	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Spd Pot/Toggle, Lights	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode &amp; on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick Module	No Charge

**CHOOSE JOYSTICK HANDLE TYPE:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	

**CHOOSE JOYSTICK MOUNT:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	
<input type="checkbox"/> ID50882-A	Attendant Control for R-net Control System Only <i>On the K450, the mounting location for the Attendant Control is off of the Push Handle Bracket. If ordering an attendant control, you must also order I104576-99-0 (Push Handles - MX Seat).</i>	
<input type="checkbox"/> I104571-99-0	Co-Pilot Attendant Control - K450	



## Joystick Options

### OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> ID51111	<b>R-net BT Mouse Module</b> <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	

## Seating System

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I104324-99-0	<b>MX Power Seating System</b> <i>Price Includes: MX Seating System, Adjustable Seat Frame, Power Seat to Floor System, Adjustable Backrest Angle, and One Piece Footplate. Weight Limit 57kg.</i>	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104529-99-0	<b>No Tilt</b> <i>Weight Limit 57kg.</i>	No Charge
<input type="checkbox"/> I104551-99-0	<b>45° Power Tilt</b> <i>Weight Limit 57kg.</i>	

## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price AUD
<input checked="" type="checkbox"/> I10585	<b>Seat Function Control Kit for R-net</b> <i>Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	

### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price AUD
<input type="checkbox"/> I104136-99-0	Toggles on Seat Function Control Kit	
<input type="checkbox"/> I104135-99-0	Push Buttons on Seat Function Control Kit	No Charge

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price AUD
<input type="checkbox"/> I10626	<b>Mount ICS Switchbox onto R-net Joystick Bracket</b> <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge
<input type="checkbox"/> I103471-99-0	<b>Mount ICS Switchbox to End of Arm Pad - Corpus</b> <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i>	No Charge
<input type="checkbox"/> I105026-99-0	Mount ICS Switchbox on Right	No Charge
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	
<input type="checkbox"/> I105027-99-0	Mount ICS Switchbox on Left	No Charge

### OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge
<input type="checkbox"/> I104196-99-0	<b>ICS Alternative Switchbox for R-net Electronics</b> <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	



## Seating System Options

### OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104769-99-0	<b>Remote Stop System - All Electronics</b> <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	
<input type="checkbox"/> I10515	<b>Memory Seat Program for K450 MX</b> <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions through the ICS Switchbox. On the K450 MX, ALL seat functions are capable of recalling a saved seat position. Please refer to the K450 Owner's Manual for more information on saving and recalling a memory seat position.</i>	No Charge

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<b>MX STYLE BACKRESTS:</b>		
<i>Use this section to select MX Style Backrests. Back heights measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Please be sure to pick both a backrest shell size and a backrest cushion style (Leatherette or Black Mesh Fabric). Note: A 16"W MX Backrest is NOT available.</i>		
<input type="checkbox"/> I104321-99-0	<b>MX Backrest - 10"W x 14"H</b> <i>Backrest can be adjusted to either 14" or 16" high. Backrest will be mounted in the 16" high position if ordering a headrest.</i>	No Charge
<input type="checkbox"/> I104322-99-0	<b>MX Backrest - 12"W x 16"H</b> <i>Backrest can be adjusted to either 16" or 18" high.</i>	No Charge
<input type="checkbox"/> I104323-99-0	<b>MX Backrest - 14"W x 20"H</b> <i>Backrest can be adjusted to either 20" or 22" high.</i>	No Charge

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<b>Please Select Cushion Style Below:</b>		
<input type="checkbox"/> I104309-99-0	MX Back Cushion 270x325, 10"W x 14"T - Leatherette	No Charge
<input type="checkbox"/> I104311-99-0	MX Back Cushion 320x375, 12"W x 16"T - Leatherette	No Charge
<input type="checkbox"/> I104313-99-0	MX Back Cushion 370x475, 14"W x 20"T - Leatherette	No Charge
<input type="checkbox"/> I104308-99-0	MX Back Cushion 270x325, 10"W x 14"T - Black Mesh	No Charge
<input type="checkbox"/> I104310-99-0	MX Back Cushion 320x375, 12"W x 16"T - Black Mesh	No Charge
<input type="checkbox"/> I104312-99-0	MX Back Cushion 370x475, 14"W x 20"T - Black Mesh	No Charge

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<b>PS JR TUBULAR STYLE BACKRESTS</b>		
<i>Use this section to select PS Jr Tubular Style Backrest Canes. Note: A 10"W Tubular Backrest is NOT available.</i>		
<input type="checkbox"/> I103717-99-0	<b>PS Jr Back Canes, Height &amp; Angle Adjustable</b> <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Solid Back, Sling Back, or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i>	No Charge

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I103719-99-0	PS Jr Solid Backrest - 12"W	No Charge
<input type="checkbox"/> I103721-99-0	PS Jr Solid Backrest - 14"W	No Charge
<input type="checkbox"/> I103723-99-0	PS Jr Solid Backrest - 16"W	No Charge

Standard headrests will not mount to sling backs.



**Backrest Options**

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I103725-99-0	PS Jr Sling Backrest - 12"W	No Charge
<input type="checkbox"/> I103727-99-0	PS Jr Sling Backrest - 14"W	No Charge
<input type="checkbox"/> I103729-99-0	PS Jr Sling Backrest - 16"W	No Charge
<input type="checkbox"/> I10449	PS Jr - Omit Backrest	No Charge

**Seat Sizing**

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104317-99-0	MX Seat Size 10"W x 10"D <i>Seat depth can be adjusted to 10", 11", or 12". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. This seat size is not compatible with the WC 19 Kit (I104553-99-0).</i>	No Charge
<input type="checkbox"/> I104318-99-0	MX Seat Size 12"W x 12"D <i>Seat depth can be adjusted to 12", 13", or 14". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan.</i>	No Charge
<input type="checkbox"/> I104319-99-0	MX Seat Size 14"W x 14"D <i>Seat depth can be adjusted to 14", 15", or 16". If ordering this seat size, with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering PS Jr Arms, the lowest allowed seat position will be 10".</i>	No Charge
<input type="checkbox"/> I104320-99-0	MX Seat Size 16"W x 16"D <i>Seat depth can be adjusted to 16", 17", or 18". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 10" from the floor to the top of the seat pan. If ordering PS Jr Arms, the lowest allowed seat position will be 16".</i>	No Charge

OPTIONAL ITEMS:

Part Number	Description
<input type="checkbox"/> I306384-99-90-0	Seat Cushions Below Are Approximately 2" Thick: Leatherette Seat Cush MX - 10"W x 10"D (250x250)
<input type="checkbox"/> I306385-99-90-0	Leatherette Seat Cush MX - 12"W x 12"D (300x300)
<input type="checkbox"/> I306386-99-90-0	Leatherette Seat Cush MX - 14"W x 14"D (350x350)
<input type="checkbox"/> I101818-99-0	Leatherette Seat Cush MX - 16"W x 16"D (400x400)
<input type="checkbox"/> I104305-99-0	Black Mesh Seat Cush MX - 10"W x 10"D (250x250)
<input type="checkbox"/> I104306-99-0	Black Mesh Seat Cush MX - 12"W x 12"D (300x300)
<input type="checkbox"/> I104307-99-0	Black Mesh Seat Cush MX - 14"W x 14"D (350x350)
<input type="checkbox"/> I104702-99-0	Black Mesh Seat Cush MX - 16"W x 16"D (400x400)

**Armrest Options**

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I10512	Height Adjustable Armrest Assembly MX Seat <i>If you are not selecting Tilt, you MUST order this item.</i>	
<input type="checkbox"/> I10514	Height Adjustable Armrest Assembly MX Seat - NC <i>No Charge if ordering Power Tilt (I104551-99-0)</i>	No Charge

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104303-99-0	MX Armrest Pad Short Leatherette 8"L x 3.5"W	No Charge
<input type="checkbox"/> I104304-99-0	MX Armrest Pad Long Leatherette 12.5"Lx3.5"W	No Charge



**Armrest Options**

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104648-99-0	Height Adj Arm Assembly - PS Jr w/ UT Arm Pads <i>Mounted off the Unitrack seat rail. Includes: 40° of angle adjustability, rotational adjustment, 2" offset adjustment both medially and laterally, and height adjustment in 1" of 7" - 11" from the seat from the seat pan. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (7" - 11"): _____</i>	
<input type="checkbox"/> I10519	Height Adj Arm Assembly - PS Jr w/ UT Arm Pads - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (7" - 11"): _____</i>	No Charge

**Legrest Options**

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price AUD
<input type="checkbox"/> I104314-99-0	MX Footplate - Small (11.5" W x 8"D)	No Charge
<input type="checkbox"/> I104315-99-0	MX Footplate - Large (14" W x 9.5"D)	No Charge

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I103427-99-0	Footstrap Kit	
<input type="checkbox"/> I103104-99-0	Footstrap Kit - Padded	
<input type="checkbox"/> I104338-99-0	Small Anti-Skid Overlay for Footplate - K450	
<input type="checkbox"/> I104337-99-0	Large Anti-Skid Overlay for Footplate - K450	
<input type="checkbox"/> IFT202SR	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	

**Accessories**

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge
<input type="checkbox"/> I10404	Bodypoint Padded Hip Belt UniTrack (C300) - Medium	
<input type="checkbox"/> I10623	Bodypoint Padded Hip Belt UniTrack - Extra Small	
<input type="checkbox"/> I104339-99-0	MF Headrest for MX Seat - Leatherette <i>Must also select PN I10552. Also includes a UT Bracket for lateral adjustment of the headrest. Please indicate below which size you need.</i> UT Brack for Lat Head Adj (10"W MX Back): ____ <i>Allows for Lateral Adjustment of the headrest.</i> UT Brack for Lat Head Adj(12-14"W MX Back): ____ <i>Allows for Lateral Adjustment of the headrest.</i>	
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	
<input type="checkbox"/> I104575-99-0	UT Trunk Suppt (B) for MX Seat - 3.5"T x 4.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	
<input type="checkbox"/> I104574-99-0	UT Trunk Suppt (F) for MX Seat - 3.5"T x 8.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i>	
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	



**Accessories**

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I104576-99-0	<b>Push Handles - K450</b> <i>Mounted off the chassis. Height and angle adjustable. Must be selected if ordering bag clips or an attendant control.</i>	
<input type="checkbox"/> I10598	<b>Medical Necessity Bag Clips - Carabiner Style</b> <i>Push Handles (I104576-99-0) must also be selected in order to mount bag clips.</i>	
<input type="checkbox"/> I100598-99-0	<b>Chest Harness</b> <i>Chest harness kit includes a standard positioning belt.</i>	
<input type="checkbox"/> ISH212	<b>BP Shoulder Harness Medium, Front Pull</b> <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	
<input type="checkbox"/> I100578-99-0	<b>Upper Extremity Support Small, Angle Adjustable</b> <i>Must also select PN I10556.</i>	
<input type="checkbox"/> I103426-99-0	<b>Upper Extremity Support Angle Adjust - UniTrack</b> <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	
<input type="checkbox"/> I10556	<b>Adj Removable - Tray Hardware (1 Each)</b>	

**Stealth Accessories**

OPTIONAL ITEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10718	<b>Stealth Comfort Plus Headrest - 8"W (CP480)</b> <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	
<input type="checkbox"/> I10721	<b>Adj Stealth Headrest Hardware - Each (TWB480)</b>	
<input type="checkbox"/> I10736	<b>Stealth Swng Away Egg Sw for Comfort Plus Headrest</b> <i>Includes Black Egg Switch and Stealth Comfort Plus Swing Away Hardware (CPS9-R). Please indicate if you want switch mounted on left or right side of headrest:                      Mount Switch on Left: _____ Mount Switch on Right: _____</i>	
<input type="checkbox"/> I10809	<b>Stealth Swing Away Facial Supprt Hdw for Cmfrt Plus</b> <i>Must also select support pad (I10813, I10814, I10815). Please indicate if you want it mounted on left or right side of headrest:                      Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	
<input type="checkbox"/> I10813	<b>Stealth #900 Support Pad (4.5" x 2") - Cmfrt Plus</b>	
<input type="checkbox"/> I10814	<b>Stealth #910 Support Pad (3" x 2") - Cmfrt Plus</b>	
<input type="checkbox"/> I10815	<b>Stealth #920 Support Pad (6.4" x 2") - Cmfrt Plus</b>	

**Alternative Drive Systems**

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price AUD
<input type="checkbox"/> I10751	<b>Mini-Joystick for Omni / IOM (9-pin connector)</b> <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	
<input type="checkbox"/> I103893-99-0	<b>Compact Joystick for Omni / IOM (9-pin connector)</b> <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	
<input type="checkbox"/> I103894-99-0	<b>Compact Joystick LITE for Omni / IOM (9-pin conn)</b> <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price AUD
<input type="checkbox"/> I104690-99-0	<b>Alt. JS Adjustable Mount Kit - Armrest</b> <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	



**Alternative Drive Systems**

**CHOOSE MOUNTING HARDWARE:**

Part Number	Description	Price AUD
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i>	
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i>	
<input type="checkbox"/> I10754	Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 10cm.</i>	
<input type="checkbox"/> I10068	Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i>	
<input type="checkbox"/> I10797	Stealth Sw Away Chin Boom for Comfort Plus Headrest <i>Includes Stealth Comfort Plus Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest: Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____</i>	

**ALT JOYSTICK ACCESSORIES:**

Part Number	Description	Price AUD
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode &amp; On/Off. Switch operation force = 130 g. Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode &amp; On/Off. Switch operation force = 130 g</i>	

**HEAD ARRAY SYSTEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10759	Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 14cm Wide x 14cm Tall. Requires Omni Input Module.</i>	
<input type="checkbox"/> I10639	Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I10760	Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 14cm Wide x 14cm Tall.</i>	
<input type="checkbox"/> I10761	Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i>	

**SIP & PUFF SYSTEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	



**Alternative Drive Systems**

*ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:*

Part Number	Description	Price AUD
<input type="checkbox"/> I10763	VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i>	
<input type="checkbox"/> I10764	VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	
<input type="checkbox"/> I10765	VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	
<input type="checkbox"/> I10317	Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 30cm Wide x 24cm Long x 4cm Thick Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	
<input type="checkbox"/> I10638	RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	
<input type="checkbox"/> I10770	Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i>	

**Alt Drive Systems Setup**

*CHOOSE OMNI OR INPUT MODULE (REQUIRED):*

Part Number	Description	Price AUD
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	
<input type="checkbox"/> I10577	Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i>	
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	

*CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:*

Part Number	Description	Price AUD
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge
<input type="checkbox"/> I10079	Mount Seat Control Switchbox on Left	No Charge
<input type="checkbox"/> I10078	Mount Seat Control Switchbox on Right	No Charge

**Alt Drive Accessories**

*OPTIONAL ITEMS:*

Part Number	Description	Price AUD
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	
<input type="checkbox"/> I10771	Touch Contact with CR2032 Battery Qty: _____	



**Alt Drive Accessories**

**OPTIONAL ITEMS:**

Part Number	Description	Price AUD
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.	

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group.