

K450 MX



Dealer Information

* **Contact:** _____

Dealer Code: _____

* **Dealer Name:** _____

Address: _____

* **City:** _____

* **Province:** _____

* **Phone#:** _____

* **Fax#:** _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* **First Name:** _____

* **Last Name:** _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* **Weight:** _____

* **Height:** _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* **I) Knee to Foot:** _____



Display ADP Codes on order/quote:

Three column pricing on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Permobil Ltd.
4 - 75 Mary St.
Aurora, Ontario
L4G 1G3

Tel: (888) 726-9330

CAN Fax: (905) 726-9990

US Fax: (800) 231-3256

www.permobil.com

Canadian prices effective **March 23, 2012** and payable in Canadian dollars.

Chair Model and Colors

| Part Number | Description | MSRP - CA | ADP Code |
|--|--|-----------|----------|
| <input checked="" type="checkbox"/> I104524-99-0 | K450 Base w/ R-net 5.2mph - K450 PS2 <i>Price Includes: K450 Power Base, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 5.0mph, Weight Capacity = 125lbs, STFH = 3"-25".</i> | 12,963.00 | |

CHOOSE DRIVE TIRES:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I104525-99-0 | Flat Free Drive Tires - K450 | No Charge | |
| <input type="checkbox"/> I10607 | Pneumatic Drive Tires & Flat Free Casters - K450 | No Charge | |

CHOOSE SHROUD COLOR:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|------------------------------------|-----------|----------|
| <input type="checkbox"/> I104534-99-0 | Shroud for K450 - Chili Pepper Red | No Charge | |
| <input type="checkbox"/> I104535-99-0 | Shroud for K450 - Onyx Black | No Charge | |
| <input type="checkbox"/> I104536-99-0 | Shroud for K450 - Platinum Silver | No Charge | |
| <input type="checkbox"/> I104537-99-0 | Shroud for K450 - Pop Star Pink | No Charge | |
| <input type="checkbox"/> I104538-99-0 | Shroud for K450 - Titan Blue | No Charge | |
| <input type="checkbox"/> I104539-99-0 | Shroud for K450 - Extreme Yellow | No Charge | |

Base Options

MUST CHOOSE SEAT ELEVATOR:

| Part Number | Description | MSRP - CA | ADP Code |
|--|--|-----------|----------|
| <input checked="" type="checkbox"/> I10510 | Power Adjustable Seat Height - 8" Travel <i>This item must be selected.</i> | 2,827.00 | |

MUST CHOOSE BATTERIES:

| Part Number | Description | MSRP - CA | ADP Code |
|---|---|-----------|----------|
| <input checked="" type="checkbox"/> IM34 SLDG | Batteries, Grp 34 (60Ah) Sealed LA/Gel, Installed | 814.00 | |

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--------------------------|-----------|----------|
| <input type="checkbox"/> I104550-99-0 | Anti-Tippers for K450 | 448.00 | |
| <input type="checkbox"/> I104540-99-0 | Light Kit - K450 R-net | 1,117.00 | |
| <input type="checkbox"/> I104553-99-0 | WC 19 Transit Kit - K450 | 1,796.00 | |

This kit includes additional anchor points for tie down straps, seat stabilizing bars and hardware, and a WC-19 Approved integrated lap and shoulder belt. This kit has been approved according to ANSI/RESNA WC19. Cannot be ordered with a 10"W Seat Size (I104317-99-0). An ICS Switchbox must be installed on the chair in order to access WC19 mode.

| | | | |
|---------------------------------------|--|--------|--|
| <input type="checkbox"/> I104554-99-0 | Fixed Ventilator Tray - K450 <i>Tray is 9.75"D x 13"W, and increases the overall length of a K450 by approximately 11". The ventilator tray is designed to accommodate an LP10 style ventilator and one U1 Battery. A custom charge will need to be applied if ordering vent tray with any of the following items: WC-19 Transit Kit (I104553-99-0), Attendant Control (ID50882), Push Handles (I104576-99-0), or Back Pack Clips (I10598).</i> | 541.00 | |
| <input type="checkbox"/> ID51113 | Programming Key for R-net - Permobil Version | 448.00 | |

Permlock

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104680-99-0 | Permlock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permlock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to: NAME: _____ ADDR: _____</i> | 1,828.00 | |

Permolock

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104652-99-0 | Permolock C3 Dual Locking Pin Kit - K450 <i>This kit includes two pins that are mounted underneath the base of the wheelchair. The pins are fixed and do not push into the wheelchair base.</i> | 613.00 | |
| <input type="checkbox"/> I10608 | Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i> | 613.00 | |

Joystick Options

CHOOSE JOYSTICK TYPE:

| Part Number | Description | MSRP - CA | ADP Code |
|--|---|-----------|----------|
| PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i> | | | |
| <input type="checkbox"/> ID51188 | R-net Remote Colour JS w/ Toggle Switches, No Lghts | 1,746.00 | |
| <input type="checkbox"/> ID51153 | R-net Remote Colour JS w/ Mono Jacks, Lights | 1,540.00 | WEX1 |
| <input type="checkbox"/> ID51149 | R-net Remote Colour JS w/ Toggle Switches, Lights | 1,700.00 | WEX1 |
| <input type="checkbox"/> ID51190 | R-net Remote Colour JS w/ Spd Pot/Toggle, Lights | 1,746.00 | WEX1 |
| <input type="checkbox"/> I103814-99-0 | Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i> | 1,700.00 | |
| <input type="checkbox"/> I10621 | Add Wrist Pad to Compact Joystick Advanced | 216.00 | |
| <input type="checkbox"/> I105246-99-0 | R-net LED Remote Joystick Module | No Charge | |

CHOOSE JOYSTICK HANDLE TYPE:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---------------------------------------|-----------|----------|
| <input type="checkbox"/> I10045 | Cone Shaped Joystick Knob - Standard. | No Charge | |
| <input type="checkbox"/> I100960-99-0 | Chin Cup For Joystick Knob | 102.00 | |
| <input type="checkbox"/> I100961-99-0 | "T" Handle For Joystick Knob | 150.00 | |
| <input type="checkbox"/> I100962-99-0 | "Large Ball" For Joystick Knob | 102.00 | |
| <input type="checkbox"/> I100963-99-0 | "Softball" For Joystick Knob | 137.00 | |
| <input type="checkbox"/> I100964-99-0 | "Mushroom" For Joystick Knob | 137.00 | |
| <input type="checkbox"/> I100965-99-0 | "Stick" For Joystick Knob | 137.00 | |

CHOOSE JOYSTICK MOUNT:

| | | | |
|---------------------------------------|--|--------|------|
| <input type="checkbox"/> I102537-99-0 | Joystick Mount Right - Swing Away, R-net & VR2 | 335.00 | WEPM |
| <input type="checkbox"/> I102538-99-0 | Joystick Mount Left - Swing Away, R-net & VR2 | 335.00 | WEPM |

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10411 | Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i> | No Charge | |
| <input type="checkbox"/> I10513 | Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i> | No Charge | |
| <input type="checkbox"/> ID50882-A | Attendant Control for R-net Control System Only <i>On the K450, the mounting location for the Attendant Control is off of the Push Handle Bracket. If ordering an attendant control, you must also order I104576-99-0 (Push Handles - MX Seat).</i> | 1,117.00 | |
| <input type="checkbox"/> I104571-99-0 | Co-Pilot Attendant Control - K450 | 3,031.00 | |

Joystick Options

| Part Number | Description | MSRP - CA | ADP Code |
|----------------------------------|--|-----------|----------|
| <input type="checkbox"/> ID51111 | R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i> | 901.00 | |

Seating System

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|--|--|-----------|----------|
| <input checked="" type="checkbox"/> I104324-99-0 | MX Power Seating System <i>Price Includes: MX Seating System, Adjustable Seat Frame, Power Seat to Floor System, Adjustable Backrest Angle, and One Piece Footplate. Weight Limit 125lbs.</i> | 7,314.00 | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104529-99-0 | No Tilt <i>Weight Limit 125lbs.</i> | No Charge | |
| <input type="checkbox"/> I104551-99-0 | 45° Power Tilt <i>Weight Limit 125lbs.</i> | 6,015.00 | |

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

| Part Number | Description | MSRP - CA | ADP Code |
|--|--|-----------|-----------|
| <input checked="" type="checkbox"/> I10585 | Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i> | 2,364.00 | SEC1P1005 |

CHOOSE PUSH BUTTONS OR TOGGLES:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104136-99-0 | Toggles on Seat Function Control Kit | 232.00 | |
| <input type="checkbox"/> I104135-99-0 | Push Buttons on Seat Function Control Kit | No Charge | |

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10626 | Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i> | No Charge | |
| <input type="checkbox"/> I103471-99-0 | Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i> | No Charge | |
| <input type="checkbox"/> I105026-99-0 | Mount ICS Switchbox on Right | No Charge | |
| <input type="checkbox"/> I103470-99-0 | Mount ICS Switchbox on Backrest | 129.00 | |
| <input type="checkbox"/> I105027-99-0 | Mount ICS Switchbox on Left | No Charge | |

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I105409-99-0 | Do Not Install ICS Switchbox, Send as Part | No Charge | |
| <input type="checkbox"/> I104196-99-0 | ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i> | 1,117.00 | |

Seating System Options

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I104769-99-0 | Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i> | 1,065.00 | |
| <input type="checkbox"/> I10515 | Memory Seat Program for K450 MX <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions through the ICS Switchbox. On the K450 MX, ALL seat functions are capable of recalling a saved seat position. Please refer to the K450 Owner's Manual for more information on saving and recalling a memory seat position.</i> | No Charge | |

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I10610 | MX STYLE BACKRESTS: <i>Use this section to select MX Style Backrests. Back heights measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Please be sure to pick both a backrest shell size and a backrest cushion style (Leatherette or Black Mesh Fabric). Note: A 16"W MX Backrest is NOT available.</i> | No Charge | |
| <input type="checkbox"/> I104321-99-0 | MX Backrest - 10"W x 14"H <i>Backrest can be adjusted to either 14" or 16" high. Backrest will be mounted in the 16" high position if ordering a headrest.</i> | No Charge | |
| <input type="checkbox"/> I104322-99-0 | MX Backrest - 12"W x 16"H <i>Backrest can be adjusted to either 16" or 18" high.</i> | No Charge | |
| <input type="checkbox"/> I104323-99-0 | MX Backrest - 14"W x 20"H <i>Backrest can be adjusted to either 20" or 22" high.</i> | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10612 | Please Select Cushion Style Below: | No Charge | |
| <input type="checkbox"/> I104309-99-0 | MX Back Cushion 270x325, 10"W x 14"T - Leatherette | No Charge | |
| <input type="checkbox"/> I104311-99-0 | MX Back Cushion 320x375, 12"W x 16"T - Leatherette | No Charge | |
| <input type="checkbox"/> I104313-99-0 | MX Back Cushion 370x475, 14"W x 20"T - Leatherette | No Charge | |
| <input type="checkbox"/> I104308-99-0 | MX Back Cushion 270x325, 10"W x 14"T - Black Mesh | No Charge | |
| <input type="checkbox"/> I104310-99-0 | MX Back Cushion 320x375, 12"W x 16"T - Black Mesh | No Charge | |
| <input type="checkbox"/> I104312-99-0 | MX Back Cushion 370x475, 14"W x 20"T - Black Mesh | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10611 | PS JR TUBULAR STYLE BACKRESTS <i>Use this section to select PS Jr Tubular Style Backrest Canes. Note: A 10"W Tubular Backrest is NOT available.</i> | No Charge | |
| <input type="checkbox"/> I103717-99-0 | PS Jr Back Canes, Height & Angle Adjustable <i>Must also select PS Jr Style Armrests if ordering this backrest. Choose PS Jr Style Backrest Type (Solid Back, Sling Back, or Omit Back) below. Backrest width should match seat width. Back height is adjustable from 16.5"T - 22.5"T, measured from the seat pan to the top of the back.</i> | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I103719-99-0 | PS Jr Solid Backrest - 12"W | No Charge | |
| <input type="checkbox"/> I103721-99-0 | PS Jr Solid Backrest - 14"W | No Charge | |
| <input type="checkbox"/> I103723-99-0 | PS Jr Solid Backrest - 16"W | No Charge | |
| <input type="checkbox"/> I10615 | Standard headrests will not mount to sling backs. | No Charge | |
| <input type="checkbox"/> I103725-99-0 | PS Jr Sling Backrest - 12"W | No Charge | |

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|-----------------------------|-----------|----------|
| <input type="checkbox"/> I103727-99-0 | PS Jr Sling Backrest - 14"W | No Charge | |
| <input type="checkbox"/> I103729-99-0 | PS Jr Sling Backrest - 16"W | No Charge | |
| <input type="checkbox"/> I10449 | PS Jr - Omit Backrest | No Charge | |

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104325-99-0 | K450 Reinforcement for MX Seat w/ PS Jr Back (12"W) | No Charge | |
| <input type="checkbox"/> I104326-99-0 | K450 Reinforcement for MX Seat w/ PS Jr Back (14"W) | No Charge | |
| <input type="checkbox"/> I104327-99-0 | K450 Reinforcement for MX Seat w/ PS Jr Back (16"W) | No Charge | |

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104317-99-0 | MX Seat Size 10"W x 10"D <i>Seat depth can be adjusted to 10", 11", or 12". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. This seat size is not compatible with the WC 19 Kit (I104553-99-0).</i> | No Charge | |
| <input type="checkbox"/> I104318-99-0 | MX Seat Size 12"W x 12"D <i>Seat depth can be adjusted to 12", 13", or 14". If ordering this seat size, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan.</i> | No Charge | |
| <input type="checkbox"/> I104319-99-0 | MX Seat Size 14"W x 14"D <i>Seat depth can be adjusted to 14", 15", or 16". If ordering this seat size, with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 3.5" from the floor to the top of the seat pan. If ordering PS Jr Arms, the lowest allowed seat position will be 10".</i> | No Charge | |
| <input type="checkbox"/> I104320-99-0 | MX Seat Size 16"W x 16"D <i>Seat depth can be adjusted to 16", 17", or 18". If ordering this seat size with MX Arms, the lowest allowed seat position when using the seat to floor function will be approximately 10" from the floor to the top of the seat pan. If ordering PS Jr Arms, the lowest allowed seat position will be 16".</i> | No Charge | |

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|--|--|-----------|----------|
| <input type="checkbox"/> I10613 | Seat Cushions Below Are Approximately 2" Thick: | No Charge | |
| <input type="checkbox"/> I306384-99-90-0 | Leatherette Seat Cush MX - 10"W x 10"D (250x250) | 539.00 | |
| <input type="checkbox"/> I306385-99-90-0 | Leatherette Seat Cush MX - 12"W x 12"D (300x300) | 539.00 | |
| <input type="checkbox"/> I306386-99-90-0 | Leatherette Seat Cush MX - 14"W x 14"D (350x350) | 539.00 | |
| <input type="checkbox"/> I101818-99-0 | Leatherette Seat Cush MX - 16"W x 16"D (400x400) | 539.00 | |
| <input type="checkbox"/> I104305-99-0 | Black Mesh Seat Cush MX - 10"W x 10"D (250x250) | 500.00 | |
| <input type="checkbox"/> I104306-99-0 | Black Mesh Seat Cush MX - 12"W x 12"D (300x300) | 500.00 | |
| <input type="checkbox"/> I104307-99-0 | Black Mesh Seat Cush MX - 14"W x 14"D (350x350) | 500.00 | |
| <input type="checkbox"/> I104702-99-0 | Black Mesh Seat Cush MX - 16"W x 16"D (400x400) | 500.00 | |

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|---|-----------|----------|
| <input type="checkbox"/> I10512 | Height Adjustable Armrest Assembly MX Seat <i>If you are not selecting Tilt, you MUST order this item.</i> | 321.00 | |
| <input type="checkbox"/> I10514 | Height Adjustable Armrest Assembly MX Seat - NC <i>No Charge if ordering Power Tilt (I104551-99-0)</i> | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|-------------|-------------|-----------|----------|
|-------------|-------------|-----------|----------|

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I104303-99-0 | MX Armrest Pad Short Leatherette 8"L x 3.5"W | No Charge | |
| <input type="checkbox"/> I104304-99-0 | MX Armrest Pad Long Leatherette 12.5"Lx3.5"W | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104648-99-0 | Height Adj Arm Assembly - PS Jr w/ UT Arm Pads <i>Mounted off the Unitrack seat rail. Includes: 40° of angle adjustability, rotational adjustment, 2" offset adjustment both medially and laterally, and height adjustment in 1" of 7" - 11" from the seat pan. In one inch increments, please indicate you would like the armrests set. SET ARMRESTS AT (7" - 11"): _____</i> | 293.00 | |
| <input type="checkbox"/> I10519 | Height Adj Arm Assembly - PS Jr w/ UT Arm Pads - NC <i>No charge when ordering power tilt. In one inch increments, please indicate where you would like the armrests set. SET ARMRESTS AT (7" - 11"): _____</i> | No Charge | |

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--------------------------------------|-----------|----------|
| <input type="checkbox"/> I104314-99-0 | MX Footplate - Small (11.5" W x 8"D) | No Charge | |
| <input type="checkbox"/> I104315-99-0 | MX Footplate - Large (14" W x 9.5"D) | No Charge | |

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I103427-99-0 | Footstrap Kit | 187.00 | |
| <input type="checkbox"/> I103104-99-0 | Footstrap Kit - Padded | 210.00 | |
| <input type="checkbox"/> I104338-99-0 | Small Anti-Skid Overlay for Footplate - K450 | 142.00 | |
| <input type="checkbox"/> I104337-99-0 | Large Anti-Skid Overlay for Footplate - K450 | 142.00 | |

Accessories

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|-----------|
| <input type="checkbox"/> I10400 | Bodypoint Non-Padded Hip Belt, UniTrack - Medium | No Charge | |
| <input type="checkbox"/> I10404 | Bodypoint Padded Hip Belt UniTrack (C300) - Medium | 123.00 | |
| <input type="checkbox"/> I10623 | Bodypoint Padded Hip Belt UniTrack - Extra Small | 123.00 | |
| <input type="checkbox"/> I104339-99-0 | MF Headrest for MX Seat - Leatherette <i>Must also select PN I10552. Also includes a UT Bracket for lateral adjustment of the headrest. Please indicate below which size you need.</i> UT Bracket for Lat Head Adj (10"W MX Back): _____ <i>Allows for Lateral Adjustment of the headrest.</i> UT Brack for Lat Head Adj (12-14"W MX Back): _____ <i>Allows for Lateral Adjustment of the headrest.</i> | 422.00 | |
| <input type="checkbox"/> IHEADADAPTER | Universal Headrest Adapter <i>For aftermarket headrests.</i> | 213.00 | |
| <input type="checkbox"/> I10552 | Adj Removable - Headrest Hardware (1 Each) | 113.00 | SEHND0020 |
| <input type="checkbox"/> I103077-99-0 | UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i> | 422.00 | SELND0015 |
| <input type="checkbox"/> I103079-99-0 | UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i> | 422.00 | SELND0015 |
| <input type="checkbox"/> I10550 | Adj Removable - Thigh Support Hardware (Pair) | 113.00 | |
| <input type="checkbox"/> I104575-99-0 | UT Trunk Suppt (B) for MX Seat - 3.5"T x 4.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i> | 604.00 | |
| <input type="checkbox"/> I104574-99-0 | UT Trunk Suppt (F) for MX Seat - 3.5"T x 8.5"D <i>Must also select the Swing Away Trunk Support Hardware (I10551).</i> | 604.00 | |
| <input type="checkbox"/> I10551 | Adj Removable - Trunk Support Hardware (Pair) | 113.00 | |

Accessories

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|-----------|
| <input type="checkbox"/> I104576-99-0 | Push Handles - K450 <i>Mounted off the chassis. Height and angle adjustable. Must be selected if ordering bag clips or an attendant control.</i> | 391.00 | |
| <input type="checkbox"/> I10598 | Medical Necessity Bag Clips - Carabiner Style <i>Push Handles (I104576-99-0) must also be selected in order to mount bag clips.</i> | 77.00 | |
| <input type="checkbox"/> I100598-99-0 | Chest Harness <i>Chest harness kit includes a standard positioning belt.</i> | 382.00 | |
| <input type="checkbox"/> I100578-99-0 | Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i> | 389.00 | SETND0010 |
| <input type="checkbox"/> I103426-99-0 | Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i> | 389.00 | SETND0010 |
| <input type="checkbox"/> I10556 | Adj Removable - Tray Hardware (1 Each) | 113.00 | |

Stealth Accessories

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10718 | Stealth Comfort Plus Headrest - 8"W (CP480) <i>Must also select I10721. Includes Universal Headrest Adapter.</i> | 358.00 | |
| <input type="checkbox"/> I10721 | Adj Stealth Headrest Hardware - Each (TWB480) | 242.00 | |
| <input type="checkbox"/> I10736 | Stealth Swng Away Egg Sw for Comfort Plus Headrest <i>Includes Black Egg Switch and Stealth Comfort Plus Swing Away Hardware (CPS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i> | 551.00 | |
| <input type="checkbox"/> I10809 | Stealth Swing Away Facial Supprt Hdw for Cmfrt Plus <i>Must also select support pad (I10813, I10814, I10815). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i> | 325.00 | |
| <input type="checkbox"/> I10813 | Stealth #900 Support Pad (4.5" x 2") - Cmfrt Plus | 127.00 | |
| <input type="checkbox"/> I10814 | Stealth #910 Support Pad (3" x 2") - Cmfrt Plus | 127.00 | |
| <input type="checkbox"/> I10815 | Stealth #920 Support Pad (6.4" x 2") - Cmfrt Plus | 127.00 | |

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I10751 | Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i> | 2,590.00 | |
| <input type="checkbox"/> I103893-99-0 | Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i> | 1,264.00 | |
| <input type="checkbox"/> I103894-99-0 | Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i> | 1,446.00 | |

CHOOSE MOUNTING HARDWARE:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> I104690-99-0 | Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> | 234.00 | |
| <input type="checkbox"/> I104695-99-0 | Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i> | 746.00 | |
| <input type="checkbox"/> I104143-99-0 | Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i> | 796.00 | |

Alternative Drive Systems

CHOOSE MOUNTING HARDWARE:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|---|-----------|----------|
| <input type="checkbox"/> I10752 | Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i> | 568.00 | |
| <input type="checkbox"/> I10753 | Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i> | 568.00 | |
| <input type="checkbox"/> I10754 | Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 10cm.</i> | 568.00 | |
| <input type="checkbox"/> I10068 | Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i> | 4,323.00 | |
| <input type="checkbox"/> I10797 | Stealth Sw Away Chin Boom for Comfort Plus Headrest <i>Includes Stealth Comfort Plus Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest: Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____</i> | 523.00 | |

ALT JOYSTICK ACCESSORIES:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|---|-----------|----------|
| <input type="checkbox"/> I104904-99-0 | Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i> | 236.00 | |
| <input type="checkbox"/> I104641-99-0 | Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i> | 428.00 | |

HEAD ARRAY SYSTEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10759 | Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 14cm Wide x 14cm Tall. Requires Omni Input Module.</i> | 4,300.00 | |
| <input type="checkbox"/> I10639 | Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below Mount on Left: _____ Mount on Right: _____</i> | 500.00 | |
| <input type="checkbox"/> I10760 | Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 14cm Wide x 14cm Tall.</i> | 4,506.00 | |
| <input type="checkbox"/> I10761 | Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i> | 236.00 | |

SIP & PUFF SYSTEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|---|-----------|----------|
| <input type="checkbox"/> I10573 | Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i> | 1,519.00 | WEX5 |
| <input type="checkbox"/> I10574 | Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i> | 2,022.00 | WEX5 |
| <input type="checkbox"/> I10762 | Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i> | 236.00 | |

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10763 | VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i> | 2,697.00 | |
| <input type="checkbox"/> I10764 | VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> | 236.00 | |

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10765 | VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i> | 746.00 | |
| <input type="checkbox"/> I10317 | Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 30cm Wide x 24cm Long x 4cm Thick Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i> | 4,731.00 | |
| <input type="checkbox"/> I10767 | Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i> | 2,667.00 | |
| <input type="checkbox"/> I10766 | Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i> | 236.00 | |
| <input type="checkbox"/> I10638 | RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i> | 5,021.00 | |
| <input type="checkbox"/> I10770 | Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i> | 236.00 | |

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

| Part Number | Description | MSRP - CA | ADP Code |
|------------------------------------|---|-----------|----------|
| <input type="checkbox"/> ID51154IR | Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i> | 2,117.00 | |
| <input type="checkbox"/> I10577 | Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i> | 113.00 | |
| <input type="checkbox"/> I10581 | Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i> | 1,494.00 | |

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------|--|-----------|----------|
| <input type="checkbox"/> I10080 | Mount Seat Control Switchbox on Backrest | No Charge | |
| <input type="checkbox"/> I10079 | Mount Seat Control Switchbox on Left | No Charge | |
| <input type="checkbox"/> I10078 | Mount Seat Control Switchbox on Right | No Charge | |

Alt Drive Accessories

OPTIONAL ITEMS:

| Part Number | Description | MSRP - CA | ADP Code |
|---------------------------------------|--|-----------|----------|
| <input type="checkbox"/> IMICROSWITCH | Microlite Switch Qty: _____ | 152.00 | |
| <input type="checkbox"/> IBUDDYBUTTON | Buddy Button Qty: _____ | 114.00 | |
| <input type="checkbox"/> I5808 | Mini Cup Switch Qty: _____ | 111.00 | |
| <input type="checkbox"/> IZTS-P | Proximity Switch Qty: _____ | 530.00 | |
| <input type="checkbox"/> IMINIJOY5SW | Mini Joystick 5 Switch OC | 901.00 | |
| <input type="checkbox"/> I10771 | Touch Contact with CR2032 Battery Qty: _____ | 255.00 | |
| <input type="checkbox"/> I10772 | Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i> | 472.00 | |

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The ADP codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*