

Chairman Koala Mini-Flex



Dealer Information

* **Contact:** _____

Dealer Code: _____

* **Dealer Name:** _____

Address: _____

* **City:** _____

* **Province:** _____

* **Phone#:** _____

* **Fax#:** _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* **First Name:** _____

* **Last Name:** _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* **Weight:** _____

* **Height:** _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* **I) Knee to Foot:** _____



Display ADP Codes on order/quote:

Three column pricing on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Permobil Ltd.

4 - 75 Mary St.
Aurora, Ontario
L4G 1G3

Tel: (888) 726-9330

CAN Fax: (905) 726-9990

US Fax: (800) 231-3256

www.permobil.com

Canadian prices effective **January 16, 2012** and payable in Canadian dollars.

Chair Model and Colors

Part Number	Description	MSRP - CA	ADP Code
<input checked="" type="checkbox"/> I100855-99-0	Koala Miniflex w/Pilot+ <i>Price Includes: FWD Koala Power Base, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	8,235.00	

CHOOSE DRIVE TIRES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I103889-99-0	Tri Spoke Split Rim w/ FF Drive Tires - 2.5 x 8"	No Charge	
<input type="checkbox"/> I103890-99-0	Tri Spoke Split Rim w/ Pneumatic Tire - 2.5 x 8"	No Charge	

CHOOSE SHROUD COLOR:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I100859-99-0	Shroud for Koala - Yellow	No Charge	
<input type="checkbox"/> I100858-99-0	Shroud for Koala - Red	No Charge	

CHOOSE HUBCAP COLOR:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver	No Charge	

Base Options

MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA	ADP Code
<input checked="" type="checkbox"/> IM40-12 SLD G	Batteries, Koala (40Ah) Sealed LA/Gel, Installed	739.00	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I100864-99-0	Light Kit - Koala	824.00	
<input type="checkbox"/> I602307-99-0	P+ PP1A Programmer - VSI, VR2 or Pilot+ Controllers	613.00	
<input type="checkbox"/> I100865-99-0	Front Bumper Koala	510.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> ID49670/3	Pilot+ series Joystick w/ 7 Buttons - No Charge	No Charge	
<input type="checkbox"/> I10207	Pilot+ Series Joystick w/ 7 Buttons	1,540.00	WEX1
<input type="checkbox"/> ID49880	Pilot+ Joystick M/J	1,643.00	WEX1
<input type="checkbox"/> ID50002	Pilot+Series Joystick T/S	1,643.00	WEX1

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	102.00	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	150.00	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	102.00	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	137.00	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	137.00	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	137.00	

CHOOSE JOYSTICK MOUNT:

<input type="checkbox"/> I100849-99-0	Joystick Mount-Right Swing Away	335.00	
<input type="checkbox"/> I100848-99-0	Joystick Mount-Left Swing Away	335.00	

Joystick Options

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net <i>Please indicate which armrest to mount midline below.</i> Mount Midline Joystick on Left Armrest: _____ Mount Midline Joystick on Right Armrest: _____	335.00	WE3D
<input type="checkbox"/> I10450	Expandable Controller - Pilot Plus <i>MUST be ordered with the expandable/remote joysticks (PN: I10207, ID49880, or ID50002).</i>	No Charge	
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	No Charge	
<input type="checkbox"/> ID49682	P+DAC Attendant Control - Pilot+ Controller Only <i>Not available with VR2 controllers. Will be sent as part if ordering sling back upholstery.</i>	1,170.00	
<input type="checkbox"/> I103182-99-0	Co-Pilot Attendant Control - MiniFlex <i>Custom charge will be applied if ordering with a tubular backrest. (PN: I10113, I10114).</i>	2,705.00	

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> IKOALA-MFII	Mini Flex Seating System <i>Price Includes: Mini Flex Seating System, Power Tilt System, Adjustable Seat Frame, Adjustable Backrest and Legrest Angle, Positioning Belt and One Piece Footrest. Weight Limit 110 lbs.</i>	6,692.00	

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10128	Single Seat Function Control Switch	No Charge	
<input type="checkbox"/> I10129	Single Seat Function Control Kit <i>For operation of a single seat function through the wheelchair electronics via the drive controls and/or membrane switch box. Order with Pilot+ or remote joystick options. Not available with VSI Controllers.</i>	1,599.00	
<input checked="" type="checkbox"/> I100579-99-0	Multiple Seat Function Control Kit <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive controls and/or membrane switch box. Order with Pilot+ or remote joystick options. Not available with VSI Controllers.</i>	2,364.00	SEC1P1005

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I100524-99-0	Toggle Switch Kit For Seat Function Control Kit <i>Add on part for either the multiple seat function control kit or single seat function control kit in order.</i>	138.00	
<input type="checkbox"/> ISW-JOY	Seat Control Joyswitch <i>Multi directional toggle to operate up to 4 power seat functions. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if selecting this part. Please specify where to mount membrane switchbox: Right Arm Pad: _____ Left Armrest Pad: _____</i>	1,597.00	
<input type="checkbox"/> I10253	Mount Joyswitch on Left	No Charge	
<input type="checkbox"/> I10254	Mount Joyswitch on Right	No Charge	
<input type="checkbox"/> ISW-ALT	Altswitch Alternating Seat Control <i>Seat function interface box which accommodates 4 single momentary switches. 1/8" mini-mom switches not included. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if selecting this part.</i>	1,179.00	

Seating System Options

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	1,015.00	

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Ergo Backs include lateral wedges and lumbar pads.</i>			
<input type="checkbox"/> I100764-99-0	Backrest Kit 10"W x 14"H Leatherette	No Charge	
<input type="checkbox"/> I100765-99-0	Backrest Kit 12"W x 16"H Leatherette	No Charge	
<input type="checkbox"/> I100766-99-0	Backrest Kit 14"W x 20"H Leatherette	No Charge	
<input type="checkbox"/> I100761-99-0	Backrest Kit 10"W x 14"H Fabric	No Charge	
<input type="checkbox"/> I100762-99-0	Backrest Kit 12"W x 16"H Fabric	No Charge	
<input type="checkbox"/> I100763-99-0	Backrest Kit 14"W x 20"H Fabric	No Charge	
<input type="checkbox"/> I10113	Tubular Backrest Frame Mini-Flex, 12" Wide <i>Backrest canes to mount aftermarket backrest systems.</i>	No Charge	
<input type="checkbox"/> I10114	Tubular Backrest Frame Mini-Flex, 14" Wide <i>Backrest canes to mount aftermarket backrest systems.</i>	No Charge	

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I100757-99-0	Seat Sizing 10"W x 10"D Leatherette	No Charge	
<input type="checkbox"/> I100758-99-0	Seat Sizing 12"W x 12"D Leatherette	No Charge	
<input type="checkbox"/> I100759-99-0	Seat Sizing 14"W x 14"D Leatherette	No Charge	
<input type="checkbox"/> I100754-99-0	Seat Sizing 10"W x 10"D Fabric	No Charge	
<input type="checkbox"/> I100755-99-0	Seat Sizing 12"W x 12"D Fabric	No Charge	
<input type="checkbox"/> I100756-99-0	Seat Sizing 14"W x 14"D Fabric	No Charge	
<input type="checkbox"/> I10115	Tubular Seat Frame Mini-Flex, 12" Wide <i>Seat frame tubes to mount aftermarket seating systems.</i>	No Charge	
<input type="checkbox"/> I10116	Tubular Seat Frame Mini-Flex, 14" Wide <i>Seat frame tubes to mount aftermarket seating systems.</i>	No Charge	

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I100770-99-0	Armrest Left, 8" Leatherette	No Charge	
<input type="checkbox"/> I100774-99-0	Armrest Left, 10" Leatherette	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I100769-99-0	Armrest Right, 8" Leatherette	No Charge	
<input type="checkbox"/> I100773-99-0	Armrest Right, 10" Leatherette	No Charge	

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10137	Height Adjustable Armrest Assembly	307.00	

Armrest Options**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10016	Height Adjustable Armrest Assembly - No Charge <i>No charge with power tilt and/or recline.</i>	No Charge	

Legrest Options**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I101088-99-0	Footplate, Small w/ 5-12" Pan to Ftplate Adj.	No Charge	
<input type="checkbox"/> I101089-99-0	Footplate, Large w/ 5-12" Pan to Ftplate Adj.	No Charge	
<input type="checkbox"/> I100599-99-0	Foot Strap Kit	77.00	

Accessories**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I100518-99-0	Positioning Belt <i>Standard positioning belt is included in the basic equipment package.</i>	No Charge	
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	254.00	
<input type="checkbox"/> I100748-99-0	Miniflex Headrest Removable - Fabric <i>Must also select PN I10552.</i>	422.00	
<input type="checkbox"/> I100749-99-0	Miniflex Headrest Removable - Leatherette <i>Must also select PN I10552.</i>	422.00	
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	213.00	
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	113.00	SEHND0020
<input type="checkbox"/> I100753-99-0	Thigh Support Kit Leatherette <i>Must also select PN I10550.</i>	422.00	
<input type="checkbox"/> I100752-99-0	Thigh Support Kit Fabric <i>Must also select PN I10550.</i>	438.00	
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	113.00	
<input type="checkbox"/> I100751-99-0	Lateral Supports Leatherette <i>Must also select PN I10551.</i>	422.00	
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	113.00	
<input type="checkbox"/> I100867-99-0	Knee Block Leatherette	481.00	
<input type="checkbox"/> I100866-99-0	Knee Block Fabric	523.00	
<input type="checkbox"/> I101473-99-0	Medical Necessity Bag Clips - Carabiner Style	77.00	
<input type="checkbox"/> I100598-99-0	Chest Harness <i>Chest harness kit includes a standard positioning belt.</i>	382.00	
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	389.00	SETND0010
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	113.00	

Stealth Accessories**OPTIONAL ITEMS:**

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10718	Stealth Comfort Plus Headrest - 8"W (CP480) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	358.00	
<input type="checkbox"/> I10721	Adj Stealth Headrest Hardware - Each (TWB485)	242.00	
<input type="checkbox"/> I10736	Stealth Swng Away Egg Sw for Comfort Plus Headrest <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the right or the left side of the head array below:</i>	551.00	
	Mount Swing Away Egg Switch on Right: _____		
	Mount Swing Away Egg Switch on Left: _____		

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	2,590.00	
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	1,264.00	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	1,446.00	

CHOOSE MOUNTING HARDWARE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104690-99-0	Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	234.00	
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	746.00	
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	796.00	
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i>	568.00	
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i>	568.00	
<input type="checkbox"/> I10754	Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 10cm.</i>	568.00	
<input type="checkbox"/> I10068	Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i>	4,323.00	

ALT JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode & On/Off. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	236.00	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 130 g</i>	428.00	

HEAD ARRAY SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10759	Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 14cm Wide x 14cm Tall. Requires Omni Input Module.</i>	4,300.00	
<input type="checkbox"/> I10639	Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below Mount on Left: _____ Mount on Right: _____</i>	500.00	
<input type="checkbox"/> I10760	Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 14cm Wide x 14cm Tall.</i>	4,506.00	
<input type="checkbox"/> I10761	Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i>	236.00	

SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,519.00	WEX5

Alternative Drive Systems

SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,022.00	WEX5
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	236.00	

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10763	VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i>	2,697.00	
<input type="checkbox"/> I10764	VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	236.00	
<input type="checkbox"/> I10765	VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	746.00	
<input type="checkbox"/> I10317	Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 30cm Wide x 24cm Long x 4cm Thick Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,731.00	
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	2,667.00	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	236.00	
<input type="checkbox"/> I10638	RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	5,021.00	
<input type="checkbox"/> I10770	Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i>	236.00	

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10073	Omni Plus Control Module for Pilot Plus <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	2,117.00	
<input type="checkbox"/> I10577	Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i>	113.00	

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Switchbox on Left	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Switchbox on Right	No Charge	

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	152.00	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	114.00	
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	111.00	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	530.00	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	901.00	
<input type="checkbox"/> I10771	Touch Contact with CR2032 Battery Qty: _____	255.00	

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting</i>	472.00	

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The ADP codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*