


# M300 PS



Display ADP Codes on order/quote:

Three column pricing on order/quote:

Please send order/quote to fax#: (800) 231-3256

Email to: [sales@permobilus.com](mailto:sales@permobilus.com)

## Dealer Information

\* Contact: \_\_\_\_\_

Dealer Code: \_\_\_\_\_

\* Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

\* City: \_\_\_\_\_

\* Province: \_\_\_\_\_

\* Phone#: \_\_\_\_\_

\* Fax#: \_\_\_\_\_

PO#: \_\_\_\_\_

Email Address: \_\_\_\_\_

## Client Information

*Permobil recommends that the client is evaluated by a certified rehab specialist.*

\* First Name: \_\_\_\_\_

\* Last Name: \_\_\_\_\_

Diagnosis: \_\_\_\_\_

Funding Source: \_\_\_\_\_

Client Age: \_\_\_\_\_

## Client Measurements

\* Weight: \_\_\_\_\_

\* Height: \_\_\_\_\_

A) Top of Shoulders: \_\_\_\_\_

B) Chest Depth: \_\_\_\_\_

C) Chest Width: \_\_\_\_\_

D) Seat Depth: \_\_\_\_\_

E) Top of Head: \_\_\_\_\_

F) Elbow to Hand: \_\_\_\_\_

G) Seat Pan to Elbow: \_\_\_\_\_

H) Hip Width: \_\_\_\_\_

\*I) Knee to Foot: \_\_\_\_\_

**Permobil Ltd.**  
4 - 75 Mary St.  
Aurora, Ontario  
L4G 1G3

Tel: (888) 726-9330

CAN Fax: (905) 726-9990

US Fax: (800) 231-3256

[www.permobil.com](http://www.permobil.com)

Canadian prices effective March 23, 2012 and payable in Canadian dollars.

## Chair Model and Colors

### CHOOSE BASE:

<input checked="" type="checkbox"/>	I105230-99-0	M300 Base PS - M300 SPO	6,895.00	WP0000455
		<i>Price Includes: M300 Power Base, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. Weight Capacity 300lbs. STFH = 17.5".</i>		
<input type="checkbox"/>	I10727	7.5mph High Speed Upgrade Package - M300	2,148.00	
		<i>Price Includes 7.5mph High Speed Motor Package. Available with R-net Electronics only.</i>		

### CHOOSE DRIVE TIRES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge	
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge	

### OPTIONAL DRIVE TIRES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8"	448.00	
	<i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels.</i>		

### CHOOSE SHROUD COLOR:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105239-99-0	M300 Shroud - Onyx Black	No Charge	
<input type="checkbox"/> I105240-99-0	M300 Shroud - Metallic Green	No Charge	
<input type="checkbox"/> I105241-99-0	M300 Shroud - Vintage Wine	No Charge	
<input type="checkbox"/> I105243-99-0	M300 Shroud - Platinum Silver	No Charge	
<input type="checkbox"/> I105242-99-0	M300 Shroud - Sublime Green	No Charge	
<input type="checkbox"/> I105244-99-0	M300 Shroud - Titan Blue	No Charge	

### CHOOSE HUBCAP COLOR:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Colour	No Charge	
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver	No Charge	

## Base Options

### MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I101994-99-0	Seat Tube Fixed	No Charge	
<input type="checkbox"/> I103779-99-0	Power Adjustable Seat Height - 8" Travel	2,570.00	
	<i>To assist in transfers and improve access for function.</i>		

### MUST CHOOSE BATTERIES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> IM24 SLD G FT	Batteries, Grp 24 (72Ah) Sealed LA/Gel, Installed	824.00	
<input type="checkbox"/> IM34 SLDG	Batteries, Grp 34 (60Ah) Sealed LA/Gel, Installed	814.00	
<input type="checkbox"/> I10828	Seat Stabilizer Option - PS Seat	60.00	
<input type="checkbox"/> I105428-99-0	Enhanced Steering Perform Unit (ESP), R-net	1,288.00	
	<i>Recommended item with switched drive inputs for improved tracking.</i>		
<input type="checkbox"/> I105297-99-0	LED Light Kit for M300 / M400 w/ R-net Electronics	985.00	
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	448.00	

## Permlock

## Permolock

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to: NAME: _____ ADDR: _____</i>	1,828.00	
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	613.00	

## Joystick Options

### CHOOSE JOYSTICK TYPE:

Part Number	Description	MSRP - CA	ADP Code
	PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input checked="" type="checkbox"/>			
<input type="checkbox"/> ID51185	R-net Remote Colour JS w/ Mono Jacks, No Lights	1,540.00	
<input type="checkbox"/> ID51188	R-net Remote Colour JS w/ Toggle Switches, No Lghts	1,746.00	
<input type="checkbox"/> ID51153	R-net Remote Colour JS w/ Mono Jacks, Lights	1,540.00	WEX1
<input type="checkbox"/> ID51149	R-net Remote Colour JS w/ Toggle Switches, Lights	1,700.00	WEX1
<input type="checkbox"/> ID51190	R-net Remote Colour JS w/ Spd Pot/Toggle, Lights	1,746.00	WEX1
<input type="checkbox"/> I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimble <i>Heavy-Duty modification enables gimble to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	1,745.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Contains two mono jacks for mode &amp; on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	1,700.00	
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	216.00	
<input type="checkbox"/> I105246-99-0	R-net LED Remote Joystick Module	No Charge	

### CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	102.00	
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	150.00	
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	102.00	
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	137.00	
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	137.00	
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	137.00	

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10424	Fixed Joystick Mount - Right - VR2	No Charge	
<input type="checkbox"/> I10425	Fixed Joystick Mount - Left VR2	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	335.00	WEPM
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	335.00	WEPM
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2	386.00	WEPM
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2	386.00	WEPM

## Joystick Options

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net <i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	335.00	WE3D
<input type="checkbox"/> I10793	Stealth Gatlin Mid Line Mount - 8 Position <i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	589.00	
<input type="checkbox"/> I10794	Stealth Gatlin Mid Line Mount - 2 Position <i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest: Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>	589.00	
<input type="checkbox"/> I10803	Stealth Gatlin Bracket for R-net Joystick	242.00	
<input type="checkbox"/> I10804	Stealth Gatlin Bracket for Compact Joystick	242.00	
<input type="checkbox"/> I10805	Stealth Gatlin Bracket for Mini Joystick	242.00	
<input type="checkbox"/> I10806	Stealth Gatlin Bracket for Omni	242.00	

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10411	Expandable Controller - R-net <i>Must be ordered with the expandable/remote joysticks (PN: I10410, I103544-99-0, ID51040, ID51108, ID51042, and I103814-99-0).</i>	No Charge	
<input type="checkbox"/> I10513	Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450)</i>	No Charge	
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i>	1,117.00	
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	901.00	

## Seating System

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104738-99-0	PS Seating System - No Tilt	No Charge	
<input type="checkbox"/> I10724	PS Seating System with 50° CG Power Tilt	6,072.00	

## Seating System Options

### CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10289	Single Seat Function Control Switch - LiberT <i>This device is a single toggle switch used to control only one seat function. If control of the seat function through the drive control is needed, then the drive electronics must be upgraded to Pilot+ or R-net Electronics, and a Single Seat Function Control Kit (Part Number I10129), must be selected.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,419.00	
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	2,364.00	SEC1P1005

## Seating System Options

### CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I102798-99-0	Push Buttons on Seat Function Control Kit	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on Seat Function Control Kit	252.00	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	

### CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	129.00	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105026-99-0	Mount ICS Switchbox on Right	No Charge	
<input type="checkbox"/> I105027-99-0	Mount ICS Switchbox on Left	No Charge	

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,117.00	

## Backrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10144	Sling Upholstery Back Support - Standard <i>Width will be the same as seat width.</i>	No Charge	
<input type="checkbox"/> I10145	Curved Metal Back Support 22" High	422.00	
<input type="checkbox"/> I10146	Curved Metal Back Support 22" High - No Charge <i>No charge if ordering tilt.</i>	No Charge	
<input type="checkbox"/> I10454	Omit Backrest	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10148	Angle Adjustable Back Canes <i>Angle adjustable from 85 to 115 degrees.</i>	131.00	
<input type="checkbox"/> I10149	Angle Adjustable Back Canes - No Charge <i>No charge if ordering tilt.</i>	No Charge	

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10150A	Seat Width 16"	No Charge	
<input type="checkbox"/> I10151A	Seat Width 18"	No Charge	
<input type="checkbox"/> I10219A	Seat Width 20"	No Charge	
<input type="checkbox"/> I10220A	Seat Width 22"	No Charge	

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10156	Seat Depth 16"	No Charge	
<input type="checkbox"/> I10157	Seat Depth 17"	No Charge	
<input type="checkbox"/> I10158	Seat Depth 18"	No Charge	
<input type="checkbox"/> I10159	Seat Depth 19"	No Charge	

## Seat Sizing

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10221	Seat Depth 20"	No Charge	
<input type="checkbox"/> I10222	Seat Depth 21"	No Charge	
<input type="checkbox"/> I10223	Seat Depth 22"	No Charge	

## Armrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10163	Height Adjustable Armrest Assembly PS <i>Adjustable in height, with a flip-back design.</i>	293.00	
<input type="checkbox"/> I10230	Height Adjustable Armrest Assembly PS - No Charge <i>Adjustable in height, with a flip-back design. No charge if ordering power tilt.</i>	No Charge	

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10290	PS Armrest Pads - 10.5"L x 2"W (Pair)	No Charge	
<input type="checkbox"/> I10291	PS Armrest Pads - 14"L x 2"W (Pair)	No Charge	
<input type="checkbox"/> I10292	Waterfall Armrest Pad for PS - 13.5"L x 3"W (Pair)	81.00	
<input type="checkbox"/> I10168	Omit Armrest Pads	No Charge	
<input type="checkbox"/> I10166	Corpus Armrest Pad Adapter Brackets <i>Must be picked if ordering Corpus armrest pads.</i>	52.00	
<input type="checkbox"/> I10208	Corpus Armrest Pad 4"x10" Left Leatherette	95.00	
<input type="checkbox"/> I10209	Corpus Armrest Pad 4"x12" Left Leatherette	95.00	
<input type="checkbox"/> I10210	Corpus Armrest Pad 4"x14" Left Leatherette	95.00	
<input type="checkbox"/> I10211	Corpus Armrest Pad 4"x16" Left Leatherette	95.00	
<input type="checkbox"/> I10212	Corpus Armrest Pad 4"x10" Right Leatherette	95.00	
<input type="checkbox"/> I10213	Corpus Armrest Pad 4"x12" Right Leatherette	95.00	
<input type="checkbox"/> I10214	Corpus Armrest Pad 4"x14" Right Leatherette	95.00	
<input type="checkbox"/> I10215	Corpus Armrest Pad 4"x16" Right Leatherette	95.00	
<input type="checkbox"/> I10169	Armrest Side Support Right Inner <i>Will limit the armrest from flipping back. For Corpus Armrest Pads only.</i>	146.00	
<input type="checkbox"/> I10170	Armrest Side Support Left Inner <i>Will limit the armrests from flipping back. For Corpus Armrest Pads only.</i>	146.00	
<input type="checkbox"/> I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	153.00	
<input type="checkbox"/> I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	153.00	

## Legrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10172	70 Degree Swing Away Legrests (Pair)	205.00	
<input type="checkbox"/> I10231	70 Degree Swing Away Legrests (Pair) <i>No charge if ordering Power Tilt.</i>	No Charge	
<input type="checkbox"/> I10625	Aeron Legrest Assembly for PS Seat - Narrow <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes Narrow, Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	335.00	
<input type="checkbox"/> I10316	Aeron Legrest Assembly for PS Seat - Standard <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes standard width, Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i>	324.00	

## Legrest Options

### CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10173	Footplates (Pair), Standard Small	No Charge	
<input type="checkbox"/> I10174	Footplates (Pair), Standard Large <i>Not available with 16" seat width.</i>	No Charge	
<input type="checkbox"/> I10175	Small Footplates, Multi - Axis Adjustable, Aluminum	201.00	
<input type="checkbox"/> I10176	Large Footplates, Multi - Axis Adjustable, Aluminum <i>Cannot be ordered with 16" or 18" seat widths.</i>	201.00	

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104684-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 &amp; I10587).</i>	No Charge	
<input type="checkbox"/> I104686-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 &amp; I10587).</i>	No Charge	
<input type="checkbox"/> I10177	Heel Loops <i>Not available with Aeron Legrest Assembly. (I10316).</i>	93.00	
<input type="checkbox"/> I10197	Calf Supports Leatherette <i>Not available with lower leg lengths of 10" or less.</i>	No Charge	

## Accessories

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10401	Bodypoint Non-Padded Hip Belt, PS - Medium	No Charge	
<input type="checkbox"/> I10406	Bodypoint Padded Hip Belt PS (C300) - Medium	123.00	
<input type="checkbox"/> I10403	Bodypoint Non-Padded Hip Belt, PS - Large	123.00	
<input type="checkbox"/> I10409	Bodypoint Padded Hip Belt PS - Large	163.00	
<input type="checkbox"/> I10199	Headrest-Removable Leatherette <i>Must also select PN I10552.</i>	422.00	
<input type="checkbox"/> I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	519.00	
<input type="checkbox"/> I10180	Universal Headrest Adapter <i>Will be sent as part if ordering sling back upholstery.</i>	203.00	
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	113.00	SEHND0020
<input type="checkbox"/> I10783	Hand Heater Kit with Buddy Button <i>Kit includes Hand Heater (2"H x 1.5"W x 2.75"D), Buddy Button, Power Cable with XLR Plug and PERMOfix Mounting Adapter. Heater has three temp/fan speed settings. For best results, a hand cone (I10784 or I10785) is recommended to trap heat around the user's hand(s). To ensure best positioning for the user, heater will NOT be mounted when chair is received.</i>	878.00	
<input type="checkbox"/> I10784	One-Hand Cone Kit with PERMOfix Mounting Set <i>Clear plastic dome (5.5"W x 13.5"L x 8"H) with forearm opening at one end. Kit includes PERMOfix mounting hardware to attach cone to end of armrest. To ensure best positioning for the user, cone will NOT be mounted when wheelchair is received.</i>	507.00	
<input type="checkbox"/> I10785	Two-Hand Cone Kit with Mounting Pins <i>Clear plastic triangular dome (15"W x 12.5"L x 8"H) with two forearm openings. Use with Upper Extremity Support is recommended. Kit includes four pins to secure cone to tray surface. Pins remove without tools for easy transfers. To ensure the best positioning for the user, cone will NOT be mounted.</i>	399.00	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	232.00	
<input type="checkbox"/> I10179	Medical Necessity Bag Clips - Carabiner Style <i>Will be sent as part if ordering sling back upholstery.</i>	77.00	

## Accessories

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I105144-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i>	160.00	
<input type="checkbox"/> I105143-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i>	160.00	
<input type="checkbox"/> I105145-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i>	160.00	
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	389.00	SETND0010
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	113.00	
<input type="checkbox"/> I101964-99-0	Oxygen Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest.</i>	242.00	WEPX
<input type="checkbox"/> I101956-99-0	Crutch Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder. This item is mounted to the ergonomic backrest.</i>	242.00	

## Stealth Accessories

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10808	Stealth Ultra Headrest <i>Includes Stealth Occipital (500 series) and Sub Occipital (600 series) Pads, bilateral Facial Support Pads (#900) with Stealth Swing Away Hardware (SUS9-B) and Flip Down Headrest Mount (TWB480-FDM). Also includes Universal Headrest Adapter.</i>	1,425.00	
<input type="checkbox"/> I10800	Stealth Swing Away Egg Switch for Ultra Headrest <i>Includes Black Egg Switch and Stealth Ultra Swing Away Hardware (SUS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i>	550.00	
<input type="checkbox"/> I10801	Stealth Swing Away Facial Supprt Hdw for Ultra <i>Must also select support pad (I10810, I10811, I10812). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	325.00	
<input type="checkbox"/> I10810	Stealth #900 Support Pad (4.5" x 2") - Ultra	127.00	
<input type="checkbox"/> I10811	Stealth #910 Support Pad (3" x 2") - Ultra	127.00	
<input type="checkbox"/> I10812	Stealth #920 Support Pad (6.4" x 2") - Ultra	127.00	
<input type="checkbox"/> I10736	Stealth Swng Away Egg Sw for Comfort Plus Headrest <i>Includes Black Egg Switch and Stealth Comfort Plus Swing Away Hardware (CPS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i>	551.00	
<input type="checkbox"/> I10809	Stealth Swing Away Facial Supprt Hdw for Cmfrt Plus <i>Must also select support pad (I10813, I10814, I10815). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	325.00	
<input type="checkbox"/> I10813	Stealth #900 Support Pad (4.5" x 2") - Cmfrt Plus	127.00	
<input type="checkbox"/> I10814	Stealth #910 Support Pad (3" x 2") - Cmfrt Plus	127.00	
<input type="checkbox"/> I10815	Stealth #920 Support Pad (6.4" x 2") - Cmfrt Plus	127.00	

## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 13 g Neutral to Full Deflection = 5mm</i>	2,590.00	

## Alternative Drive Systems

### ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 220 g Neutral to Full Deflection = 28mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	1,264.00	
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 50 g Neutral to Full Deflection = 17mm Joystick includes two 3.5mm mono jacks for Mode and On/Off.</i>	1,446.00	

### CHOOSE MOUNTING HARDWARE:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104690-99-0	Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	234.00	
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	746.00	
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	796.00	
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 19cm.</i>	568.00	
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 14cm.</i>	568.00	
<input type="checkbox"/> I10754	Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 10cm.</i>	568.00	
<input type="checkbox"/> I10068	Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i>	4,323.00	
<input type="checkbox"/> I10795	Stealth Swing Away Chin Boom for Ultra Headrest <i>Includes Stealth Ultra Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest: Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____</i>	523.00	
<input type="checkbox"/> I10797	Stealth Sw Away Chin Boom for Comfort Plus Headrest <i>Includes Stealth Comfort Plus Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest: Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____</i>	523.00	
<input type="checkbox"/> I10807	Stealth Gatlin Eclipse Tray w/ Mini JS Mount <i>Designed to mount on Stealth Gatlin Mid Line Mounts. Includes bracket for mounting a Mini Joystick under tray. Through hole in tray surface provides joystick access. Please indicate side of tray for hole position: Position Hole on Left: _____ Position Hole on Right: _____</i>	633.00	

### ALT JOYSTICK ACCESSORIES:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid Switches can be used for Mode &amp; On/Off. Switch operation force = 130 g Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	236.00	
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode &amp; On/Off. Switch operation force = 130 g</i>	428.00	

### HEAD ARRAY SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10759	Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 14cm Wide x 14cm Tall. Requires Omni Input Module.</i>	4,300.00	

## Alternative Drive Systems

### HEAD ARRAY SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10639	Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below Mount on Left: _____ Mount on Right: _____</i>	500.00	
<input type="checkbox"/> I10760	Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 14cm Wide x 14cm Tall.</i>	4,506.00	
<input type="checkbox"/> I10761	Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i>	236.00	

### SIP & PUFF SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,519.00	WEX5
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,022.00	WEX5
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	236.00	

### ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10763	VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i>	2,697.00	
<input type="checkbox"/> I10764	VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	236.00	
<input type="checkbox"/> I10765	VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	746.00	
<input type="checkbox"/> I10317	Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 30cm Wide x 24cm Long x 4cm Thick Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,731.00	
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	2,667.00	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	236.00	
<input type="checkbox"/> I10638	RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	5,021.00	
<input type="checkbox"/> I10770	Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i>	236.00	

## Alt Drive Systems Setup

### CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	2,117.00	
<input type="checkbox"/> I10577	Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i>	113.00	

## Alt Drive Systems Setup

### CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,494.00	

### CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Switchbox on Left	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Switchbox on Right	No Charge	

## Alt Drive Accessories

### OPTIONAL ITEMS:

Part Number	Description	MSRP - CA	ADP Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	152.00	
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	114.00	
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	111.00	
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	530.00	
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	901.00	
<input type="checkbox"/> I10771	Touch Contact with CR2032 Battery Qty: _____	255.00	
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights. Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 150mm rod, one 200mm rod and one 225mm gooseneck switch with 1.2m cable and mounting clamp.</i>	472.00	

Notes & Additional Instructions:

**Chair Order Policy:**

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, LTD. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

**Order Acknowledgement:**

I, \_\_\_\_\_, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: \_\_\_\_\_

Printed Name: \_\_\_\_\_

Title: \_\_\_\_\_

Date: \_\_\_\_\_

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The ADP codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*