


Permolock



**Permolock C
(C400 & C500)**



**Permolock C3
(C300, C350, K300, K450,
M300, M400)**

*** Indicates a Required Field**

Dealer Information

*Contact: _____
* Permobil Dealer Code: _____
*NMEDA Certification #: _____
*Dealer Name: _____
Address: _____
City: _____
Province: _____
*Phone#: _____
Fax#: _____
PO#: _____

Client Information

Permobil requires that a Permolock Base be installed by a Certified NMEDA Dealer.

First Name: _____
Last Name: _____
Serial Number of Wheelchair: _____
*Wheelchair Base Type: _____
*Seating System Type: _____

Memo on Transporting Wheelchairs in Vehicles Using Permolock

A user may ride in a Permobil C300, C350, K450, C400, C500, M300 or M400 of appropriate configuration in a motor vehicle, when properly secured by both a Permolock wheelchair lockdown system AND a crash tested lap and shoulder belt, which is secured to the motor vehicle. A "crash tested lap belt" does NOT include the positioning belt attached to the wheelchair. The wheelchair belt is designed to position the user in wheelchair during ordinary use. It is not designed to protect the user against injury while riding in the wheelchair in a vehicle. Users who pursue riding in or driving vehicles in their wheelchairs assume a great risk. We urge them to give careful consideration to this pursuit.

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Prices effective **January 16, 2012**. All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if the order is for a new Permolock. Any parts ordered for service or future alterations will carry different pricing and warranty.

Permobil Ltd.

4 - 75 Mary St.

Aurora, Ontario

L4G 1G3

Tel: (888) 726-9330

CAN Fax: (905) 726-9990

US Fax: (800) 231-3256

www.permobil.com



Part Number	Description	Price CAD
PERMOLOCK C ITEMS (USE WITH C400 & C500)		
<i>The items in this section are compatible with the C400 and C500 bases. Choose the Permolock C when ordering one of these power wheelchairs.</i>		
<input type="checkbox"/> I102843-99-0	Permolock C Kit - COMPLETE <i>Includes: Locking Base, Plunger, Manual Release Handle, and Vehicle Mounted Release Button. The Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer Name & Address where the locking Base will ship to NAME: _____ ADDR: _____</i>	2,560.00
<input type="checkbox"/> I102840-99-0	Base ONLY for Permolock C <i>Price Includes: Locking Base only. This item must be installed by a NMEDA Dealer.</i>	2,266.00
<input type="checkbox"/> I102842-99-0	Plunger ONLY for Permolock C <i>Price Includes: Plunger only. This item may be purchased by any Permobil Inc or NMEDA certified dealer. The standard mounting position for the button to control the plunger is at the end of the armrest pad, opposite the joystick.</i>	1,205.00
<input type="checkbox"/> I103242-99-0	3.6m Plunger Sw Ext Cable for VS with Permolock C <i>This Cable MUST be selected when a Permolock C is ordered with a VS seating system.</i>	
<input type="checkbox"/> I102841-99-0	Manual Release Handle ONLY for Permolock C <i>Price Includes: Manual Release Handle only. This item may be purchased by any Permobil Inc or NMEDA Certified Dealer.</i>	314.00
<input type="checkbox"/> I1824545	Permolock C OEM Seat Install Bracket <i>Allows for installation of original automotive seat.</i>	695.00
PERMOLOCK C3 ITEMS		
<i>The items in this section are compatible with the C300, C350, K300, M300, M400 and K450 bases. Choose the Permolock C3 when ordering one of these power wheelchairs.</i>		
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base C300, C350, K300, K450 <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to: NAME: _____ ADDR: _____</i>	1,828.00
<input type="checkbox"/> I10609	Permolock C3 Dual Locking Pin Kit - C350 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes a set of reinforced braces must be installed whenever the spring loaded pins are mounted in a C350 Base.</i>	613.00
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit - C300/K300 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	613.00
<input type="checkbox"/> I104652-99-0	Permolock C3 Dual Locking Pin Kit - K450 <i>This kit includes two pins that are mounted underneath the base of the wheelchair. The pins are fixed and do not push into the wheelchair base.</i>	613.00
<input type="checkbox"/> I1824275	Permolock C3 OEM Seat Install Bracket <i>Allows for installation of original automotive seat.</i>	541.00
PERMOLOCK ACCESSORIES - FIT ANY PERMOLOCK SYSTEM		
<i>The accessories below are compatible with both the Permolock C and Permolock C3 systems.</i>		
<input type="checkbox"/> I1821551	Permolock Height Adapter for OEM Seat Bracket <i>Used in combination with the OEM Seat Install Bracket to provide height adjustablty of original automotive seat. Range of height adjustability is 8" - 12", in 1" increments.</i>	530.00
<input type="checkbox"/> I104232-99-0	2m Extension Cable for Permolock Base (C or C3) <i>This 2m Extension Cable increases the length of the "power" cable for the Permolock Locking Base, making it possible to position it farther from the dashboard or in a rear passenger position.</i>	70.00

Notes & Additional Information: