



**Permolock C  
(C400 & C500)**



**Permolock C3  
(C300, C350, K300, K450,  
M300, M400)**

**\* Indicates a Required Field**  
Dealer Information

\*Contact: \_\_\_\_\_

\* Permobil Dealer Code: \_\_\_\_\_

\*NMEDA Certification #: \_\_\_\_\_

\*Dealer Name: \_\_\_\_\_

Address: \_\_\_\_\_

City: \_\_\_\_\_

State/Zip: \_\_\_\_\_

\*Phone#: \_\_\_\_\_

Fax#: \_\_\_\_\_

PO#: \_\_\_\_\_

Client Information

*Permobil requires that a Permolock Base be installed by a NMEDA Dealer.*

First Name: \_\_\_\_\_

Last Name: \_\_\_\_\_

Serial Number of Wheelchair: \_\_\_\_\_

\*Wheelchair Base Type: \_\_\_\_\_

\*Seating System Type: \_\_\_\_\_

**Memo on Transporting Wheelchairs in Vehicles Using Permolock**

A user may ride in a Permobil C300, C350, K450, C400, C500, M300 or M400 of appropriate configuration in a motor vehicle, when properly secured by both a Permolock wheelchair lockdown system AND a crash tested lap and shoulder belt, which is secured to the motor vehicle. A "crash tested lap belt" does NOT include the positioning belt attached to the wheelchair. The wheelchair belt is designed to position the user in wheelchair during ordinary use. It is not designed to protect the user against injury while riding in the wheelchair in a vehicle. Users who pursue riding in or driving vehicles in their wheelchairs assume a great risk. We urge them to give careful consideration to this pursuit.

Please send order/quote to fax#: (800) 231-3256  
Email to: [sales@permobilus.com](mailto:sales@permobilus.com)

Prices effective **January 16, 2012**. All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if the order is for a new Permolock. Any parts ordered for service or future alterations will carry different pricing and warranty.

**Permobil Inc.**  
300 Duke Dr.  
Lebanon, TN 37090  
Tel: (800) 736-0925  
Fax: (800) 231-3256

Part Number	Description	Price
	<b>PERMOLOCK C ITEMS (USE WITH C400 &amp; C500)</b> <i>The items in this section are compatible with the C400 and C500 bases. Choose the Permolock C when ordering one of these power wheelchairs.</i>	
<input type="checkbox"/> I103049-99-0	Permolock C Locking Base Kit <i>Includes: Locking Base, Manual Release Handle and Vehicle Mounted Release Button. PN I102842-99-0 must be selected below. The Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer Name &amp; Address where the locking Base will ship to</i> NAME: _____ ADDR: _____	1,600.00
<input type="checkbox"/> I102842-99-0	Plunger ONLY for Permolock C <i>Price Includes: Plunger only. This item may be purchased by any Permobil Inc or NMEDA certified dealer. The standard mounting position for the button to control the plunger is at the end of the armrest pad, opposite the joystick.</i>	1,057.00
<input type="checkbox"/> I103242-99-0	3.6m Plunger Sw Ext Cable for VS with Permolock C <i>This Cable MUST be selected when a Permolock C is ordered with a VS seating system.</i>	No Charge
<input type="checkbox"/> I102839-99-0	Permolock C OEM Seat Install Bracket <i>Allows for installation of original automotive seat onto Permolock C Locking Base.</i>	666.00
	<b>PERMOLOCK C3 ITEMS</b> <i>The items in this section are compatible with the C300, C350, K300 M300, M400 and K450 bases. Choose the Permolock C3 when ordering one of these power wheelchairs.</i>	
<input type="checkbox"/> I104680-99-0	Permolock C3 Locking Base Kit <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to: NAME: _____</i> ADDR: _____	1,860.00
<input type="checkbox"/> I10608	Permolock C3 Dual Locking Pin Kit- C/K300, M-series <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permolock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i>	623.00
<input type="checkbox"/> I10609	Permolock C3 Dual Locking Pin Kit - C350 <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes a set of reinforced braces must be installed whenever the spring loaded pins are mounted in a C350 Base.</i>	623.00
<input type="checkbox"/> I104652-99-0	Permolock C3 Dual Locking Pin Kit - K450 <i>This kit includes two pins that are mounted underneath the base of the wheelchair. The pins are fixed and do not push into the wheelchair base.</i>	623.00
<input type="checkbox"/> I104708-99-0	Permolock C3 OEM Seat Install Bracket <i>Allows for installation of original automotive seat onto Permolock C3 Locking Base.</i>	550.00
	<b>PERMOLOCK ACCESSORIES - FIT ANY PERMOLOCK SYSTEM</b> <i>The accessories below are compatible with both the Permolock C and Permolock C3 systems.</i>	
<input type="checkbox"/> I1821551	Permolock Height Adapter for OEM Seat Bracket <i>Used in combination with the OEM Seat Install Bracket to provide height adjustablty of original automotive seat. Range of height adjustability is 8" - 12", in 1" increments.</i>	540.00
<input type="checkbox"/> I104232-99-0	2m Extension Cable for Permolock Base (C or C3) <i>This 2m Extension Cable increases the length of the "power" cable for the Permolock Locking Base, making it possible to position it farther from the dashboard or in a rear passenger position.</i>	68.00
<input type="checkbox"/> I102859-99-0	Quick Release Kit for Permolock Manual Release <i>The Locking Base's Manual Release Handle mounts to this kit. Use when quick removal of the Manual Release is desired.</i>	127.00

Notes & Additional Information: