


C500 Corpus 3G



*** Indicates a Required Field**

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State/Zip: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____

E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* I) Knee to Foot: _____

Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Permobil Inc.
300 Duke Dr.
Lebanon, TN 37090
Tel: (800) 736-0925
Fax: (800) 231-3256
www.permobil.com

Prices effective **January 16, 2012.**

C500 Corpus 3G

Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105036-99-0	C500 Base Corpus 3G 5mph - C500 PS2 <i>Price Includes: FWD C500 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. Weight Capacity = 300lbs, STFH = 18.75".</i>	11,943.00	K0884/K0014
<input type="checkbox"/> I105037-99-0	C500 Base Corpus 3G 6.5mph - C500 PS2 <i>Price Includes: FWD C500 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Anti-tippers, ESP (Enhanced Steering Performance), High Speed Motors, 120 Amp Electronics, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. Weight Capacity = 300lbs, STFH = 18.75". Anti-tippers can not be removed.</i>	13,127.00	K0884/K0014

CHOOSE DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103499-99-0	Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8"	No Charge	
<input type="checkbox"/> I103498-99-0	Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8"	No Charge	

OPTIONAL DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10782	Tri-Spoke Split Rim w/ Snow Tire Pneumatic 3x8" <i>Snow tires are designed for use in cold weather conditions. The nylon compound in these tires remains flexible at low temperatures, but wears more easily. Therefore, the snow tires (sold as wheel assemblies) are not recommended for year-long use and must be ordered along with standard drive wheels.</i>	407.00	

CHOOSE SHROUD COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103816-99-0	C500 Shroud Onyx Black	No Charge	
<input type="checkbox"/> I102260-99-0	C500 Shroud Vintage Wine	No Charge	
<input type="checkbox"/> I102261-99-0	C500 Shroud Crystal Blue	No Charge	
<input type="checkbox"/> I102262-99-0	C500 Shroud Metallic Charcoal Gray	No Charge	
<input type="checkbox"/> I102263-99-0	C500 Shroud Green Metallic	No Charge	
<input type="checkbox"/> I103833-99-0	C500 Shroud Platinum Silver	No Charge	

CHOOSE HUBCAP COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Color	No Charge	
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver	No Charge	

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105140-99-0	Power Adjustable Seat Height - 8" Travel	2,724.00	E2300
<input type="checkbox"/> I105405-99-0	Seat Tube Fixed	No Charge	

MUST CHOOSE BATTERIES:

Part Number	Description	Price	HCPCS Code
<input checked="" type="checkbox"/> IM24 SLD G FT	Batteries, Grp 24 (72Ah) Sealed LA/Gel, Installed	824.00	E2363

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103635-99-0	Light Kit, R-net, C400/C500	1,001.00	
<input checked="" type="checkbox"/> I104082-99-0	R-net ESP w C400/C500 Mounting Bracket <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	No Charge	K0108
<input type="checkbox"/> I10715	Compact Laptop Ventilator Tray - Corpus 3G Seat <i>Mounts to rear of Corpus 3G backrest hinge. Accommodates Pulmonetic Systems / CareFusion LTV Series and similar sized Laptop Ventilators. The tray is 13"T, adjustable in width from 8" - 13", and adjustable in depth from 2" - 6.5". Overall backrest angle limited to 150° when mounting this ventilator tray. Custom Option.</i>	1,385.00	E1029

C500 Corpus 3G

Base Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10787	Battery Tray for Laptop Ventilator Tray <i>Measures 11 3/8" W x 6 1/2" H x 1 5/8" D. Accomodates Pulmonetic Systems Lithium-Ion Battery Pack, CareFusion SprintPack Battery System, or Richardson PowerTech Vent Power Center. Custom Option.</i>	208.00	K0108
<input type="checkbox"/> I10788	O2 Holder for Laptop Ventilator Tray <i>Mounts to Left or Right Side of Laptop Vent Tray. O2 Holder is 7" Tall. Width infinitely adjustable 3.5" - 7". Depth infinitely adjustable 3.75" - 5". Custom Option.</i>	295.00	K0108
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	456.00	
<input type="checkbox"/> I105304-99-0	Anti-Tippers for C400 / C500 <i>Anti-Tippers are required for any user weight above 220lbs. They are included in the High Speed Package for the C500 VS Base (I104060-99-0).</i>	418.00	

Permolock

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102843-99-0	Permolock C Kit - COMPLETE <i>Includes: Locking Base, Plunger, Manual Release Handle, and Vehicle Mounted Release Button. The Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer Name & Address where the locking Base will ship to NAME: _____ ADDR: _____</i>	2,604.00	
<input type="checkbox"/> I102842-99-0	Plunger ONLY for Permolock C <i>Price Includes: Plunger only. This item may be purchased by any Permobil Inc or NMEDA certified dealer. The standard mounting position for the button to control the plunger is at the end of the armrest pad, opposite the joystick.</i>	1,057.00	

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>For additional input devices, refer to Alternative Drive Systems Section.</i>		
<input type="checkbox"/> ID51185	R-net Remote Color JS w/ Mono Jacks, No Lights	716.00	
<input type="checkbox"/> ID51188	R-net Remote Color JS w/ Toggle Switches, No Lghts	917.00	
<input type="checkbox"/> ID51153	R-net Remote Color JS w/ Mono Jacks, Lights	716.00	
<input type="checkbox"/> ID51149	R-net Remote Color JS w/ Toggle Switches, Lights	917.00	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Spd Pot/Toggle, Lights	917.00	
<input type="checkbox"/> I105700-99-0	R-net Remote Color JS w/ Lights, Heavy Duty Gimble <i>Heavy-Duty modification enables gimble to hold up to extreme operating forces. Reinforced housing is attached to a standard Color LCD R-net Joystick. 3" long joystick knob has 3/4" diameter. Not compatible with alternative joystick handles.</i>	972.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>	634.00	E2373
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	204.00	

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	99.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	145.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	99.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	133.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	133.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	133.00	E2323
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	133.00	E2323
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	133.00	E2323

C500 Corpus 3G

Joystick Options

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	133.00	E2323

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10664	Joystick Mount Right - Fixed	No Charge	
<input type="checkbox"/> I10665	Joystick Mount Left - Fixed	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	306.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	306.00	E1028
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2	510.00	E1028
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2	510.00	E1028
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net	276.00	E1028
	<i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>		
<input type="checkbox"/> I10793	Stealth Gatlin Mid Line Mount - 8 Position	535.00	E1028
	<i>Features eight locking positions each 45° apart. Adjustments are made by pulling out on mounting arm and rotating to desired angle. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest:</i>		
	<i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>		
<input type="checkbox"/> I10794	Stealth Gatlin Mid Line Mount - 2 Position	535.00	E1028
	<i>Designed for easy swing away. Mount can be set to swing from horizontal to 95° or from horizontal to 180°. Joystick bracket sold separately. Please indicate if you want it mounted on left or right armrest:</i>		
	<i>Mount Midline on Left Armrest: _____ Mount Midline on Right Armrest: _____</i>		
<input type="checkbox"/> I10803	Stealth Gatlin Bracket for R-net Joystick	220.00	E1028
<input type="checkbox"/> I10804	Stealth Gatlin Bracket for Compact Joystick	220.00	E1028
<input type="checkbox"/> I10805	Stealth Gatlin Bracket for Mini Joystick	220.00	E1028
<input type="checkbox"/> I10806	Stealth Gatlin Bracket for Omni	220.00	E1028
<input type="checkbox"/> I10802	Stealth Height Adjustable Mnt (ARM250-P)	195.00	E1028
	<i>Provides 1.75" of height adjustability to all Permobil joystick mounts and Stealth Gatlin Mid Line Mounts.</i>		

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10411	Expandable Controller - R-net	561.00	E2377
	<i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i>		
<input type="checkbox"/> I10513	Harness for Expandable Controller	489.00	E2313
	<i>Must be ordered with any Expandable Controller (I10411 or I10450).</i>		
<input type="checkbox"/> I105225-99-0	Co-Pilot Attendant Control R-net - Corpus 3G	4,375.00	E2331
<input type="checkbox"/> I105408-99-0	Attendant Control R-net - Corpus 3G	1,136.00	E2331
<input type="checkbox"/> ID51111	R-net BT Mouse Module	833.00	
	<i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>		
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System	1,042.00	E2351
	<i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>		

Seating System

CHOOSE SEAT MOUNTING PLATE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105056-99-0	Standard (0°) Tilt Mounting Plate	No Charge	
	<i>Range of Tilt will be 0° to +50°.</i>		

C500 Corpus 3G

Seating System

CHOOSE SEAT MOUNTING PLATE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105068-99-0	5° Anterior Tilt Mounting Plate <i>Range of Tilt will be -5° to +45°. This option lowers the front Seat to Floor Height by 1".</i>	232.00	

CHOOSE SEATING SYSTEM:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE: <i>Weight Limit 300lbs. For the Corpus 3G Seat on the C400/C500, there is a half speed limitation when the Total Backrest Angle (TBA) > 150°. TBA = Tilt Angle + Backrest Recline Angle. Example 1: 45° Tilt + 100° Recline = 145° TBA (Full Speed) Example 2: 35° Tilt + 120° Recline = 155° TBA (Half Speed)</i>			
<input type="checkbox"/> I10647	C3G Seat w/ No Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>	5,995.00	K0108
<input type="checkbox"/> I10732	C3G Seat w/ No Tilt & 175° Power Recline	7,850.00	E1004
<input type="checkbox"/> I10733	C3G Seat w/ 50° DuraCore Tilt & Manual Recline <i>Manual Recline angle limited to 120°.</i>	6,640.00	E1002
<input type="checkbox"/> I10734	C3G Seat w/ 50° DuraCore Tilt & 175° Power Recline	11,395.00	E1007

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,540.00	E2310
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	3,013.00	E2311

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox	256.00	
<input type="checkbox"/> I105409-99-0	Do Not Install ICS Switchbox, Send as Part	No Charge	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I104651-99-0	Mount ICS Switchbox onto Armrest <i>Mount on Left: _____ Mount on Right: _____</i>	No Charge	
<input type="checkbox"/> I105026-99-0	Mount ICS Switchbox on Right	No Charge	
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest <i>Mount on Left: _____ Mount on Right: _____</i>	131.00	
<input type="checkbox"/> I105027-99-0	Mount ICS Switchbox on Left	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,032.00	K0108

C500 Corpus 3G

Seating System Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104769-99-0	Remote Stop System - All Electronics <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	1,032.00	
<input type="checkbox"/> I10750	Memory Seat Program, Corpus 3G w R-net	103.00	

Backrest Options

CHOOSE FIXED OR SLIDING BACKREST KIT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105067-99-0	Sliding Backrest Kit - 5" Mechanical Shear	No Charge	
<input type="checkbox"/> I105066-99-0	Fixed Backrest Kit - No Mechanical Shear	30.00	

CHOOSE BACKREST TYPE:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

BACKREST HEIGHT NOTE:

Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Ergo Backs include lateral wedges and lumbar pads.

<input type="checkbox"/> I105113-99-0	Corpus 3G Ergo Backrest - 14"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	898.00	E2620
<input type="checkbox"/> I105114-99-0	Corpus 3G Ergo Backrest - 16"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	898.00	E2620
<input type="checkbox"/> I105109-99-0	Corpus 3G Ergo Backrest - 18"W x 20"T (C3GERGBCA) <i>Power Recline will be limited to 150° when ordering this backrest.</i>	898.00	E2620
<input type="checkbox"/> I10667	Solid Back 16"W X 20"T - Corpus 3G Seat <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	
<input type="checkbox"/> I10669	Solid Back 18"W x 20"T - Corpus 3G Seat <i>Power Recline will be limited to 150° when ordering this backrest.</i>	No Charge	

PLEASE NOTE:

All taller backrests listed below include: lateral wedges, lumbar pads, a UT rail mounted at the top for lateral adjustment of the headrest, and split backrest feature with height adjustment in 1" increments. Please indicate below where you want the height set at so that a properly sized backrest cushion can be provided.

<input type="checkbox"/> I105110-99-0	Corpus 3G Ergo Back -14"W x 23 - 28"T (C3GERGBCA)	898.00	E2620
<input type="checkbox"/> I105111-99-0	Corpus 3G Ergo Back -16"W x 23 - 28"T (C3GERGBCA)	898.00	E2620
<input type="checkbox"/> I105112-99-0	Corpus 3G Ergo Back -18"W x 23 - 28"T (C3GERGBCA)	898.00	E2620
<input type="checkbox"/> I10668	Solid Back 16"W x 23 - 28"T - Corpus 3G Seat	No Charge	
<input type="checkbox"/> I10670	Solid Back 18"W x 23 - 28"T - Corpus 3G Seat	No Charge	
	Set Backrest at 23" _____		
	Set Backrest at 24" _____		
	Set Backrest at 25" _____		
	Set Backrest at 26" _____		
	Set Backrest at 27" _____		
	Set Backrest at 28" _____		

CHOOSE BACKREST CUSHION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10659	Leatherette Cushion - Corpus 3G Ergo Back	No Charge	
<input type="checkbox"/> I10660	Black Mesh Cushion - Corpus 3G Ergo Back	30.00	
<input type="checkbox"/> I10661	Gray Mesh Cushion - Corpus 3G Ergo Back	30.00	
<input type="checkbox"/> I105117-99-0	Omit Backrest Cushion	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
-------------	-------------	-------	------------

C500 Corpus 3G

Backrest Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10726	Tubular Backrest Frame - Corpus 3G (1" Canes) <i>Provides 1" diameter, vertical canes for mounting aftermarket backrests. Infinitely adjustable in width from 14" to 20". Height of canes is 21" - 23" measured from the seat pan to the top of the cane. Provides an infinitely adjustable 1"-4" gap between the backrest hinge and the cane. Depending on the aftermarket back, some seat depth may be lost when mounting this system. This is a custom option. Please also indicate the length of the UniTrack linkage disks below: 1" UniTrack Link Kit, C3G Tubular Back: _____ 2" UniTrack Link Kit, C3G Tubular Back: _____ 3" UniTrack Link Kit, C3G Tubular Back: _____</i>	397.00	
<input type="checkbox"/> I10725	Split Backrest Cushion - Corpus 3G <i>Custom Option.</i>	283.00	
<input type="checkbox"/> I10786	Tailored Backrest Cushion - Corpus 3G <i>Sides of cushion are trimmed to provide additional clearance for lateral trunk supports. You must also select a backrest cushion (I10659, I10660 or I10661) when ordering this custom option.</i>	95.00	

Seat Sizing

CHOOSE SEAT FRAME AND WIDTH SETTING:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>Narrow Seat Frames can be adjusted to 17"W, 19"W, or 21"W without any additional parts. Wide Seat Frames can be adjusted to 19"W, 21"W, or 23"W without any additional parts. Please select both the Frame Kit and Width Setting below.</i>		
<input type="checkbox"/> I10673	Narrow Corpus 3G Seat Frame Kit (17", 19", 21")	No Charge	
<input type="checkbox"/> I10674	Wide Corpus 3G Seat Frame Kit (19", 21", 23") Set Seat Width at 17": _____ Set Seat Width at 19": _____ Set Seat Width at 21": _____ Set Seat Width at 23": _____	No Charge	

CHOOSE SEAT DEPTH:

Part Number	Description	Price	HCPCS Code
	PLEASE NOTE: <i>Odd numbered seat depths can be made 1" shorter without any additional parts. For example, a 19"D Seat can also be set at 18"D.</i>		
<input type="checkbox"/> I105040-99-0	Seat Depth 14"	No Charge	
<input type="checkbox"/> I105457-99-0	Seat Depth 15"	No Charge	
<input type="checkbox"/> I105095-99-0	Seat Depth 16"	No Charge	
<input type="checkbox"/> I105458-99-0	Seat Depth 17"	No Charge	
<input type="checkbox"/> I105042-99-0	Seat Depth 18"	No Charge	
<input type="checkbox"/> I105459-99-0	Seat Depth 19"	No Charge	
<input type="checkbox"/> I105045-99-0	Seat Depth 20"	No Charge	
<input type="checkbox"/> I105460-99-0	Seat Depth 21"	No Charge	
<input type="checkbox"/> I105047-99-0	Seat Depth 22"	No Charge	

CHOOSE SEAT CUSHION:

Part Number	Description	Price	HCPCS Code
	NOTE: <i>For seat widths of 23" and seat depths of 22", a Corpus II Ergonomic Seat Cushion will be provided.</i>		
<input type="checkbox"/> I10651	Corpus 3G Ergo Seat Cush - Leatherette (C3GERGSCA)	497.00	E2605
<input type="checkbox"/> I10652	Corpus 3G Ergo Seat Cush - Black Mesh (C3GERGSCA)	508.00	E2605
<input type="checkbox"/> I10653	Corpus 3G Ergo Seat Cush - Gray Mesh (C3GERGSCA)	508.00	E2605
<input type="checkbox"/> I105118-99-0	Omit Seat Cushion	No Charge	

Armrest Options

CHOOSE ARMBAR KIT AND MOUNTING POSITION:

C500 Corpus 3G

Armrest Options

CHOOSE ARMBAR KIT AND MOUNTING POSITION:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE:			
<i>Narrow Armbar Kit can be adjusted to 15"W, 17"W or 19"W between the armrest pads without any additional parts. Wide Armbar Kit can be adjusted to 19"W, 21"W, or 23"W between the armrest pads without any additional parts. Please select both Armbar Kit and Width Setting below.</i>			
<input type="checkbox"/> I105440-99-0	Narrow Armbar Kit, Adjustable 15", 17", 19"W	No Charge	
<input type="checkbox"/> I105441-99-0	Wide Armbar Kit, Adjustable 19", 21", 23"W	No Charge	
	Set Armbar at 15"W Between Armrests: _____		
	Set Armbar at 17"W Between Armrests: _____		
	Set Armbar at 19"W Between Armrests: _____		
	Set Armbar at 21"W Between Armrests: _____		
	Set Armbar at 23"W Between Armrests: _____		

CHOOSE ARMREST HEIGHT MOUNTING POSITION:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105081-99-0	Height Adjustable Amrests	298.00	
<input type="checkbox"/> I105069-99-0	Height Adjustable Amrests - NC	No Charge	
<i>No charge when ordering power tilt and/or recline.</i>			
SET ARMREST HEIGHT AT: _____			
<i>Corpus 3G Armrests are infinitely adjustable from 7" - 15", measured from the seat pan to the top of the armrest pad. No additional parts are needed to get this full range of adjustment. In one inch increments (ex: 8", 9", etc.), please indicate where you would like the armrests set.</i>			

CHOOSE ARMREST SWIVEL LOCKED OR ENABLED:

Part Number	Description	Price	HCPCS Code
ARMREST MEDIAL SWIVEL NOTE:			
<i>Every Corpus 3G Armrest can flip backwards for transfers. The armrests also include a feature that allows the armrest pads to swivel medially towards the center of the seating system. This swivel feature can be locked out by simply moving a bolt located under the armrest bracket. Please indicate below whether each armrest should have the swivel feature enabled or locked.</i>			
<input type="checkbox"/> I10645	Left Armrest Bracket - Medial Swivel Locked	No Charge	
<input type="checkbox"/> I10646	Right Armrest Bracket - Medial Swivel Locked	No Charge	
<input type="checkbox"/> I105073-99-0	Left Armrest Bracket - Medial Swivel Enabled	No Charge	
<input type="checkbox"/> I105072-99-0	Right Armrest Bracket - Medial Swivel Enabled	No Charge	

CHOOSE ARMREST PAD RIGHT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105100-99-0	4"x10" Corpus 3G Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105104-99-0	4"x13" Corpus 3G Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105105-99-0	4"x16" Corpus 3G Armrest Pad Right - Leatherette	No Charge	
<input type="checkbox"/> I105099-99-0	4"x10" Corpus 3G Armrest Pad Right - Black Mesh	No Charge	
<input type="checkbox"/> I105103-99-0	4"x13" Corpus 3G Armrest Pad Right - Black Mesh	No Charge	
<input type="checkbox"/> I105108-99-0	4"x16" Corpus 3G Armrest Pad Right - Black Mesh	No Charge	
<input type="checkbox"/> I105199-99-0	4"x10" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge	
<input type="checkbox"/> I105201-99-0	4"x13" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge	
<input type="checkbox"/> I105203-99-0	4"x16" Corpus 3G Armrest Pad Right - Gray Mesh	No Charge	

CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105098-99-0	4"x10" Corpus 3G Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105102-99-0	4"x13" Corpus 3G Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105107-99-0	4"x16" Corpus 3G Armrest Pad Left - Leatherette	No Charge	
<input type="checkbox"/> I105097-99-0	4"x10" Corpus 3G Armrest Pad Left - Black Mesh	No Charge	
<input type="checkbox"/> I105101-99-0	4"x13" Corpus 3G Armrest Pad Left - Black Mesh	No Charge	
<input type="checkbox"/> I105106-99-0	4"x16" Corpus 3G Armrest Pad Left - Black Mesh	No Charge	
<input type="checkbox"/> I105200-99-0	4"x10" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge	

C500 Corpus 3G

Armrest Options

CHOOSE ARMREST PAD LEFT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105202-99-0	4"x13" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge	
<input type="checkbox"/> I105204-99-0	4"x16" Corpus 3G Armrest Pad Left - Gray Mesh	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105463-99-0	Corpus 3G Arm Trough Adapter 1" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts.</i>	287.00	
<input type="checkbox"/> I105462-99-0	Corpus 3G Arm Trough Adpater 7/8" - Pair <i>For aftermarket armrest pads. Please indicate below if they should be installed or on chair sent as parts.</i>	287.00	
Install Arm Trough Adapters: _____ Send Arm Trough Adapters as Parts: _____			
<input type="checkbox"/> I100671-99-0	Armrest Pouch-Right	199.00	
<input type="checkbox"/> I100672-99-0	Armrest Pouch-Left	199.00	
<input type="checkbox"/> I10705	Corpus 3G Armrest Side Support - Right Outer <i>Custom Option. Leatherette Pad (2"T x 8.5"L).</i>	221.00	
<input type="checkbox"/> I10706	Corpus 3G Armrest Side Support - Left Outer <i>Custom Option. Leatherette Pad (2"T x 8.5"L).</i>	221.00	
<input type="checkbox"/> I10707	Corpus 3G Armrest Side Support - Right Inner <i>Custom Option. Leatherette Pad (2"T x 8.5"L). Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	221.00	
<input type="checkbox"/> I10708	Corpus 3G Armrest Side Support - Left Inner <i>Custom Option. Leatherette Pad (2"T x 8.5"L). Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	221.00	

Legrest Options

CHOOSE LEGREST TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105080-99-0	Corpus 3G Manual Adj Legrest Elevation	No Charge	
<input type="checkbox"/> I105079-99-0	Power Elevating Center Mount Legrests - 170° <i>Includes mechanical articulation.</i>	2,895.00	

FOOTPLATE HEIGHT POSITION:

Part Number	Description	Price	HCPCS Code
	SET FOOTPLATE HEIGHT AT: _____ <i>The height of the footplates are infinitely adjustable from 6" - 17". This is measured from the SEAT PAN to the top of the FOOTPLATE, when the legrest assembly is at 90°. DO NOT INLCUDE CUSHION THICKNESS IN YOUR MEASUREMENT. In one inch increments (ex: 12", 15", etc.), please indicate where you would like the footplates set. Calf Pads not available at footplate heights less than 9".</i>		

CHOOSE FOOTPLATE STYLE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105076-99-0	2pc Corpus 3G Footplates - 9"D x 7.5"W	No Charge	
<input type="checkbox"/> I105077-99-0	2pc Narrow Corpus 3G Footplates - 9"D x 6"W	221.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105252-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad.</i>	No Charge	
<input type="checkbox"/> I105258-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge	
<input type="checkbox"/> I105250-99-0	UT Calf Support Kit (A) w Mesh Cover & Adj Hardware <i>5.5"W x 6"T Pad.</i>	No Charge	
<input type="checkbox"/> I105256-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	327.00	

C500 Corpus 3G

Legrest Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105254-99-0	UT Calf Support Kit (D) w Mesh Cover & Adj Hardware 8"W x 8"T Pad.	327.00	

Accessories

MUST CHOOSE POSITIONING BELT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge	
<input type="checkbox"/> I10405	Bodypoint Padded Hip Belt UniTrack (N/C) - Medium	No Charge	
<input type="checkbox"/> I10402	Bodypoint Non-Padded Hip Belt, UniTrack - Large	97.00	E0978
<input type="checkbox"/> I10408	Bodypoint Padded Hip Belt UniTrack - Large	163.00	E0978
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	246.00	K0108
<input type="checkbox"/> I100600-99-0	Positioning Belt Retractable Long	267.00	K0108

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I105291-99-0	Headrest Corpus 3G Back, 10"W x 5"T, Leatherette <i>Must also select PN I10552.</i>	286.00	E0955
<input type="checkbox"/> I105290-99-0	Headrest Corpus 3G Back, 10"W x 5"T, Black Fabric <i>Must also select PN I10552.</i>	306.00	E0955
<input type="checkbox"/> I105212-99-0	Headrest Corpus 3G Back, Ergo Shape, Black Fabric <i>Must also select PN I10552.</i>	306.00	E0955
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	211.00	K0108
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	222.00	E1028
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) 7"W x 7"T Curved Pad. <i>Must also select PN I10553.</i>	131.00	E1020
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) 8"W x 8"T Flat Pad. <i>Must also select PN I10553.</i>	131.00	E1020
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Left: _____ Mount Amp Suppt on Right: _____</i>	222.00	E1028
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) 4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. <i>Must also select PN I10550.</i>	236.00	E0956
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) 6"W x 5"T Pad w/ 2" Tall Mounting Bracket. <i>Must also select PN I10550.</i>	236.00	E0956
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) 6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. <i>Must also select PN I10550.</i>	236.00	E0956
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) 8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. <i>Must also select PN I10550.</i>	236.00	E0956
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) 12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. <i>Must also select PN I10550.</i>	236.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	444.00	E1028
<input type="checkbox"/> I100978-99-0	Abductor Kit - Leatherette <i>Must also select PN I10555.</i>	320.00	E0957
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	222.00	E1028
<input type="checkbox"/> I10783	Hand Heater Kit with Buddy Button	822.00	
<input type="checkbox"/> I10784	One-Hand Cone Kit with PERMOfix Mounting Set <i>Kit includes Hand Heater (2"H x 1.5"W x 2.75"D), Buddy Button, Power Cable with XLR Plug and PERMOfix Mounting Adapter. Heater has three temp/fan speed settings. For best results, a hand cone (I10784 or I10785) is recommended to trap heat around the user's hand(s). To ensure best positioning for the user, heater will NOT be mounted when chair is received.</i>	507.00	
	<i>Clear plastic dome (5.5"W x 13.5"L x 8"H) with forearm opening at one end. Kit includes PERMOfix mounting hardware to attach cone to end of armrest. To ensure best positioning for the user, cone will NOT be mounted when wheelchair is received.</i>		

C500 Corpus 3G

Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10785	Two-Hand Cone Kit with Mounting Pins <i>Clear plastic triangular dome (15"W x 12.5"L x 8"H) with two forearm openings. Use with Upper Extremity Support is recommended. Kit includes four pins to secure cone to tray surface. Pins remove without tools for easy transfers. To ensure the best positioning for the user, cone will NOT be mounted.</i>	399.00	
<input type="checkbox"/> I105260-99-0	Transfer Handles, Corpus 3G Seat	409.00	
<input type="checkbox"/> I105292-99-0	Push Handles, Corpus 3G Seat	355.00	
<input type="checkbox"/> I103044-99-0	Cup/Cell Phone Holder	47.00	
<input type="checkbox"/> I105563-99-0	Medical Necess Bag Clips - Carabiner Style (C3G)	78.00	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	236.00	
<input type="checkbox"/> I101500-99-0	Essentials Carrier for Chassis <i>It is recommended that you also select Rear Tie Down Points Extension Kit (I101515-99-0) below when selecting this item.</i>	398.00	
<input type="checkbox"/> I101515-99-0	Rear Tie Down Points Extension Kit <i>Recommended when ordering Essentials Carrier for Chassis (I101500-99-0) in order to access rear tiedown points.</i>	118.00	
<input type="checkbox"/> I105144-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i>	149.00	E0960
<input type="checkbox"/> I105143-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i>	149.00	E0960
<input type="checkbox"/> I105145-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i>	149.00	E0960
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	249.00	E0950
<input type="checkbox"/> I105221-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	280.00	E0950
<input type="checkbox"/> I105222-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	293.00	E0950
<input type="checkbox"/> I105220-99-0	Upper Extremity Support Angle Adjustable. 18" W <i>Must also select PN I10556.</i>	320.00	E0950
<input type="checkbox"/> I105262-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	338.00	E0950
<input type="checkbox"/> I105263-99-0	Upper Extremity Support Joystick Cutout 18" W <i>Must also select PN I10556.</i>	338.00	E0950
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	310.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	222.00	E1028
<input type="checkbox"/> I105224-99-0	Oxygen Holder, Corpus 3G	235.00	E2208
<input type="checkbox"/> I105223-99-0	Crutch Holder, Corpus 3G	235.00	E2207

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10808	Stealth Ultra Headrest <i>Includes Stealth Occipital (500 series) and Sub Occipital (600 series) Pads, bilateral Facial Support Pads (#900) with Stealth Swing Away Hardware (SUS9-B) and Flip Down Headrest Mount (TWB480-FDM). Also includes Universal Headrest Adapter.</i>	1,295.00	E0955
<input type="checkbox"/> I10800	Stealth Swing Away Egg Switch for Ultra Headrest <i>Includes Black Egg Switch and Stealth Ultra Swing Away Hardware (SUS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i>	500.00	E1028

C500 Corpus 3G

Stealth Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10801	Stealth Swing Away Facial Supprt Hdw for Ultra <i>Must also select support pad (I10810, I10811, I10812). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	295.00	E1028
<input type="checkbox"/> I10810	Stealth #900 Support Pad (4.5" x 2") - Ultra	115.00	K0108
<input type="checkbox"/> I10811	Stealth #910 Support Pad (3" x 2") - Ultra	115.00	K0108
<input type="checkbox"/> I10812	Stealth #920 Support Pad (6.4" x 2") - Ultra	115.00	K0108
<input type="checkbox"/> I10718	Stealth Comfort Plus Headrest - 8"W (CP480) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	335.00	E0955
<input type="checkbox"/> I10719	Stealth Comfort Plus Headrest - 10"W (CP280) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	335.00	E0955
<input type="checkbox"/> I10720	Stealth Comfort Plus Headrest - 14"W (CP180) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	335.00	E0955
<input type="checkbox"/> I10721	Adj Stealth Headrest Hardware - Each (TWB480)	227.00	E1028
<input type="checkbox"/> I10736	Stealth Swng Away Egg Sw for Comfort Plus Headrest <i>Includes Black Egg Switch and Stealth Comfort Plus Swing Away Hardware (CPS9-R). Please indicate if you want switch mounted on left or right side of headrest: Mount Switch on Left: _____ Mount Switch on Right: _____</i>	500.00	E1028
<input type="checkbox"/> I10809	Stealth Swing Away Facial Supprt Hdw for Cmfrt Plus <i>Must also select support pad (I10813, I10814, I10815). Please indicate if you want it mounted on left or right side of headrest: Mount Hardware on Left: _____ Mount Hardware on Right: _____</i>	295.00	E1028
<input type="checkbox"/> I10813	Stealth #900 Support Pad (4.5" x 2") - Cmfrt Plus	115.00	K0108
<input type="checkbox"/> I10814	Stealth #910 Support Pad (3" x 2") - Cmfrt Plus	115.00	K0108
<input type="checkbox"/> I10815	Stealth #920 Support Pad (6.4" x 2") - Cmfrt Plus	115.00	K0108
<input type="checkbox"/> I10778	Stealth Lateral Supports 5"D x 3"T (TWBSTL-CS3G35) <i>Must also select I10723. Contains the THIN Profile Stealth Mounting Hardware, which allows for flipping back the armrest for transfers. Works with Corpus 3G backrest shells only.</i>	330.00	E0956
<input type="checkbox"/> I10735	Stealth Lateral Supports - 6"D x 4"T (TWBSTL-CS3G) <i>Must also select I10723. Contains the THIN Profile Stealth Mounting Hardware, which allows for flipping back the armrest for transfers. Works with Corpus 3G backrest shells only.</i>	330.00	E0956
<input type="checkbox"/> I10723	Adjustable Stealth Lateral Hardware (Pair)	469.00	E1028

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	2,590.00	E2312
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.5" (38mm) Foam Ball Knob.</i>	1,264.00	E2373
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off. Includes Standard Joystick Knob and 1.5" (38mm) Foam Ball Knob.</i>	1,446.00	E2373

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104690-99-0	Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	222.00	E1028
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	711.00	E1028
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	758.00	E1028

C500 Corpus 3G

Alternative Drive Systems

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 7½".</i>	539.00	E1028
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 5½".</i>	539.00	E1028
<input type="checkbox"/> I10754	Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 4".</i>	539.00	E1028
<input type="checkbox"/> I10068	Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i>	4,117.00	K0108
<input type="checkbox"/> I10795	Stealth Swing Away Chin Boom for Ultra Headrest <i>Includes Stealth Ultra Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest: Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____</i>	475.00	E1028
<input type="checkbox"/> I10797	Stealth Sw Away Chin Boom for Comfort Plus Headrest <i>Includes Stealth Comfort Plus Swing Away Hardware, 18" boom and bracket for attaching a Mini Joystick or Compact Joystick. Please indicate which bracket to include and if you want boom mounted on left or right side of headrest: Mini Joystick Bracket: _____ Compact Joystick Bracket: _____ Mount Boom on Left: _____ Mount Boom on Right: _____</i>	475.00	E1028
<input type="checkbox"/> I10807	Stealth Gatlin Eclipse Tray w/ Mini JS Mount <i>Designed to mount on Stealth Gatlin Mid Line Mounts. Includes bracket for mounting a Mini Joystick under tray. Through hole in tray surface provides joystick access. Please indicate side of tray for hole position: Position Hole on Left: _____ Position Hole on Right: _____</i>	575.00	E1028

ALT JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g) Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i>	226.00	K0108
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g)</i>	407.00	K0108

HEAD ARRAY SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10759	Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 5½" W x 5½" H. Requires Omni Input Module.</i>	4,300.00	E2330
<input type="checkbox"/> I10639	Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below: Mount on Left: _____ Mount on Right: _____</i>	500.00	E1028
<input type="checkbox"/> I10760	Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 5½" W x 5½" H.</i>	4,506.00	E2330
<input type="checkbox"/> I10761	Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i>	222.00	E1028

SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,519.00	E2325,E2326
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	2,022.00	E2325,E2326

C500 Corpus 3G

Alternative Drive Systems

SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	222.00	E1028

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10763	VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i>	2,568.00	E2328
<input type="checkbox"/> I10764	VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	222.00	E1028
<input type="checkbox"/> I10765	VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	711.00	E1028
<input type="checkbox"/> I10317	Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5" Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,506.00	E2328
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	2,724.00	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	222.00	E1028
<input type="checkbox"/> I10638	RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	5,021.00	E2328
<input type="checkbox"/> I10770	Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i>	222.00	E1028

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> <i>Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	1,541.00	
<input type="checkbox"/> I10577	Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i>	222.00	E1028
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,042.00	E2351

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Switchbox on Left	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Switchbox on Right	No Charge	

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	156.00	K0108
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	116.00	K0108
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	113.00	K0108
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	540.00	K0108
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	621.00	K0108
<input type="checkbox"/> I10771	Touch Contact with CR2032 Battery Qty: _____	243.00	K0108

C500 Corpus 3G

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____</i> <i>Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i>	430.00	K0108

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.