


C400 Corpus Junior



*** Indicates a Required Field**

Dealer Information

* **Contact:** _____
 Dealer Code: _____
 * **Dealer Name:** _____
 Address: _____
 * **City:** _____
 * **State/Zip:** _____
 * **Phone#:** _____
 * **Fax#:** _____
 PO#: _____
 Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* **First Name:** _____
 * **Last Name:** _____
 Diagnosis: _____
 Funding Source: _____
 Client Age: _____

Client Measurements

* **Weight:** _____
 * **Height:** _____
 A) Top of Shoulders: _____
 B) Chest Depth: _____
 C) Chest Width: _____
 D) Seat Depth: _____
 E) Top of Head: _____
 F) Elbow to Hand: _____
 G) Seat Pan to Elbow: _____
 H) Hip Width: _____
 * **I) Knee to Foot:** _____

Three column pricing on order/quote:
 Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**
 Email to: sales@permobilus.com

Permobil Inc.
 300 Duke Dr.
 Lebanon, TN 37090
 Tel: (800) 736-0925
 Fax: (800) 231-3256
www.permobil.com

Prices effective **August 29, 2011.**

C400 Corpus Junior

Chair Model and Colors

CHOOSE BASE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10475	C400 Base Corpus R-net - C400 PS0 <i>Price Includes: FWD C400 Power Base Set Up For No Power Seat Functions, Shock Absorbing Suspension System, R-net Basic Controller, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	9,395.00	K0848
<input type="checkbox"/> I10476	C400 Base Corpus R-net - C400 PS1 <i>Price Includes: FWD C400 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, R-net Controller, Tie Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	9,395.00	K0856
<input type="checkbox"/> I103259-99-0	C400 Base Corpus R-net - C400 PS2 <i>Price Includes: FWD C400 Power Base Set Up For Multiple Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below.</i>	9,395.00	K0861

CHOOSE DRIVE TIRES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103889-99-0	Tri Spoke Split Rim w/ FF Drive Tires - 2.5 x 8"	No Charge	
<input type="checkbox"/> I103890-99-0	Tri Spoke Split Rim w/ Pneumatic Tire - 2.5 x 8"	No Charge	

CHOOSE SHROUD COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102190-99-0	C400 Shroud Vintage Wine	No Charge	
<input type="checkbox"/> I103957-99-0	C400 Shroud Onyx Black	No Charge	
<input type="checkbox"/> I102194-99-0	C400 Shroud Crystal Blue	No Charge	
<input type="checkbox"/> I102191-99-0	C400 Shroud Metallic Charcoal Gray	No Charge	
<input type="checkbox"/> I102192-99-0	C400 Shroud Green Metallic	No Charge	
<input type="checkbox"/> I103958-99-0	C400 Shroud Extreme Yellow	No Charge	

CHOOSE HUBCAP COLOR:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10451	Match Hub Cap to Shroud Color	No Charge	
<input type="checkbox"/> I103537-99-0	Hub Cap Platinum Silver	No Charge	

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102085-99-0	Seat Tube Fixed	No Charge	
<input type="checkbox"/> I103865-99-0	Power Adjustable Seat Height, 8" Travel (R-net)	2,645.00	E2300

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IM24 SLD G FT	Batteries, Grp 24 (72Ah) Sealed LA/Gel, Installed	824.00	E2363
<input type="checkbox"/> I102210-99-0	Anti-Tippers for C400/C500	406.00	
<input type="checkbox"/> I103635-99-0	Light Kit, R-net, C400/C500	972.00	
<input checked="" type="checkbox"/> I100570-99-0	Seat Stabilizer	No Charge	
<input type="checkbox"/> I103241-99-0	Articulating Vent Tray for C400/C500 <i>If selecting this part, you MUST also pick: PN I104082-99-0 - ESP Unit. Fits an LP10 Sized Vent and a Group 34 Battery. Cannot be ordered with backrest height less than 26.5", a tubular back, or a 13" seat depth. Customization may be required if ordered with any of the following items Co-Pilot, any standard attendant control, push handles, or an oxygen holder. Only available with expandable controllers (R-net Electronics.)</i>	1,318.00	E1030
<input checked="" type="checkbox"/> I104082-99-0	R-net ESP w C400/C500 Mounting Bracket <i>MUST be selected if ordering any ventilator tray (I102769-99-0, I103241-99-0, or I102768-99-0). Recommended item with switched drive inputs for improved tracking.</i>	No Charge	K0108
<input type="checkbox"/> ID51113	Programming Key for R-net - Permobil Version	443.00	

Joystick Options

C400 Corpus Junior

Joystick Options

CHOOSE JOYSTICK TYPE:

Part Number	Description	Price	HCPCS Code
PLEASE NOTE:			
<i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i>			
<input type="checkbox"/> I10746	R-net Remote Color JS w/ Mono Jacks, No Lights, NC	No Charge	
<input type="checkbox"/> ID51185	R-net Remote Color JS w/ Mono Jacks, No Lights	695.00	
<input type="checkbox"/> ID51188	R-net Remote Color JS w/ Toggle Switches, No Lghts	890.00	
<input type="checkbox"/> I10747	R-net Remote Color JS w/ Mono Jacks, Lights, NC	No Charge	
<input type="checkbox"/> ID51153	R-net Remote Color JS w/ Mono Jacks, Lights	695.00	
<input type="checkbox"/> ID51149	R-net Remote Color JS w/ Toggle Switches, Lights	890.00	
<input type="checkbox"/> ID51190	R-net Remote Color JS w/ Spd Pot/Toggle, Lights	890.00	
<input type="checkbox"/> I103814-99-0	Compact Advanced Joystick for R-net	616.00	E2373
<i>Compact sized joystick that does not contain an LCD Display. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i>			
<input type="checkbox"/> I10621	Add Wrist Pad to Compact Joystick Advanced	198.00	

CHOOSE JOYSTICK HANDLE TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10045	Cone Shaped Joystick Knob - Standard.	No Charge	
<input type="checkbox"/> I100960-99-0	Chin Cup For Joystick Knob	96.00	E2324
<input type="checkbox"/> I100961-99-0	"T" Handle For Joystick Knob	141.00	E2323
<input type="checkbox"/> I100962-99-0	"Large Ball" For Joystick Knob	96.00	E2323
<input type="checkbox"/> I100963-99-0	"Softball" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100964-99-0	"Mushroom" For Joystick Knob	129.00	E2323
<input type="checkbox"/> I100965-99-0	"Stick" For Joystick Knob	129.00	E2323
<input type="checkbox"/> IPC101B	Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC102B	Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft	129.00	E2323
<input type="checkbox"/> IPC107B	Bodypoint J/S Handle Dome Shaped, Rubber	129.00	E2323

CHOOSE JOYSTICK MOUNT:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10412	Joystick Mount Right - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I10413	Joystick Mount Left - Fixed, R-net Only	No Charge	
<input type="checkbox"/> I102537-99-0	Joystick Mount Right - Swing Away, R-net & VR2	297.00	E1028
<input type="checkbox"/> I102538-99-0	Joystick Mount Left - Swing Away, R-net & VR2	297.00	E1028
<input type="checkbox"/> I103820-99-0	Retractable Joystick Mnt Right - R-net, VR2	495.00	E1028
<input type="checkbox"/> I103819-99-0	Retractable Joystick Mnt Left - R-net, VR2	495.00	E1028
<input type="checkbox"/> I101260-99-0	Midline Control Mount - Pilot+ & R-net	268.00	E1028

Please indicate which armrest to mount midline below.

Mount Midline Joystick on Left Armrest: _____

Mount Midline Joystick on Right Armrest: _____

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10411	Expandable Controller - R-net	545.00	E2377
<i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i>			
<input type="checkbox"/> I10513	Harness for Expandable Controller	475.00	E2313
<i>Must be ordered with any Expandable Controller (I10411 or I10450).</i>			
<input type="checkbox"/> I104003-99-0	Co-Pilot Attendant Control, Corpus R-net	4,248.00	E2331
<i>Custom charge may apply if ordering a Corpus Tubular Backrest OR a backrest height of 21" or Lower. Not available with VR2 Controllers.</i>			
<input type="checkbox"/> ID50882	Attendant Control for R-net Control System Only	1,103.00	E2331
<i>Not available with VR2 Controllers.</i>			

C400 Corpus Junior

Joystick Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51111	R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i>	809.00	
<input type="checkbox"/> I10582	Output Module (IOM) for R-net Control System <i>This allows a user to control an external device (ex: Communication Device, ECU, Door Opener, etc.) through the drive input (JSM, Head Array, etc.). One output module is required for each separate external device to be controlled. Up to three output modules can be installed on one wheelchair.</i>	1,012.00	E2351

Seating System

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
	The Seating Systems Listed Below Include: <i>Corpus Seat Frame with 1" Depth Adjustments from 13" - 18" (& up to 22" Deep with Standard Corpus Seat Frame) Back Frame with Adjustable Recline Angle, Reclining Armrest Assembly, Solid Seat Pan, Legrest Assembly with Adjustable Legrest Elevation, Weight Limit 165 lbs for Seat Depths Less Than 16". Weight Limit of 250 lbs for Seat Depths of 16" or Greater.</i>		
<input type="checkbox"/> I10491	Corpus Junior Seat w/o Pwr Tilt or Pwr Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	5,689.00	K0108
<input type="checkbox"/> I10492	Corpus Junior Seat w/ 45° Pwr Tilt & Manual Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	6,995.00	E1002
<input type="checkbox"/> I10493	Corpus Junior Seat w 150° Power Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	7,356.00	E1004
<input type="checkbox"/> I10494	Corpus Jr Seat w 45° Pwr Tilt & 150° Pwr Recline <i>Includes:Solid Backrest, Height Adjustable Armrest and Seat Sizes Adjustable in 1" Increments</i>	10,895.00	E1007

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10418	Single Seat Function Control Switch for R-net <i>Please choose either Push Buttons or Toggle Switches below.</i>	No Charge	
<input type="checkbox"/> I10419	Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	1,495.00	E2310
<input type="checkbox"/> I10585	Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i>	2,925.00	E2311

CHOOSE PUSH BUTTONS OR TOGGLES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I102798-99-0	Push Buttons on ICS Switchbox	No Charge	
<input type="checkbox"/> I102799-99-0	Toggle Switches on ICS Switchbox	249.00	
<input type="checkbox"/> I104621-99-0	Omit ICS Switchbox	No Charge	

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10626	Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i>	No Charge	
<input type="checkbox"/> I103471-99-0	Mount ICS Switchbox to End of Arm Pad - Corpus <i>ICS Switchbox will be mounted to the end of the armrest pad. If this item is selected, the ICS Switchbox will NOT swing away with the Joystick Module. Please indicate below if you would like it mounted onto the left or right arm pad.</i>	No Charge	

C400 Corpus Junior

Seating System Options

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

Part Number	Description	Price	HCPCS Code
	Mount ICS Switchbox on Left Armpad _____		
	Mount ICS Switchbox on Right Armpad _____		
<input type="checkbox"/> I103470-99-0	Mount ICS Switchbox on Backrest	127.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104196-99-0	ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i>	1,002.00	K0108
<input type="checkbox"/> I104746-99-0	Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i>	1,002.00	
<input type="checkbox"/> I10497	Memory Seat Program, C400 or C500 Corpus w R-net <i>If this is selected, the chair will be programmed to allow the user to store and save up to 3 different memory seat positions. On the C400/C500 Corpus and Lowrider, the following seat functions are capable of remembering a seat position: Power Tilt, Power Recline, & Power Legs. The Power Adj Seat Height Function is NOT capable of recalling a saved memory seat position. Please refer to the Corpus Seat Owner's Manual for more information on saving and recalling memory seat positions.</i>	100.00	

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
	BACKREST HEIGHT NOTE: <i>Backrest heights are measured from the seat pan to the top of the backrest shell with the backrest set at 90°. Ergo Backs include lateral wedges and lumbar pads.</i>		
<input type="checkbox"/> I101853-99-0	Ergo Back Fixed 14" W x 21" H - Leatherette	852.00	E2620
<input type="checkbox"/> I101863-99-0	Ergo Back Fixed 14"W x 26.5" High - Leatherette	852.00	E2620
<input type="checkbox"/> I101854-99-0	Ergo Back Sliding 16"W x 21"H - Leatherette	852.00	E2620
<input type="checkbox"/> I101861-99-0	Ergo Back Sliding 16"W x 26.5" H - Leatherette	852.00	E2620
<input type="checkbox"/> I10455	Ergo Back Fixed 14" W x 21" H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10456	Ergo Back Fixed 14"W x 26.5" High - Black Mesh	874.00	E2620
<input type="checkbox"/> I10457	Ergo Back Sliding 16"W x 21"H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10458	Ergo Back Sliding 16"W x 26.5" H - Black Mesh	874.00	E2620
<input type="checkbox"/> I10004	Rectangular Backrest 16"W - Plastic, Fixed <i>Choose backrest height from seat pan (19"-29"): _____ "</i>	1,139.00	
<input type="checkbox"/> I10240	Solid Back 16"W X 21"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>	No Charge	
<input type="checkbox"/> I10130	Solid Back 16"W x 26.5"H <i>Standard Option with Tilt and/or Recline Systems. Includes Leatherette Back Cushion ONLY.</i>	No Charge	
<input type="checkbox"/> I10111	Tubular Backrest Frame Corpus, 16" Wide <i>Backrest canes to mount aftermarket backrest systems. For tilt only systems. Measurement between the armrests is 17.5". Seat pan to top of canes is 22".</i>	363.00	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10014	Custom Backrest Height _____ " <i>Backrest height measured from seat pan and cannot be lower than 19". Price in addition to above backrest.</i>	636.00	
<input type="checkbox"/> I100972-99-0	Split Backrest <i>Split Backrest 8" from top to allow for transfers and/or transportation. This modification cannot be made on back heights lower than 26.5". Must select either PN: I10429 or I10430 below.</i>	544.00	
<input type="checkbox"/> I10429	Split Backrest Cushion	No Charge	

C400 Corpus Junior

Backrest Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10430	Do NOT Split Backrest Cushion	No Charge	
<input type="checkbox"/> I10284	Flat Backrest Cushion for Ergo Back - Leatherette <i>This style of cushion is ONLY available in Leatherette. Recommended cushion when ordering lateral supports. Lateral wedges will be omitted if ordering this part.</i>	91.00	
<input type="checkbox"/> I10318	ROHO Backrest Cushion For Ergo Back <i>Price Includes: ROHO MID PROFILE Backrest Cushion which will be matched to fit the Ergonomic Backrest that you select. Approximate height of MID PROFILE cells is 3" in their inflated state. Also includes a Mesh Cover for the ROHO Backrest Cushion, Hand Pump, and Patch Kit. This will NOT fit these part numbers: (I10004, I10005, I10111, I10112, I305510-00-90-0). Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>	636.00	
<input type="checkbox"/> I10319	ROHO Backrest Cushion For Ergo Back - Notched <i>Recommended ROHO cushion when ordering Lateral Trunk Supports. Please indicate below if you would like to OMIT the Ergo Back Cushion.</i>	636.00	
<input type="checkbox"/> I10320	Omit Ergonomic Backrest Cushion	No Charge	
<input type="checkbox"/> I10278	Notch Backrest Cushion for Lateral Supports <i>Backrest Cushion will be notched 2" in on either side to allow for closer adjustment of trunk supports.</i>	188.00	

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10264	Corpus Junior Seat Width - Standard <i>Corpus seat frame is 17" wide at the outside and 15" wide between the seat rails. Actual hip width can be accommodated with size appropriate cushion and/or hip guides.</i>	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10265	Corpus Junior Seat 13" Deep	No Charge	
<input type="checkbox"/> I10266	Corpus Junior Seat 14" Deep	No Charge	
<input type="checkbox"/> I10267	Corpus Junior Seat 15" Deep	No Charge	
<input type="checkbox"/> I10268	Corpus Junior Seat 16" Deep	No Charge	
<input type="checkbox"/> I10269	Corpus Junior Seat 17" Deep	No Charge	
<input type="checkbox"/> I10270	Corpus Junior Seat 18" Deep	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10276	Miniflex Seat Cushion Leatherette <i>Approximately 2" thick. Not available with seat widths or depths greater than 14".</i>	483.00	
<input type="checkbox"/> I10038	Ergonomic Seat Cushion Leatherette <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	483.00	E2601,E2602
<input type="checkbox"/> I10461	Ergonomic Seat Cushion - Black Mesh <i>Not available with seat widths less than 16" or seat depths less than 14".</i>	493.00	E2601

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103339-99-0	4"x10" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103341-99-0	4"x12" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103343-99-0	4"x14" Corpus Arm Right - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103331-99-0	4"x10" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103333-99-0	4"x12" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103335-99-0	4"x14" Corpus Arm Right - Dual Taper, Black Fabric	No Charge	

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103338-99-0	4"x10" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103340-99-0	4"x12" Corpus Arm Left - Dual Taper, Leatherette	No Charge	
<input type="checkbox"/> I103342-99-0	4"x14" Corpus Arm Left - Dual Taper, Leatherette	No Charge	

C400 Corpus Junior

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103330-99-0	4"x10" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103332-99-0	4"x12" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	
<input type="checkbox"/> I103334-99-0	4"x14" Corpus Arm Left - Dual Taper, Black Fabric	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10137	Height Adjustable Armrest Assembly	289.00	E0973
<input type="checkbox"/> I10016	Height Adjustable Armrest Assembly - No Charge <i>No charge with power tilt and/or recline.</i>	No Charge	
<input type="checkbox"/> I101003-99-0	Armrest Mount High <i>Range of adjustment in high position is 11" - 13.5" measured from the seat pan to the top of the armrest pad.</i>	No Charge	
<input type="checkbox"/> I101002-99-0	Armrest Mount Low (High Position Standard) <i>Range of adjustment in low position is 9" - 11.5" measured from the seat pan to the top of the armrest pad.</i>	No Charge	
<input type="checkbox"/> IARMADAP1"	Arm Trough Adaptor 1" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed or on chair sent as parts.</i>	279.00	
<input type="checkbox"/> IARMADAP7/8"	Arm Trough Adaptor 7/8" Pair <i>For aftermarket arm troughs. Please indicate below if they should be installed on chair or sent as parts.</i>	279.00	
Install Arm Trough Adapters: _____			
Send Arm Trough Adapters as Parts: _____			
<input type="checkbox"/> I100671-99-0	Armrest Pouch-Right	193.00	
<input type="checkbox"/> I100672-99-0	Armrest Pouch-Left	193.00	
<input type="checkbox"/> I101261-99-0	Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i>	144.00	
<input type="checkbox"/> I101262-99-0	Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i>	144.00	
<input type="checkbox"/> I101263-99-0	Armrest Side Support Right Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	144.00	
<input type="checkbox"/> I101264-99-0	Armrest Side Support Left Inner <i>Recommend 2" wider width between the armrest to ensure armrests will flip-up.</i>	144.00	
<input type="checkbox"/> I10285	Mount Corpus Armrest Pads On Inside of Bracket <i>Used to achieve 14"W between the arm measurement when ordering a 14"W Backrest.</i>	No Charge	

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103867-99-0	Power Articulating Center Mount Elevating Legrest	2,995.00	K0108
<input type="checkbox"/> I102212-99-0	Manual Adjustment Legrest Elevation <i>Includes Calf Supports.</i>	No Charge	

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I103055-99-0	UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less.</i>	No Charge	
<input type="checkbox"/> I102982-99-0	UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad.</i>	No Charge	
<input type="checkbox"/> I102978-99-0	UT Calf Support Kit (H) w Mesh Cover & Adj Hardware <i>7"W x 7"T Curved Pad.</i>	317.00	E1028
<input type="checkbox"/> I100561-99-0	Legrest Assm. 1pc Footplate (Complete)	402.00	
<input type="checkbox"/> I102051-99-0	Legrest Assm. 2pc Narrow Footplates (Complete),CII	No Charge	
<input type="checkbox"/> I10006	3" Width Between Two- Piece Footplates <i>Not available if chair has anti-tips and lower leg length is greater than 14" (measured from seat pan to top of footplate.)</i>	504.00	

C400 Corpus Junior

Legrest Options

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10271	Corpus Junior Lower Leg Length <i>Modification to allow for lower leg lengths of 8" - 10". Standard length is 11"-15". Measurement is taken from seat pan to the top of the footplates. Please specify lower leg length measurement in inches: _____</i>	No Charge	
<input type="checkbox"/> IFT202SR	BP Ankle Huggers, Med Side Release, Adj (Pair) <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	229.00	K0108

Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10400	Bodypoint Non-Padded Hip Belt, UniTrack - Medium	No Charge	
<input type="checkbox"/> I10405	Bodypoint Padded Hip Belt UniTrack (N/C) - Medium	No Charge	
<input type="checkbox"/> I10623	Bodypoint Padded Hip Belt UniTrack - Extra Small	115.00	
<input type="checkbox"/> I10402	Bodypoint Non-Padded Hip Belt, UniTrack - Large	94.00	K0108
<input type="checkbox"/> I10408	Bodypoint Padded Hip Belt UniTrack - Large	158.00	K0108
<input type="checkbox"/> I100519-99-0	Positioning Belt Retractable	239.00	K0108
<input type="checkbox"/> I100749-99-0	Miniflex Headrest Removable - Leatherette <i>Must also select PN I10552.</i>	268.00	E0955
<input type="checkbox"/> I104131-99-0	Headrest for Ergo Back w/ Curved Mount Brkt - Leath <i>Must also select PN I10552. This part consists of a leatherette headrest pad (10"W x 5"T) with a 10"L slotted, curved mounting bracket.</i>	268.00	E0955
<input type="checkbox"/> I103460-99-0	UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i>	297.00	E0955
<input type="checkbox"/> I1823249	UT Lateral Adjustment Bar for UniTrack Headrest <i>Allows for approximately 4" of lateral adjustment to the UT Ergo Headrest either right or left. NOTE: Due to the width of this item (approx 9"), it may be necessary to send other items such as back pack clips unmounted. Custom charge may apply with lower back heights.</i>	56.00	
<input type="checkbox"/> IHEADADAPTER	Universal Headrest Adapter <i>For aftermarket headrests.</i>	205.00	K0108
<input type="checkbox"/> I10552	Adj Removable - Headrest Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I103070-99-0	UT Amp. Support Pad (H) w/ Mesh Cover - (1 Each) <i>7"W x 7"T Curved Pad. Must also select PN I10553.</i>	127.00	E1020
<input type="checkbox"/> I103071-99-0	UT Amp. Support Pad (D) w/ Mesh Cover - (1 Each) <i>8"W x 8"T Flat Pad. Must also select PN I10553.</i>	127.00	E1020
<input type="checkbox"/> I10553	Adj Removable - Amp Support Hardware (1 Each) <i>Please indicate if you would like the amputation support mounted on the right or left side of the seat frame: Mount Amp Suppt on Right: _____ Mount Amp Suppt on Left: _____</i>	216.00	E1028
<input type="checkbox"/> I103077-99-0	UT Hip Suppt (B) w/ Mesh Cover - Low (Pair) <i>4.5"W x 3.5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I103079-99-0	UT Hip Suppt (A) w/ Mesh Cover - Low (Pair) <i>6"W x 5"T Pad w/ 2" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I103094-99-0	UT Thigh/Hip Suppt (J) w/ Mesh Cover - High (Pair) <i>6.5"W x 5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I103095-99-0	UT Thigh/Hip Suppt (C) w/ Mesh Cover - High (Pair) <i>8"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I102977-99-0	UT Thigh/Hip Suppt (G) w/ Mesh Cover - High (Pair) <i>12"W x 3.5"T Pad w/ 4" Tall Mounting Bracket. Must also select PN I10550.</i>	229.00	E0956
<input type="checkbox"/> I10550	Adj Removable - Thigh Support Hardware (Pair)	431.00	E1028
<input type="checkbox"/> I10272	Fixed Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	250.00	E0956
<input type="checkbox"/> I10274	Swing Away Trunk Supports Jr Corpus - Leatherette <i>Must also select PN I10551.</i>	270.00	E0956
<input type="checkbox"/> I10551	Adj Removable - Trunk Support Hardware (Pair)	431.00	E1028

C400 Corpus Junior

Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I100978-99-0	Abductor Kit - Leatherette <i>Must also select PN I10555.</i>	311.00	E0957
<input type="checkbox"/> I10555	Adj Removable - Abductor Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I10015	Knee Block Assembly Corpus	641.00	K0108
<input type="checkbox"/> I103651-99-0	Transfer Handles for UniTrack Accessory Rail	397.00	
<input type="checkbox"/> I100975-99-0	Push Handles for Ergo Backrest <i>Order with Ergonomic or Rectangular Backs. Custom charge may apply if ordering push handles with an articulating vent tray.</i>	345.00	
<input type="checkbox"/> I307103-00-0	Push Handles for Tubular Backrest <i>Order with Corpus Tubular Backs (PN: I10111 & I10112).</i>	343.00	
<input type="checkbox"/> I102059-99-0	Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i>	229.00	
<input type="checkbox"/> I101473-99-0	Medical Necessity Bag Clips - Carabiner Style	76.00	
<input type="checkbox"/> I101500-99-0	Essentials Carrier for Chassis <i>It is recommended that you also select Rear Tie Down Points Extension Kit (I101515-99-0) below when selecting this item.</i>	386.00	
<input type="checkbox"/> I101515-99-0	Rear Tie Down Points Extension Kit <i>Recommended when ordering Essentials Carrier for Chassis (I101500-99-0) in order to access rear tiedown points.</i>	115.00	
<input type="checkbox"/> I105144-99-0	Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> I105143-99-0	Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> I105145-99-0	Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i>	145.00	E0960
<input type="checkbox"/> ISH212	BP Shoulder Harness Medium, Front Pull <i>This item will be sent as a part with the wheelchair. It will not be installed.</i>	199.00	E0960
<input type="checkbox"/> I100578-99-0	Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i>	242.00	E0950
<input type="checkbox"/> I100575-99-0	Upper Extremity Support Angle Adjustable. 14" W <i>Must also select PN I10556.</i>	272.00	E0950
<input type="checkbox"/> I100576-99-0	Upper Extremity Support Angle Adjustable. 16" W <i>Must also select PN I10556.</i>	284.00	E0950
<input type="checkbox"/> I100571-99-0	Upper Extremity Support Joystick Cutout 16" W <i>Must also select PN I10556.</i>	328.00	E0950
<input type="checkbox"/> I100573-99-0	Upper Extremity Support Mid Joystick Cutout 16"W <i>Must also select PN I10556.</i>	402.00	E0950
<input type="checkbox"/> I103426-99-0	Upper Extremity Support Angle Adjust - UniTrack <i>Side mounted tray that attaches to the UniTrack Accessory Rail. Dimensions are: 19.5"W x 13"D. Can only be ordered with chairs that offer the UniTrack System. Must also select PN I10556.</i>	301.00	E0950
<input type="checkbox"/> I10556	Adj Removable - Tray Hardware (1 Each)	216.00	E1028
<input type="checkbox"/> I101964-99-0	Oxygen Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest. Customization may be required if ordering on a Street.</i>	228.00	E2208
<input type="checkbox"/> I101956-99-0	Crutch Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder. This item is mounted to the ergonomic backrest. Customization may be required if ordering on a Street.</i>	228.00	E2207
STEALTH OPTIONS			
<input type="checkbox"/> I10718	Stealth Comfort Plus Headrest - 8"W (CP480) <i>Must also select I10721. Includes Universal Headrest Adapter.</i>	325.00	E0955

C400 Corpus Junior

Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10721	Adj Stealth Headrest Hardware - Each (TWB485)	220.00	E1028
<input type="checkbox"/> I10736	Stealth Swng Away Egg Sw for Comfort Plus Headrest <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the right or the left side of the head array below:</i> Mount Swing Away Egg Switch on Right: _____ Mount Swing Away Egg Switch on Left: _____	485.00	E1028

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10751	Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i>	2,515.00	E2312
<input type="checkbox"/> I103893-99-0	Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off.</i>	1,227.00	E2373
<input type="checkbox"/> I103894-99-0	Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off.</i>	1,404.00	E2373

CHOOSE MOUNTING HARDWARE:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104690-99-0	Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	216.00	E1028
<input type="checkbox"/> I104695-99-0	Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	690.00	E1028
<input type="checkbox"/> I104143-99-0	Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i>	736.00	E1028
<input type="checkbox"/> I10752	Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 7½".</i>	523.00	E1028
<input type="checkbox"/> I10753	Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 5½".</i>	523.00	E1028
<input type="checkbox"/> I10754	Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 4".</i>	523.00	E1028
<input type="checkbox"/> I10068	Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i>	3,997.00	E1028

ALT JOYSTICK ACCESSORIES:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I104904-99-0	Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g)</i> Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____	219.00	K0108
<input type="checkbox"/> I104641-99-0	Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g)</i>	395.00	K0108

HEAD ARRAY SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10759	Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 5½" W x 5½" H. Requires Omni Input Module.</i>	4,175.00	E2330
<input type="checkbox"/> I10639	Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below:</i> Mount on Left: _____ Mount on Right: _____	485.00	E1028

C400 Corpus Junior

Alternative Drive Systems

HEAD ARRAY SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10760	Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 5½" W x 5½" H.</i>	4,375.00	E2330
<input type="checkbox"/> I10761	Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i>	216.00	E1028

SIP & PUFF SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10573	Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,475.00	E2325,E2326
<input type="checkbox"/> I10574	Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i>	1,963.00	E2325,E2326
<input type="checkbox"/> I10762	Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i>	216.00	E1028

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10763	VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i>	2,493.00	E2328
<input type="checkbox"/> I10764	VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i>	216.00	E1028
<input type="checkbox"/> I10765	VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i>	216.00	E1028
<input type="checkbox"/> I10317	Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5" Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i>	4,375.00	E2328
<input type="checkbox"/> I10767	Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i>	2,645.00	
<input type="checkbox"/> I10766	Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i>	690.00	E1028
<input type="checkbox"/> I10638	RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i>	4,875.00	E2328
<input type="checkbox"/> I10770	Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i>	216.00	E1028

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> ID51154IR	Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i>	1,496.00	
<input type="checkbox"/> I10577	Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i>	216.00	E1028
<input type="checkbox"/> I10581	Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i>	1,012.00	E2351

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> I10080	Mount Seat Control Switchbox on Backrest	No Charge	
<input type="checkbox"/> I10079	Mount Seat Control Switchbox on Left	No Charge	
<input type="checkbox"/> I10078	Mount Seat Control Switchbox on Right	No Charge	

C400 Corpus Junior

Alt Drive Accessories

OPTIONAL ITEMS:

Part Number	Description	Price	HCPCS Code
<input type="checkbox"/> IMICROSWITCH	Microlite Switch Qty: _____	151.00	K0108
<input type="checkbox"/> IBUDDYBUTTON	Buddy Button Qty: _____	113.00	K0108
<input type="checkbox"/> I5808	Mini Cup Switch Qty: _____	110.00	K0108
<input type="checkbox"/> IZTS-P	Proximity Switch Qty: _____	524.00	K0108
<input type="checkbox"/> IMINIJOY5SW	Mini Joystick 5 Switch OC	603.00	K0108
<input type="checkbox"/> I10771	Touch Contact with CR2032 Battery Qty: _____	236.00	K0108
<input type="checkbox"/> I10772	Dome Switch, Backrest Mnt w Gooseneck Qty: _____	417.00	K0108

Mount on Left: _____ Mount on Right: _____
Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.