

M300 PS



*** Indicates a Required Field**

Dealer Information

* Contact: _____

Dealer Code: _____

* Dealer Name: _____

Address: _____

* City: _____

* State/Zip: _____

* Phone#: _____

* Fax#: _____

PO#: _____

Email Address: _____

Client Information

Permobil recommends that the client is evaluated by a certified rehab specialist.

* First Name: _____

* Last Name: _____

Diagnosis: _____

Funding Source: _____

Client Age: _____

Client Measurements

* Weight: _____

* Height: _____

A) Top of Shoulders: _____

B) Chest Depth: _____

C) Chest Width: _____

D) Seat Depth: _____


E) Top of Head: _____

F) Elbow to Hand: _____

G) Seat Pan to Elbow: _____

H) Hip Width: _____

* I) Knee to Foot: _____



Three column pricing on order/quote:

Display HCPCS Codes on order/quote:

Please send order/quote to fax#: **(800) 231-3256**

Email to: sales@permobilus.com

Permobil Inc.

300 Duke Dr.

Lebanon, TN 37090

Tel: (800) 736-0925

Fax: (800) 231-3256

www.permobil.com

Prices effective **August 29, 2011.**

M300 PS

Chair Model and Colors

CHOOSE BASE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|----------|------------|
| <input type="checkbox"/> I10643 | M300 Base PS - M300 NPO <i>Price Includes: M300 Power Base Set Up For No Power Seat Functions, Shock Absorbing Suspension System, Tie-Down Hardware For Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 6.0mph, Weight Capacity = 300lbs, STFH = 17.5".</i> | 6,645.00 | K0848 |
| <input type="checkbox"/> I105230-99-0 | M300 Base PS - M300 SPO <i>Price Includes: M300 Power Base Set Up For One Power Seat Function, Shock Absorbing Suspension System, Tie-Down Hardware for Strap Systems, Charger, Flat Free Tires, and Standard Positioning Belt. Must Choose Flat Free or Pneumatic Tires Below. 6.0mph, Weight Capacity = 300lbs, STFH = 17.5".</i> | 6,645.00 | K0856 |
| <input type="checkbox"/> I10727 | 7.5mph High Speed Upgrade Package - M300 <i>Price Includes 7.5mph High Speed Motor Package. Available with R-net Electronics only.</i> | 1,895.00 | |

CHOOSE DRIVE TIRES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I103499-99-0 | Tri-Spoke Split Rim w/ Flat Free Drive Tire 3x8" | No Charge | |
| <input type="checkbox"/> I103498-99-0 | Tri-Spoke Split Rim w/ Pneumatic Tires - 3x8" | No Charge | |

CHOOSE SHROUD COLOR:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|-------------------------------|-----------|------------|
| <input type="checkbox"/> I105239-99-0 | M300 Shroud - Onyx Black | No Charge | |
| <input type="checkbox"/> I105240-99-0 | M300 Shroud - Metallic Green | No Charge | |
| <input type="checkbox"/> I105241-99-0 | M300 Shroud - Vintage Wine | No Charge | |
| <input type="checkbox"/> I105243-99-0 | M300 Shroud - Platinum Silver | No Charge | |
| <input type="checkbox"/> I105242-99-0 | M300 Shroud - Sublime Green | No Charge | |
| <input type="checkbox"/> I105244-99-0 | M300 Shroud - Titan Blue | No Charge | |

CHOOSE HUBCAP COLOR:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|-------------------------------|-----------|------------|
| <input type="checkbox"/> I10451 | Match Hub Cap to Shroud Color | No Charge | |
| <input type="checkbox"/> I103537-99-0 | Hub Cap Platinum Silver | No Charge | |

Base Options

MUST CHOOSE EITHER SEAT ELEVATOR OR FIXED SEAT TUBE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I101994-99-0 | Seat Tube Fixed | No Charge | |
| <input type="checkbox"/> I103779-99-0 | Power Adjustable Seat Height - 8" Travel | 2,445.00 | E2300 |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|--|--|-----------|------------|
| <input type="checkbox"/> IM24 SLD G FT | Batteries, Grp 24 (72Ah) Sealed LA/Gel, Installed | 824.00 | E2363 |
| <input type="checkbox"/> IM34 SLDG | Batteries, Grp 34 (60Ah) Sealed LA/Gel, Installed | 814.00 | K0108 |
| <input type="checkbox"/> I105428-99-0 | Enhanced Steering Perform Unit (ESP), R-net <i>Recommended item with switched drive inputs for improved tracking.</i> | 1,170.00 | K0108 |
| <input type="checkbox"/> I105297-99-0 | LED Light Kit for M300 / M400 w/ R-net Electronics | 995.00 | |
| <input type="checkbox"/> I105295-99-0 | Seat Stabilizer M300 PS with Tilt | No Charge | |
| <input type="checkbox"/> ID51113 | Programming Key for R-net - Permobil Version | 443.00 | |

Permolock

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|----------|------------|
| <input type="checkbox"/> I104680-99-0 | Permolock C3 Locking Base <i>Kit Includes: Permolock C3 Locking Base, Owner's Manual, Manual Release Handle, and Vehicle Mounted Release Button. Does NOT include Dual Locking Pin Kit (I10608). Locking Base must be installed by a NMEDA Dealer. Please provide the NMEDA Dealer and Address where the locking Base will ship to:</i> NAME: _____ ADDR: _____ | 1,806.00 | |

M300 PS

Permlock

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|--------|------------|
| <input type="checkbox"/> I10608 | Permlock C3 Dual Locking Pin Kit- C/K300, M-Series <i>This kit includes two spring loaded pins that will be mounted underneath the base of the wheelchair. These pins are spring loaded so they push into the wheelchair base when an obstacle is encountered. The kit also includes spacers that may be used if the Permlock C3 Base is being mounted in a vehicle with a sloped floor. The spacers are only needed for a C300 or K300.</i> | 605.00 | |

Joystick Options

CHOOSE JOYSTICK TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| | PLEASE NOTE: <i>For all remote joystick modules shown below, you must also select an expandable controller. For additional input devices, refer to Alternative Drive Systems Section.</i> | | |
| <input type="checkbox"/> ID51185 | R-net Remote Color JS w/ Mono Jacks, No Lights | 695.00 | |
| <input type="checkbox"/> ID51188 | R-net Remote Color JS w/ Toggle Switches, No Lghts | 890.00 | |
| <input type="checkbox"/> ID51153 | R-net Remote Color JS w/ Mono Jacks, Lights | 695.00 | |
| <input type="checkbox"/> ID51149 | R-net Remote Color JS w/ Toggle Switches, Lights | 890.00 | |
| <input type="checkbox"/> ID51190 | R-net Remote Color JS w/ Spd Pot/Toggle, Lights | 890.00 | |
| <input type="checkbox"/> I103814-99-0 | Compact Advanced Joystick for R-net <i>Compact sized joystick that does not contain an LCD Display. Must order expandable controller (PN: I10411) with this part. Contains two mono jacks for mode & on/off. Does NOT contain a wrist support pad. If you want a wrist support pad, please select I10621, below.</i> | 616.00 | E2373 |
| <input type="checkbox"/> ID51108 | R-net LED Remote Joystick Module <i>Standard Joystick Module. This Joystick does NOT have an LCD Display and does NOT contain mono jacks. Limited to five profiles with this Joystick Module. Please note if ordering this item, the ICS Switchbox must be mounted onto either the armrest pad or the backrest. Cannot order PN: I1822845 or I103469-99-0 with this style of joystick.</i> | No Charge | |

CHOOSE JOYSTICK HANDLE TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10045 | Cone Shaped Joystick Knob - Standard. | No Charge | |
| <input type="checkbox"/> I100960-99-0 | Chin Cup For Joystick Knob | 96.00 | E2324 |
| <input type="checkbox"/> I100961-99-0 | "T" Handle For Joystick Knob | 141.00 | E2323 |
| <input type="checkbox"/> I100962-99-0 | "Large Ball" For Joystick Knob | 96.00 | E2323 |
| <input type="checkbox"/> I100963-99-0 | "Softball" For Joystick Knob | 129.00 | E2323 |
| <input type="checkbox"/> I100964-99-0 | "Mushroom" For Joystick Knob | 129.00 | E2323 |
| <input type="checkbox"/> I100965-99-0 | "Stick" For Joystick Knob | 129.00 | E2323 |
| <input type="checkbox"/> IPC101B | Bodypoint J/S Handle 3" U-Shaped w/ Flex-Shaft | 129.00 | E2323 |
| <input type="checkbox"/> IPC102B | Bodypoint J/S Handle 4" U-Shaped w/ Flex-Shaft | 129.00 | E2323 |
| <input type="checkbox"/> IPC107B | Bodypoint J/S Handle Dome Shaped, Rubber | 129.00 | E2323 |
| <input type="checkbox"/> I10412 | Joystick Mount Right - Fixed, R-net Only | No Charge | |
| <input type="checkbox"/> I10413 | Joystick Mount Left - Fixed, R-net Only | No Charge | |
| <input type="checkbox"/> I102537-99-0 | Joystick Mount Right - Swing Away, R-net & VR2 | 297.00 | E1028 |
| <input type="checkbox"/> I102538-99-0 | Joystick Mount Left - Swing Away, R-net & VR2 | 297.00 | E1028 |
| <input type="checkbox"/> I103820-99-0 | Retractable Joystick Mnt Right - R-net, VR2 | 495.00 | E1028 |
| <input type="checkbox"/> I103819-99-0 | Retractable Joystick Mnt Left - R-net, VR2 | 495.00 | E1028 |
| <input type="checkbox"/> I101260-99-0 | Midline Control Mount - Pilot+ & R-net <i>Please indicate which armrest to mount midline below.</i> | 268.00 | E1028 |

Please indicate which armrest to mount midline below.

Mount Midline Joystick on Left Armrest: _____

Mount Midline Joystick on Right Armrest: _____

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|-------------|-------------|-------|------------|
|-------------|-------------|-------|------------|

M300 PS

Joystick Options

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|----------------------------------|--|----------|------------|
| <input type="checkbox"/> I10411 | Expandable Controller - R-net <i>This item is required when ordering power tilt, power recline, and power elevating legrests together OR when ordering an Alternative Drive Control.</i> | 545.00 | E2377 |
| <input type="checkbox"/> I10513 | Harness for Expandable Controller <i>Must be ordered with any Expandable Controller (I10411 or I10450).</i> | 475.00 | E2313 |
| <input type="checkbox"/> ID50882 | Attendant Control for R-net Control System Only <i>Not available with VR2 Controllers.</i> | 1,103.00 | E2331 |
| <input type="checkbox"/> ID51111 | R-net BT Mouse Module <i>The R-net BT Mouse Module allows for mouse control through a standard R-net Joystick or other type of input device (ex: Mini Joystick w/ Omni). It is only usable as a mouse on a computer. It CANNOT be paired to work with other Bluetooth devices such as cell phones, wireless headsets, etc. This kit includes: an R-net BT Mouse Module, a 1.5m R-net Bus Cable, and a USB Style Bluetooth Adapter for the user's computer. Note: The Bluetooth Adapter may not work with all computers.</i> | 809.00 | |

Seating System

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I104738-99-0 | PS Seating System - No Tilt | No Charge | |
| <input type="checkbox"/> I10724 | PS Seating System with 50° CG Power Tilt <i>Weight Limit 300lbs.</i> | 5,165.00 | E1002 |

Seating System Options

CHOOSE SEAT FUNCTION CONTROL TYPE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10289 | Single Seat Function Control Switch - LiberT <i>This device is a single toggle switch used to control only one seat function. If control of the seat function through the drive control is needed, then the drive electronics must be upgraded to Pilot+ or R-net Electronics, and a Single Seat Function Control Kit (Part Number I10129), must be selected.</i> | No Charge | |
| <input type="checkbox"/> I10419 | Single Seat Function Control Kit for R-net <i>For operation of a single seat function through the wheelchair electronics via the drive control and/or ICS Switchbox. Order with R-net Controller or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i> | 1,495.00 | E2310 |
| <input type="checkbox"/> I10585 | Multiple Seat Function Control Kit for R-net <i>MUST be selected if ordering two or more seat functions. This device is used to operate two or more seat functions through the wheelchair electronics via the drive control and/or ICS switch box. Order with R-net electronics or remote joystick options. Please indicate below if you want push buttons or toggles on the ICS Switchbox.</i> | 2,925.00 | E2311 |

CHOOSE PUSH BUTTONS OR TOGGLES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|----------------------------------|-----------|------------|
| <input type="checkbox"/> I102798-99-0 | Push Buttons on ICS Switchbox | No Charge | |
| <input type="checkbox"/> I102799-99-0 | Toggle Switches on ICS Switchbox | 249.00 | |
| <input type="checkbox"/> I104621-99-0 | Omit ICS Switchbox | No Charge | |

CHOOSE MOUNTING POSITION OF ICS SWITCHBOX:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10626 | Mount ICS Switchbox onto R-net Joystick Bracket <i>ICS Switchbox will be mounted to the same bracket as the R-net Joystick Module. If this item is selected, the ICS Switchbox will swing away with the Joystick Module.</i> | No Charge | |
| <input type="checkbox"/> I103470-99-0 | Mount ICS Switchbox on Backrest | 127.00 | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|----------|------------|
| <input type="checkbox"/> I104196-99-0 | ICS Alternative Switchbox for R-net Electronics <i>Not available with VR2 Electronics. Seat function interface box which accommodates up to 8 momentary switches. Provides an additional method to accessing the seat functions. You must still order the Multiple Seat Function Control Kit if ordering this part. 1/8" mini-mono switches not included.</i> | 1,002.00 | K0108 |

M300 PS

Seating System Options

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|----------|------------|
| <input type="checkbox"/> I104746-99-0 | Remote Stop System - R-net, Pilot+ <i>Consists of a Radio Frequency (RF) handheld transmitter and an RF receiver mounted on the wheelchair. Each transmitter and receiver is a serialized matched pair, meaning they must be used together. Custom charge may be required if ordering with a low backrest height or a pediatric wheelchair. If ordering Remote Stop System with a K450, push handles must also be selected.</i> | 1,002.00 | |

Backrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10144 | Sling Upholstery Back Support - Standard <i>Width will be the same as seat width.</i> | No Charge | |
| <input type="checkbox"/> I10145 | Curved Metal Back Support 22" High | 417.00 | |
| <input type="checkbox"/> I10146 | Curved Metal Back Support 22" High - No Charge <i>No charge if ordering tilt.</i> | No Charge | |
| <input type="checkbox"/> I10454 | Omit Backrest | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10148 | Angle Adjustable Back Canes <i>Angle adjustable from 85 to 115 degrees.</i> | 129.00 | |
| <input type="checkbox"/> I10149 | Angle Adjustable Back Canes - No Charge <i>No charge if ordering tilt.</i> | No Charge | |

Seat Sizing

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|----------------------------------|----------------|-----------|------------|
| <input type="checkbox"/> I10150A | Seat Width 16" | No Charge | |
| <input type="checkbox"/> I10151A | Seat Width 18" | No Charge | |
| <input type="checkbox"/> I10219A | Seat Width 20" | No Charge | |
| <input type="checkbox"/> I10220A | Seat Width 22" | No Charge | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|----------------|-----------|------------|
| <input type="checkbox"/> I10156 | Seat Depth 16" | No Charge | |
| <input type="checkbox"/> I10157 | Seat Depth 17" | No Charge | |
| <input type="checkbox"/> I10158 | Seat Depth 18" | No Charge | |
| <input type="checkbox"/> I10159 | Seat Depth 19" | No Charge | |
| <input type="checkbox"/> I10221 | Seat Depth 20" | No Charge | |
| <input type="checkbox"/> I10222 | Seat Depth 21" | No Charge | |
| <input type="checkbox"/> I10223 | Seat Depth 22" | No Charge | |

Armrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10163 | Height Adjustable Armrest Assembly PS <i>Adjustable in height, with a flip-back design.</i> | 289.00 | E0973 |
| <input type="checkbox"/> I10230 | Height Adjustable Armrest Assembly PS - No Charge <i>Adjustable in height, with a flip-back design. No charge if ordering power tilt.</i> | No Charge | |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10290 | PS Armrest Pads - 10.5"L x 2"W (Pair) | No Charge | |
| <input type="checkbox"/> I10291 | PS Armrest Pads - 14"L x 2"W (Pair) | No Charge | |
| <input type="checkbox"/> I10292 | Waterfall Armrest Pad for PS - 13.5"L x 3"W (Pair) | 80.00 | |
| <input type="checkbox"/> I10168 | Omit Armrest Pads | No Charge | |

M300 PS

Armrest Options

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> I10166 | Corpus Armrest Pad Adapter Brackets <i>Must be picked if ordering Corpus armrest pads.</i> | 51.00 | |
| <input type="checkbox"/> I10208 | Corpus Armrest Pad 4"x10" Left Leatherette | 94.00 | |
| <input type="checkbox"/> I10209 | Corpus Armrest Pad 4"x12" Left Leatherette | 94.00 | |
| <input type="checkbox"/> I10210 | Corpus Armrest Pad 4"x14" Left Leatherette | 94.00 | |
| <input type="checkbox"/> I10211 | Corpus Armrest Pad 4"x16" Left Leatherette | 94.00 | |
| <input type="checkbox"/> I10212 | Corpus Armrest Pad 4"x10" Right Leatherette | 94.00 | |
| <input type="checkbox"/> I10213 | Corpus Armrest Pad 4"x12" Right Leatherette | 94.00 | |
| <input type="checkbox"/> I10214 | Corpus Armrest Pad 4"x14" Right Leatherette | 94.00 | |
| <input type="checkbox"/> I10215 | Corpus Armrest Pad 4"x16" Right Leatherette | 94.00 | |
| <input type="checkbox"/> I10169 | Armrest Side Support Right Inner <i>Will limit the armrest from flipping back. For Corpus Armrest Pads only.</i> | 142.00 | |
| <input type="checkbox"/> I10170 | Armrest Side Support Left Inner <i>Will limit the armrests from flipping back. For Corpus Armrest Pads only.</i> | 142.00 | |
| <input type="checkbox"/> I101261-99-0 | Armrest Side Support Right Outer <i>For Corpus Armrest Pads only.</i> | 144.00 | |
| <input type="checkbox"/> I101262-99-0 | Armrest Side Support Left Outer <i>For Corpus Armrest Pads only.</i> | 144.00 | |

Legrest Options

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10172 | 70 Degree Swing Away Legrests (Pair) | 202.00 | |
| <input type="checkbox"/> I10231 | 70 Degree Swing Away Legrests (Pair) <i>No charge if ordering Power Tilt.</i> | No Charge | |
| <input type="checkbox"/> I10316 | Aeron Legrest Assembly for PS Seat <i>90° Legrest Hanger with Manual Adjustment Legrest Elevation. Price includes Corpus Style, Angle Adjustable, Two Piece Footplates. Range of adjustment from the seat pan to the top of the footplate is: 12 - 17" (with calf pads) and 6.5 - 17" (without calf pads).</i> | 321.00 | |

CHOOSE ONE OF THE FOLLOWING:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10173 | Footplates (Pair), Standard Small | No Charge | |
| <input type="checkbox"/> I10174 | Footplates (Pair), Standard Large <i>Not available with 16" seat width.</i> | No Charge | |
| <input type="checkbox"/> I10175 | Small Footplates, Multi - Axis Adjustable, Aluminum | 198.00 | K0108 |
| <input type="checkbox"/> I10176 | Large Footplates, Multi - Axis Adjustable, Aluminum <i>Cannot be ordered with 16" or 18" seat widths.</i> | 198.00 | K0108 |

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I104684-99-0 | UT Calf Support Kit (B) w Mesh Cover & Adj Hardware <i>3.5"W x 4.5"T Pad. Recommended pad to be used with lower leg lengths of 12" or less. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i> | No Charge | |
| <input type="checkbox"/> I104686-99-0 | UT Calf Support Kit (J) w Mesh Cover & Adj Hardware <i>5"W x 6.5"T Pad. Bracket designed for use with Aeron Center Mount Legrest Assembly. Cannot be ordered with Swing Away Legs (I10586 & I10587).</i> | No Charge | |
| <input type="checkbox"/> I10177 | Heel Loops <i>Not available with Aeron Legrest Assembly. (I10316).</i> | 92.00 | E0951 |
| <input type="checkbox"/> I10197 | Calf Supports Leatherette <i>Not available with lower leg lengths of 10" or less.</i> | No Charge | |

Accessories

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|-------------|-------------|-------|------------|
|-------------|-------------|-------|------------|

M300 PS

Accessories

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|-----------|------------|
| <input type="checkbox"/> I10401 | Bodypoint Non-Padded Hip Belt, PS - Medium | No Charge | |
| <input type="checkbox"/> I10406 | Bodypoint Padded Hip Belt PS (C300) - Medium | 126.00 | K0108 |
| <input type="checkbox"/> I10403 | Bodypoint Non-Padded Hip Belt, PS - Large | 95.00 | K0108 |
| <input type="checkbox"/> I10409 | Bodypoint Padded Hip Belt PS - Large | 158.00 | K0108 |
| <input type="checkbox"/> I10199 | Headrest-Removable Leatherette <i>Must also select PN I10552.</i> | 268.00 | E0955 |
| <input type="checkbox"/> I103460-99-0 | UT Ergo Headrest w/ Black Fabric Cover <i>Must also select PN I10552.</i> | 297.00 | E0955 |
| <input type="checkbox"/> I10180 | Universal Headrest Adapter <i>Will be sent as part if ordering sling back upholstery.</i> | 205.00 | K0108 |
| <input type="checkbox"/> I10552 | Adj Removable - Headrest Hardware (1 Each) | 216.00 | E1028 |
| <input type="checkbox"/> I102059-99-0 | Medical Necessities Bag <i>Includes Medical Necessity Bag Clips (Carabiner Style).</i> | 229.00 | |
| <input type="checkbox"/> I10179 | Medical Necessity Bag Clips - Carabiner Style <i>Will be sent as part if ordering sling back upholstery.</i> | 76.00 | |
| <input type="checkbox"/> I105144-99-0 | Bodypoint Monoflex Large Chest Support Belt Kit <i>21.5"L x 4.5"W Pad Size. Item will be shipped as parts.</i> | 145.00 | E0960 |
| <input type="checkbox"/> I105143-99-0 | Bodypoint Monoflex Medium Chest Support Belt Kit <i>18.25"L x 4"W Pad Size. Item will be shipped as parts.</i> | 145.00 | E0960 |
| <input type="checkbox"/> I105145-99-0 | Bodypoint Monoflex X-Lrge Chest Support Belt Kit <i>24.5"L x 5"W Pad Size. Item will be shipped as parts.</i> | 145.00 | E0960 |
| <input type="checkbox"/> I100578-99-0 | Upper Extremity Support Small, Angle Adjustable <i>Must also select PN I10556.</i> | 242.00 | E0950 |
| <input type="checkbox"/> I10556 | Adj Removable - Tray Hardware (1 Each) | 216.00 | E1028 |
| <input type="checkbox"/> I101964-99-0 | Oxygen Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Custom charge may be applied if ordering in combination with accessories such as medical necessities bags, attendant controls, push handles, back pack clips, etc. This item is mounted to the Ergonomic Backrest. Customization may be required if ordering on a Street.</i> | 228.00 | E2208 |
| <input type="checkbox"/> I101956-99-0 | Crutch Holder for PS or Corpus Seat <i>Cannot be ordered with a Corpus Tubular Backrest. (I10111 or I10112). Accessories such as medical necessities bags and push handles may need to be sent as parts if ordering a crutch holder. This item is mounted to the ergonomic backrest. Customization may be required if ordering on a Street.</i> | 228.00 | E2207 |

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL JOYSTICK SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|----------|------------|
| <input type="checkbox"/> I10751 | Mini-Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 0.5 oz (13 g) Neutral to Full Deflection = 0.2" (5mm)</i> | 2,515.00 | E2312 |
| <input type="checkbox"/> I103893-99-0 | Compact Joystick for Omni / IOM (9-pin connector) <i>Operating Force = 7.8 oz (220 g) Neutral to Full Deflection = 1.1" (28mm) Joystick includes two 1/8" mono jacks for Mode and On/Off.</i> | 1,227.00 | E2373 |
| <input type="checkbox"/> I103894-99-0 | Compact Joystick LITE for Omni / IOM (9-pin conn) <i>Operating Force = 1.8 oz (50 g) Neutral to Full Deflection = 0.7" (17mm) Joystick includes two 1/8" mono jacks for Mode and On/Off.</i> | 1,404.00 | E2373 |

CHOOSE MOUNTING HARDWARE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|--|--------|------------|
| <input type="checkbox"/> I104690-99-0 | Alt. JS Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> | 216.00 | E1028 |
| <input type="checkbox"/> I104695-99-0 | Alt. JS Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i> | 690.00 | E1028 |
| <input type="checkbox"/> I104143-99-0 | Alt. JS Adjustable Mount Kit - Chin Boom <i>Swing-Away / Removable Chin Boom mounted on right of backrest shell. Custom charge may be applied if ordering lower backrest heights.</i> | 736.00 | E1028 |

M300 PS

Alternative Drive Systems

CHOOSE MOUNTING HARDWARE:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|----------|------------|
| <input type="checkbox"/> I10752 | Alt. JS Adjustable Mount Kit - Haims Harness, Lrg <i>For neck width up to 7½".</i> | 523.00 | E1028 |
| <input type="checkbox"/> I10753 | Alt. JS Adjustable Mount Kit - Haims Harness, Sml <i>For neck width up to 5½".</i> | 523.00 | E1028 |
| <input type="checkbox"/> I10754 | Alt. JS Adjustable Mount Kit - Haims Harness, Peds <i>For neck width up to 4".</i> | 523.00 | E1028 |
| <input type="checkbox"/> I10068 | Alt. JS Adj. Mount Kit - Power SwingAway Chin Boom <i>Includes power swing away mount off the backrest with adjustable mounting hardware. Custom charge may be applied if ordering lower backrest heights.</i> | 3,997.00 | E1028 |

ALT JOYSTICK ACCESSORIES:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> I104904-99-0 | Dual Switches for Compact Joystick - Fixed Mount <i>Includes two pushbutton switches mounted on rigid arms. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g) Mount at front of Joystick: _____ Mount at rear of Joystick: _____ Mount at left of Joystick: _____ Mount at right of Joystick: _____</i> | 219.00 | K0108 |
| <input type="checkbox"/> I104641-99-0 | Dual Switches for Compact Joystick - Gooseneck Mnt <i>Includes two pushbutton switches mounted on adjustable gooseneck arms. Switches are positioned on Left and Right sides of Compact Joystick. Switches can be used for Mode & On/Off. Switch operation force = 4.5oz (130 g)</i> | 395.00 | K0108 |

HEAD ARRAY SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|----------|------------|
| <input type="checkbox"/> I10759 | Head Array - 3 Switch w Curved Wings for Omni <i>Center Pad measures 5½" W x 5½" H. Requires Omni Input Module.</i> | 4,175.00 | E2330 |
| <input type="checkbox"/> I10639 | Stealth Swing Away Egg Switch for 3SW Head Array <i>Includes: Black Egg Switch with Swing Away Stealth Mount. Please indicate if you want it mounted on the left or the right side of the head array below: Mount on Left: _____ Mount on Right: _____</i> | 485.00 | E1028 |
| <input type="checkbox"/> I10760 | Head Array - 5 Switch w Flat Wings for Omni / IOM <i>Center Pad measures 5½" W x 5½" H.</i> | 4,375.00 | E2330 |
| <input type="checkbox"/> I10761 | Adjustable Mounting Hrdw - Head Array <i>Required when selecting PN I10759 or I10760.</i> | 216.00 | E1028 |

SIP & PUFF SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|---|----------|-------------|
| <input type="checkbox"/> I10573 | Sip & Puff System for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i> | 1,475.00 | E2325,E2326 |
| <input type="checkbox"/> I10574 | Sip & Puff System with Fiber Optic Opti-Stop <i>Requires Omni Control Module. Cannot use this part with an Input Module. Includes breath tube kit.</i> | 1,963.00 | E2325,E2326 |
| <input type="checkbox"/> I10762 | Adjustable Mounting Hrdw - Sip & Puff <i>Required when selecting PN I10573 or I10574.</i> | 216.00 | E1028 |

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|----------|------------|
| <input type="checkbox"/> I10763 | VIC Proportional Finger Control for Omni / IOM <i>The VIC CCD translates small finger movements into proportional driving commands.</i> | 2,493.00 | E2328 |
| <input type="checkbox"/> I10764 | VIC Adjustable Mount Kit - Armrest <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____</i> | 216.00 | E1028 |
| <input type="checkbox"/> I10765 | VIC Adjustable Mount Kit - Midline <i>Mount on Left Arm: _____ Mount on Right Arm: _____</i> | 216.00 | E1028 |
| <input type="checkbox"/> I10317 | Proportional Tablet Control for Omni / IOM <i>The Tablet Control is a unique drive input which enables a user to control a power wheelchair in a very similar manner to a mouse being controlled by a touchpad on a laptop computer. The dimensions of the Tablet Control Pad are: 12"W x 9.5"D x 1.5" Kit also includes: Buddy Button Switch, and Tablet Control Mounting Bracket.</i> | 4,375.00 | E2328 |

M300 PS

Alternative Drive Systems

ALTERNATIVE PROPORTIONAL CONTROL SYSTEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|----------|------------|
| <input type="checkbox"/> I10767 | Proportional Foot Control <i>Mount on Left: _____ Mount on Right: _____</i> | 2,645.00 | |
| <input type="checkbox"/> I10766 | Adjustable Mounting Hrdw - Foot Control <i>Required when selecting PN I10767.</i> | 690.00 | E1028 |
| <input type="checkbox"/> I10638 | RIM Control - SwitchIT for Omni <i>Requires Omni Control Module. Cannot use this part with an Input Module.</i> | 4,875.00 | E2328 |
| <input type="checkbox"/> I10770 | Adjustable/Removable Hardware - RIM Control <i>Required when selecting PN I10638.</i> | 216.00 | E1028 |

Alt Drive Systems Setup

CHOOSE OMNI OR INPUT MODULE (REQUIRED):

| Part Number | Description | Price | HCPCS Code |
|------------------------------------|---|----------|------------|
| <input type="checkbox"/> ID51154IR | Omni Control Module for R-net w/ Infrared <i>Mount in front of Left Arm: _____ Mount in front of Right Arm: _____ Please specify which type of alternative drive input you will be using if not ordering the input device from Permobil. Input Device: _____</i> | 1,496.00 | |
| <input type="checkbox"/> I10577 | Adjustable Swing-Away Hardware - Omni <i>Required when selecting Part Number ID51154IR.</i> | 216.00 | E1028 |
| <input type="checkbox"/> I10581 | Input Module (IOM) for R-net Control System <i>This device allows a proportional or switched drive input to be connected to an R-net Control System without the use of an Omni. It must be ordered in conjunction with an R-net Remote Joystick. It does NOT support any of the following input devices: Sip-n-Puff, 3-Switch Head Array or RIM control.</i> | 1,012.00 | E2351 |

CHOOSE MOUNTING POSITION OF SEAT CONTROL SWITCHBOX:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------|--|-----------|------------|
| <input type="checkbox"/> I10080 | Mount Seat Control Switchbox on Backrest | No Charge | |
| <input type="checkbox"/> I10079 | Mount Seat Control Switchbox on Left | No Charge | |
| <input type="checkbox"/> I10078 | Mount Seat Control Switchbox on Right | No Charge | |

Alt Drive Accessories

OPTIONAL ITEMS:

| Part Number | Description | Price | HCPCS Code |
|---------------------------------------|---|--------|------------|
| <input type="checkbox"/> IMICROSWITCH | Microlite Switch Qty: _____ | 151.00 | K0108 |
| <input type="checkbox"/> IBUDDYBUTTON | Buddy Button Qty: _____ | 113.00 | K0108 |
| <input type="checkbox"/> I5808 | Mini Cup Switch Qty: _____ | 110.00 | K0108 |
| <input type="checkbox"/> IZTS-P | Proximity Switch Qty: _____ | 524.00 | K0108 |
| <input type="checkbox"/> IMINIJOY5SW | Mini Joystick 5 Switch OC | 603.00 | K0108 |
| <input type="checkbox"/> I10771 | Touch Contact with CR2032 Battery Qty: _____ | 236.00 | K0108 |
| <input type="checkbox"/> I10772 | Dome Switch, Backrest Mnt w Gooseneck Qty: _____ <i>Mount on Left: _____ Mount on Right: _____ Mounts at top of backrest shell. Custom charge may be applied if ordering with short backrest heights Kit Includes: One PERMOfix joint with Flat Plate, two PERMOfix joints with 45° swivel adjustment, one 6" rod, one 8" rod and one 9" gooseneck switch with 50" cable and mounting clamp.</i> | 417.00 | K0108 |

Notes & Additional Instructions:

Chair Order Policy:

You are highly encouraged to contact Permobil Sales with ANY questions regarding proper completion of our order forms at 1-800-736-0925. The order cannot be processed without the following:

- A client name or code, the client's measurements, the client's weight
- The order form must also be signed by a person authorized by your company to acknowledge the items selected. Permobil, Inc. is not responsible for configuration or size discrepancies resulting from customer errors on the order form.
- Order confirmations will be provided by Permobil summarizing the items selected. Please review this carefully.
- Once the order is shipped, any changes are subject to the returns policy and may be prohibited.

Order Acknowledgement:

I, _____, am an agent of the medical equipment provider named on this order form and I have the authority to contract for the purchase of powered wheelchairs and related parts on behalf of said provider. I acknowledge that I have reviewed this order and that it is complete and accurate to the best of my knowledge. I further acknowledge that any changes to this order after submission of this order form are subject to the returns policy and may carry additional charges.

Signed: _____

Printed Name: _____

Title: _____

Date: _____

NOTE: All specifications and prices are subject to change without notice. Please note that prices displayed are only valid if ordered with the wheelchair. Any parts ordered for service or future alterations will carry different pricing and warranty. MID PROFILE® Cushion is a registered trademark of the The ROHO Group. *The HCPCS codes provided are not intended to be billing or legal advice, rather our interpretation of the code definitions. Use of the codes does not ensure coverage or payment for the item. For coverage information, verify the policy of the appropriate payer.*