

permobil



## KNOWLEDGE IS POWER

The rear-wheel drive C350 PS has two different speed packages of up to 6.5 MPH, ideal for quick class changes. The C350 PS offers easy access to the battery as well as a seat-to-floor height of 18.25." Available in Black Onyx, Crystal Blue, Extreme Yellow, Platinum Silver and Vintage Wine, the C350 PS comes with standard matching tri-spoke wheel covers to help complete the look. Plus, the PS seating system features a sliding, telescoping frame and expanding seat pan to allow for growth. The adjustability of the PS seating system is unmatched with 5° at the footplate, armrests that both flip up and move out of the way, and footplate hangers that swing away or inside. The C350 PS is designed to support aftermarket options and accessories, a smart choice for people with hectic schedules and active lives.

**C350**  
PS

**POWER FUNCTIONS**



Power adjustable tilt



VR2 electronics  
standard with PS seating system



4-wheel adjustable independent suspension




Power seat elevator

\*shown with optional after-market seat cushion




Simple flip-up and removeable armrest


SHROUD COLORS



Crystal Blue



Vintage Wine



Platinum Silver



Xtreme Yellow



Black Onyx

**SPECIFICATIONS**

Maximum Speed	6.5 MPH	Power Seat Elevator	8"
Driving Range	15-22 Miles	Seat Widths	16" - 22" (Adjustable)
Max User Weight	300 lbs	Seat Depths	16" - 22" (Adjustable)
Ground Clearance	3.1"	Backrest Heights	22"
Battery Capacity	Group 34	Seat To Floor Height	18.25"
Width	25"	Electronics (standard)	VR2 Basic RNet electronics optional
Length w/o Front Riggings	36.5"	Alternative Drive Controls	Yes w/R-Net
Recommended Battery	Group 34	Effect of Driving	1/2 Speed @ 135° back angle
Gradient Capacity	10° Up 6° Dn		



THE POWER OF MOBILITY

**Permobil Inc.**  
**300 Duke Drive**  
**Lebanon, Tennessee 37090**  
**800.736.0925**  
**www.permobil.com**