



READY. SET. ROLL!

With Permobility's K300 PS Junior, cool kids of all sizes can get in on the action. The K300 comes with a standard adult power base, and the expandable seat pan allows for a custom fit no matter how fast a child grows. Features such as EZ adjust armrests, the power seat elevator and power tilt functions, R-Net electronics and the Uni-track system provide for maximum comfort and allow for an infinite number of adjustments and add-on accessories. Plus, the hot new colors—Chili Pepper Red, Black Onyx, Platinum Silver, Popstar Pink, Metallic Green and Titan Blue—and the tri-spoke wheel covers make any child the envy of the block.

K300
PS Jr.



POWER FUNCTIONS



Power seat elevator and power tilt



R-Net electronics for improved programmability

SHROUD COLORS



Crystal Blue



Platinum Silver



Black Onyx



Chili Pepper Red



Metallic Green



Popstar Pink



Titan Blue



Adjustable seat pan fits kids of all sizes



Push handles



Uni-track system allows for extensive customization



One-step removable armrest



EZ adjust armrests easily adjust height angle and rotation



Easy touch swing-away footrests with adjustable high strength multi-axis footplates. Available with 70° or 90° front hangers.



Optional stroller handle and manual tilt feature with tool-less height adjustability that fold down for transport



SPECIFICATIONS

Maximum Speed	5 MPH	Power Seat Elevator	8"
Max User Weight	165 lbs	Seat Widths	11" - 16" (Adjustable) in 1" increments
Ground Clearance	3"	Seat Depths	10" - 18" (Adjustable)
Battery Capacity	Group 34	Backrest Heights	16" - 22"
Minimum Width	24.25"	Backrest Angle Adjustments	85° - 120°
Base Length	34.5"	Seat To Floor Height w/ Tilt	16.5"
Battery Charger Capacity	8A	Tilt Angle	45° Power/Manual
Weight of Wheelchair w/Batteries	265 lbs	Electronics (standard)	R-Net
Uphill Angle	10°	Alternative Drive Controls	Yes w/R-Net
Downhill Angle	6°	Effect of Driving When Tilted	1/2 Speed @ 20°
Side Slope	10°		



THE POWER OF MOBILITY

Permobil Inc.
300 Duke Drive
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com