



READY OR NOT...

Designed to grow with a child for years to come, the Koala Miniflex is the ideal pediatric wheelchair for toddlers. While it is scaled to fit the needs of children, the Koala is built with the same rugged motor and front-wheel-drive technology as Permobil's full-size adult chairs. The low seat-to-floor height allows little users to interact at the same height as their peers, and the integrated elevation provides them with easy access to a variety of surface heights such as tables and desks. Plus, the tilt feature offers improved stability and rest, a necessity after a long day of play.

Koala
MINIFLEX



POWER FUNCTIONS




Power tilt

SHROUD COLORS



Yellow



Red

* Custom shroud colors available for additional charge.



Pilot+™ Series programmable electronics are standard



8" power seat elevator



Remote stopping device allows for parental control



Optional bumper — an ideal safety accommodation for beginner drivers



Available with leatherette or fabric upholstery

SPECIFICATIONS

Maximum Speed	4 MPH	Gradient Capacity	10° Up 6° Dn
Driving Range	15-24 Miles	Power Seat Elevator	8"
Max User Weight	110 lbs	Seat sizes available	10" x 10" / 12" x 12" / 14" x 14"
Ground Clearance	2.5"	Seat To Floor Height	15"
Width	22"	Electronics	Pilot+
Length	28"-32"	Alternative Drive Controls	Yes
Recommended Battery	M40 Gell Cell M40-12 SLDG		



THE POWER OF MOBILITY

Permobil Inc.
300 Duke Drive
Lebanon, Tennessee 37090
800.736.0925
www.permobil.com